



PolyStain DS Kit - for 2 Mouse antibody on Human tissue

(DAB/Fast-Red)

NB-23-00100- 3(120 ml)

NB-23-00100- 2(36 ml)

NB-23-00100- 1(12 ml)



PolyStain DS Kit - for 2 Mouse antibody on Human tissue (DAB/Fast-Red)

NB-23-00100-1; NB-23-00100-2; NB-23-00100-3

Storage: 2-8°C

INTENDED USE:

The PolyStain DS Kit is designed to use with user supplied two mouse antibodies to detect two distinct antigens on human tissue or cell samples. This kit has been tested in paraffin tissue. However, this kit can be used on frozen specimen and freshly prepared monolayer cell smears. Double staining is one of most common methods used in immunohistostaining that allow revealing two distinct antigens in a single tissue. PolyStain DS Kit from NeoBiotech Labs supplies two polymer enzyme conjugates: HRP polymer anti-Mouse IgG and AP polymer anti-Mouse IgG with two distinct substrates/chromogens, DAB (brown color, react with HRP polymer anti-Mouse IgG) and Fast Red (red color, react with AP polymer anti-Mouse IgG). PolyStain DS Kit is a non-biotin system that avoids endogenous biotin non-specific binding.

KIT COMPONENTS:

Component No.	Content	12mL Kit	36mL Kit	120mL Kit
Reagent 1	HRP Polymer anti-Mouse (RTU)	6ml	18ml	60ml
Reagent 2A	DAB Substrate (RTU)	12ml	18ml	70ml
Reagent 2B	DAB Chromogen (20x)	1.5ml	2ml	3ml
Reagent 3A	DS-MM Blocker A (RTU)	6ml	18ml	60ml
Reagent 3B	DS-MM Blocker B (RTU)	6ml	18ml	60ml
Reagent 4	AP Polymer anti-Mouse (RTU)	6ml	18ml	60ml
Reagent 5A	Fast Red Tablet	6 tablets	18 tablets	60 tablets
Reagent 5B	Fast Red Substrate (RTU)	5ml x 6	5ml x 18	5ml x 60
Reagent 6	NeoMount Universal (RTU)	7ml	18ml	70ml

RECOMMENDED PROTOCOL:

1. Fixation: To ensure the quality of the staining and obtain reproducible performance, user needs to supply appropriately fixed tissue and well prepared slides.
2. Tissue needs to be properly adhered to the slide tightly to avoid falling off.
3. Paraffin embedded section must be deparaffinized with xylene and rehydrated with a graded series of ethanol before staining.
4. Cell smear samples should be made as much monolayer as possible to obtain satisfactory results.
5. Three control slides will aid the interpretation of the result: positive tissue control, reagent control (slides treated with Isotype control reagent), and negative control.
6. It takes about 30 minutes to dissolve Fast Red tablet into the substrate buffer. Make sure to start preparing Fast Red solution near the end of the secondary antibody incubation.
7. Proceed IHC staining: **DO NOT** let specimen or tissue dry from this point on.

Reagent	Staining Procedure	Incubation Time (Min.)
1. Peroxidase Blocking Reagent Not provided	<ol style="list-style-type: none"> a. Incubate slides in peroxidase blocking reagent (Ready-to-use 3% H₂O₂ solution) for 10 minutes. b. Rinse the slide using distilled water. 	10 min.
2. HIER Pretreatment: Refer to antibody data sheet.	<ol style="list-style-type: none"> a. Heat Induced Epitope Retrieval (HIER) may be required for primary antibody suggested by vendor. b. Wash with PBS for 2 min., 3 times. 	
3. Preblock (optional)	For paraffin section, Improved formula saves the need for a preblock step. For frozen tissue, preblock may or may not be required depending on fixative. (Preblock catalogue No.NB-23-00169 was Recommended.)	
4. Mouse Antibody 1: Supplied by user	<p>Notes: Investigator needs to optimize dilution and incubation times prior to double staining.</p> <ol style="list-style-type: none"> a. Apply 2 drops or enough volume of mouse primary antibody 1 to cover the tissue completely. Incubate in moist chamber for 30-60 min. b. Rinse with PBS containing 0.05% Tween-20 for 2 min., 3 times. 	30 - 60 min.

<p>5. Reagent 1: HRP Polymer anti-Mouse (RTU)</p>	<p>a. Apply 1drop (50µl) of Reagent 1 HRP Polymer anti-Mouse to cover each section. b. Incubate in moist chamber for 15-30 min. c. Rinse with PBS containing 0.05% Tween-20 for 2 min.,3 times</p>	<p>15 – 30 min</p>
<p>6. Reagents 2A, 2B: Reagents 2A: DAB Substrate(RTU) Reagents 2B: DAB Chromogen (20x)</p>	<p>a. Add 1 drop or 2 drops (for higher sensitivity and contrast) of Reagent 2B to 1 ml Reagent 2A. Mix well. Protect from light and use within 5 hours. b. Apply 2 drops or enough volume of DAB chromogen working solution to completely cover tissue. Incubate for 3-10 min. c. Rinse thoroughly with distilled water 4 times, 2 minutes each time.</p>	<p>3 - 10 min</p>
<p>7. Reagent 3A: DS-MM Blocker A</p>	<p>a. Apply 2 drops or enough volume of Reagent 3A DS-MM Blocker A to cover the tissue completely. Mix well on the slide and Incubate in moist chamber for 30 min. b. Rinse with PBS containing 0.05% Tween-20 for 2 min., 3 times.</p>	<p>30 min</p>
<p>8. Reagent 3B: DS-MM Blocker B</p>	<p>a. Apply 2 drops or enough volume of Reagent 3B DS-MM Blocker B to cover the tissue completely. Mix well on the slide and Incubate in moist chamber for 5 min. b. Rinse with PBS containing 0.05% Tween-20 for 2 min., 3 times.</p>	<p>5 min</p>
<p>9. Mouse antibody 2: Supplied by user</p>	<p>Notes: Investigator needs to optimize dilution and incubation times prior to double staining. a. Apply 2 drops or enough volume of mouse primary antibody 2 to cover the tissue completely. b. Rinse with PBS containing 0.05% Tween-20 for 2 min., 3 times.</p>	<p>30 - 60 min.</p>
<p>10. Reagent 4: AP Polymer anti-Mouse (RTU)</p>	<p>a. Apply 1drop (50µl) of Reagent 4 AP Polymer anti-Mouse to cover each section. b. Incubate in moist chamber for 15-30 min. c. Rinse with PBS containing 0.05% Tween-20 for 2 min., 3 times. d. Rinse well with tap water.</p>	<p>15 - 30 min.</p>
<p>11. Reagent 5A, 5B: Reagent 5A: Fast Red Tablet Reagent 5B: Fast Red Substrate(RTU)</p>	<p>Notes: Tablet dissolves in about 30 minutes. Prepare in advance. a. Dissolve 1 Reagent 5A Fast Red tablet in 5ml Reagent 5B Fast Red Substrate, vortex until the tablet dissolved completely. Use within 1 hour. b. Apply 2 drops (100µl) or enough volume of Fast -Red working solution to completely cover the tissue. Incubate for 10-20 min, observe appropriate color development. c. Rinse well with distilled water. (Fast Red is alcohol soluble; do not dehydrate.)</p>	<p>10 - 20 min</p>

12. HEMATOXYLIN Not provided	a. Counterstain with 2 drops (100µl) or enough volume of hematoxylin to completely cover tissue. Incubate for 10-15 seconds. b. Rinse thoroughly with tap water for 2-3 min. c. Put slides in PBS until show blue color (about ½ - 1 min.) d. Rinse well in distilled water.	
13. Reagent 6 NeoMount Universal (RTU)	a. Apply 2 drops (100µl) or enough volume of Reagent 6 NeoMount Universal to cover tissue when tissue is wet. Rotate the slides to allow NeoMount Universal spread evenly. DO NOT coverslip. b. Place slides horizontally in an oven at 40-50°C for at least 30 minutes or leave at room temperature until slides are thoroughly dried. Hardened NeoMount Universal forms an impervious polymer barrier to organic solvent. Do not use oil directly on the top of dried NeoMount Universal.	30 min. in 40- 50°C oven Or overnight at room temperature

PROTOCOL NOTES:

1. The fixation, tissue slide thickness, antigen retrieval and primary antibody dilution and incubation time effect results significantly. Investigator needs to consider all factors and determine optimal conditions when interpreting the results.
2. NeoMount Universal is water-based mounting medium for immunohistochemistry. It may be used as the permanent mounting media. It does not need coverslip. However, if you need to coverslip, after Simpo-Mount dried, dip the slide in xylene and take out immediately. Apply O-Mount (NeoMount Perm, Catalog No. NB-23-00156) on the tissue and place cover glass on the slide. Store it after dry completely. Step not recommended.

PRECAUTIONS:

DAB may be carcinogenic. Appropriate personal protective equipment should be worn at all times while handling.

FOR RESEARCH USE



