

UltraMarathonRT

20,000 U/mL | Store at -20°C

TEMPLATE SWITCHING KIT

UltraMarathonRT (uMRT) undergoes an efficient and highly selective form of template switching that can be used to introduce a specific primer binding site to the 5' end of RNAs within a library. Using the Template Switching Oligonucleotides (TSO) as its new template, uMRT then adds a section of known sequence to the 3' end of the resulting cDNA. In this way, cDNA products are tagged at each end, enabling efficient PCR amplification and RNA-seq library preparation.

➔ To complete the template switching experiment, additional reagents not contained in this kit are necessary.

Components Not Provided

dNTP mix (10 mM each) • dATP (10 mM) • uMRT dT₁₈ & AmpPCR Primer
RNaseOUT • KAPA HiFi DNA polymerase (Roche, Cat# KK2102)

uMRT dT₁₈: RT primer:

CCCTCTCTCTCTTTCTCTCTCTTTTTTTTTTTTTTTTTT

AmpPCR: PCR primer for cDNA amplification:

CCCTCTCTCTCTTTCTCTCTC



Distributed by:

CliniSciences Group

Template Switching Quick Start Protocol

ANNEALING RT PRIMERS TO RNA TEMPLATES

1. Gently mix the following components in a nuclease-free microcentrifuge tube and collect the liquid with a quick spin.

Components		Volume (3 μ L total)	Final
Primer	uMRT-dT18 (5 μ M)	0.5 μ L	2.5 pmol
Template RNA	Total RNA	variable	1 ng-20 ng
	poly(A)-RNA	variable	0.1 ng-20 ng
dNTP mix (10 mM each)		0.5 μ L	0.5 mM
Nuclease-free Water		added to total 3 μ L	

2. Incubate at 95°C for 30 sec and then snap cool on ice to anneal the primer to the template.

REVERSE TRANSCRIPTION

3. Gently mix the following components in a nuclease-free microcentrifuge tube and collect the liquid with a quick spin.

Components	Volume (7 μ L total)
2x UltraMarathonRT Buffer	5 μ L
UltraMarathonRT Enzyme	0.5 μ L
UltraMarathonRT Boost	0.5 μ L
RNase Out™ (40 U/ μ L) (optional)	0.5 μ L
Nuclease-free Water	0.5 μ L

4. Add RT reaction mix to the annealed RNA and mix gently by tapping the tube (for a total volume of 10 μ L).
5. Incubate the mixture at 42°C for 45 min and then hold the tubes at 4°C in a PCR thermocycler.

Distributed by:

CliniSciences Group



Above: uMRT Template Switching Kit Components

TEMPLATE SWITCHING

- Gently mix the following components in a nuclease-free microcentrifuge tube and collect the liquid with a quick spin.

Note: Use a heated lid (55°C) to avoid water evaporation during the reaction.

Components	Volume (10 μ L total)
5x Template Switching Buffer	4 μ L
UltraMarathonRT Enzyme	2 μ L
uMRT-TSO (10 μ M)	2 μ L
dATP (10 mM)	2 μ L

- Add the template switching reaction mix to the RT reaction tube from step 5 (for a total volume of 20 μ L) and gently mix by tapping the tube.
- Incubate at 42°C for 45 min to catalyze the template switching reaction.
- Inactivate the enzyme by heating at 95°C for 1 min. then hold the tube at 4°C in a PCR thermocycler.

Note: Use a heated lid (55°C) to avoid water evaporation during the reaction.

Note: The resulting cDNA contains adapters at both ends and can be stored at -20°C until future use.

Distributed by:

CliniSciences Group

PRE-AMPLIFY CDNA USING KAPA HiFi

10. Gently mix the following components in a nuclease-free microcentrifuge tube and collect the liquid with a quick spin.

Components	Volume (100 µL total)
Unpurified cDNA product (from Step 9)	20 µL
PCR primer AmpPCR (20 uM)	10 µL
5x GC enhanced buffer	20 µL
dNTP mix (10 mM each)	3 µL
KAPA HiFi DNA Polymerase (1 U/µL)	2 µL
Nuclease-free Water	45 µL

11. Run the below PCR program in a thermocycler.

KAPA HiFi PCR Conditions				
Cycle	Denature	Anneal	Extend	Hold
1	98°C, 2 min	–	–	–
11-15*	98°C, 15 s	62°C, 30 s	72°C, 6 min	–
1	–	–	72°C, 5 min	4°C
∞	–	–	–	4°C

Note: 12 cycles are suggested for 10 ng total RNA or 1 ng mRNA input. 15 cycles are suggested for 1 ng total RNA or 0.1 ng of mRNA.

Note: The resulting pre-amplified cDNA contains adapters at both ends and can be stored at -20°C.

This product and its labels and packaging may be covered by one or more patents, trademarks or copyrights owned or licensed by RNAConnect, Inc. The use of this product may require additional third-party intellectual property rights for certain applications. This product is for research use only and is not intended for use in humans.

CliniSciences Group

Austria

Company: ClniSciences GmbH
Address: Sternwartestrasse 76, A-1180
Wien - Austria
Telephone: +43 720 115 580
Fax: +43 720 115 577
Email: oesterreich@clinisciences.com
Web: <https://www.clinisciences.com>



Belgium

Company: ClniSciences S.R.L.
Address: Avenue Stalingrad 52, 1000
Brussels - Belgium
Telephone: +32 2 31 50 800
Fax: +32 2 31 50 801
Email: belgium@clinisciences.com
Web: <https://www.clinisciences.com>



Denmark

Company: ClniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: denmark@clinisciences.com
Web: <https://www.clinisciences.com>



Finland

Company: ClniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: suomi@clinisciences.com
Web: <https://www.clinisciences.com>



France

Company: ClniSciences S.A.S
Address: 74 Rue des Salses, 92000
Nanterre - France
Telephone: +33 9 77 40 09 09
Fax: +33 9 77 40 10 11
Email: info@clinisciences.com
Web: <https://www.clinisciences.com>



Germany

Company: Biotrend Chemikalien GmbH
Address: Wilhelm-Mauser-Str. 41-43,
50827 Köln - Germany
Telephone: +49 221 9498 320
Fax: +49 221 9498 325
Email: info@biotrend.com
Web: <https://www.biotrend.com>



Iceland

Company: ClniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: island@clinisciences.com
Web: <https://www.clinisciences.com>



Ireland

Company: ClniSciences Limited
Address: Ground Floor, 71 lower Baggot street
Dublin D02 P593 - Ireland
Telephone: +353 1 6971 148
Fax: +353 1 6971 147
Email: ireland@clinisciences.com
Web: <https://www.clinisciences.com>



Italy

Company: ClniSciences S.r.l.
Address: Via Maremmana inferiore 378
Roma D0012 Galdonia Montacello - Italy
Telephone: +39 06 94 80 56 71
Fax: +39 06 94 80 00 21
Email: italia@clinisciences.com
Web: <https://www.clinisciences.com>



Netherlands

Company: ClniSciences B.V.
Address: Kralenhoftstraat 137A,
1018RG Amsterdam, Netherlands
Telephone: +31 85 2082 351
Fax: +31 85 2082 353
Email: nederland@clinisciences.com
Web: <https://www.clinisciences.com>



Norway

Company: ClniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: norge@clinisciences.com
Web: <https://www.clinisciences.com>



Poland

Company: ClniSciences sp. z o.o.
Address: ul. Rolmistwa Witolda Pileckiego 67
lok. 200 - 02-781 Warszawa - Poland
Telephone: +48 22 307 0535
Fax: +48 22 307 0532
Email: polska@clinisciences.com
Web: <https://www.clinisciences.com>



Portugal

Company: Quimigen Unipessoal LDA
Address: Rua Almada Negreiros, Lote 5, Loja 14,
2615-275 Alverca Do Ribatejo - Portugal
Telephone: +351 30 8808 050
Fax: +351 30 8808 052
Email: info@quimigen.pt
Web: <https://www.quimigen.pt>



Spain

Company: ClniSciences Lab Solutions
Address: C/ Hermanos del Moral 13
(Bajo E), 28019, Madrid - Spain
Telephone: +34 91 269 40 65
Fax: +34 91 269 40 74
Email: espana@clinisciences.com
Web: <https://www.clinisciences.com>



Sweden

Company: ClniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: sverige@clinisciences.com
Web: <https://www.clinisciences.com>



Switzerland

Company: ClniSciences Limited
Address: Marktgasse 18 8302 Kloten -
Switzerland
Telephone: +41 (0)44 805 76 81
Fax: +41 (0)44 805 76 75
Email: switzerland@clinisciences.com
Web: <https://www.clinisciences.com>



UK

Company: ClniSciences Limited
Address: 11 Progress Business center, White
Parkway, SL1 6DQ Slough - United Kingdom
Telephone: +44 (0)1753 866 511
or +44 (0) 330 684 0982
Fax: +44 (0)1753 208 899
Email: uk@clinisciences.com
Web: <https://www.clinisciences.com>



USA

Company: Biotrend Chemicals LLC
Address: c/o Carr Riggs Ingram,
500 Grand Boulevard, Suite 210 Miramar
Beach, FL 32550- USA
Telephone: +1 850 650 7790
Fax: +1 850 650 4383
Email: info@biotrend-usa.com
Web: <https://www.biotrend-usa.com>

