

## EZgrip™ HDPE Carboy

- **VersaCap®** technology provides the flexibility of interchangeable connectors
- **EZgrip™** body features advanced ergonomic design
- **EZfil™** wide-mouth design for easy filling and cleaning
- **VersaBarb®** technology provides interchangeable spigot options

### Made in the USA



### Caution

- Do not use carboys under pressure or full vacuum
- Do not use carboys at temperatures greater than the rated maximum
- Do not store oxidizing agents in carboys
- Do not flame carboy
- Do not use carboy spigot as a handle
- Do not use for long-term storage of strong organic compounds

Foxx Life Sciences' **EZgrip™** Carboy is an advanced fluid handling system developed with the end-user in mind. The ergonomically designed **EZgrip™** body, coupled with our **VersaCap®** technology provides the user with an easy-to-use and easy-to-clean carboy.

The **EZgrip™** HDPE Carboy is made from a translucent USP Class VI, FDA Grade resin, with excellent chemical resistance, making it ideal for most buffers and reagents as well as powders and solids.


The **EZgrip™** Carboy's body includes design features that allow the user to easily hold and manipulate during lifting, filling, and pouring, while **VersaCap®** technology provides interchangeable, leak-proof top connections. **EZgrip™** Carboy features include large top handles, molded hand grips in the base, easy-to-read metric graduation marks certified to  $\pm 5\%$  accuracy, and the carboy material molded into the side for easy identification.

The **EZgrip™** HDPE Carboy offers an **EZfil™** 120mm neck for easier filling and cleaning. This is available on 20L, 40L, and 75L models.

ISO 13485

- 1 Novel shape maximizes bench and shelf space
- 2 Bottom handles provide a secure grip while lifting and transporting
- 3 Large top handles for a full hand grip
- 4 Molded carboy identification marks
- 5 Large, clear graduation marks with  $\pm 5\%$  accuracy
- 6 Autoclavable spigot and **VersaBarb®** technology
- 7 **VersaCap®** technology with versatile adapters and replaceable gasket
- 8 **EZgrip™** body design allows for maximum solution recovery
- 9 **EZfil™** wide-mouth shape for easy filling and cleaning





Product Number	Description	Capacity (L)	Material	Spigot	Cap Size (mm)	Qty/Pack
152-0111-OEM	Carboy, HDPE, 2.5L, <b>VersaCap®</b> 83B	2.5	HDPE	No	83	1
152-1111-OEM	Carboy, HDPE, 5L, <b>VersaCap®</b> 83B	5	HDPE	No	83	1
152-1211-OEM	Carboy, HDPE, 5L, Spigot, <b>VersaCap®</b> 83B	5	HDPE	Yes	83	1
152-2111-OEM	Carboy, HDPE, 10L, <b>VersaCap®</b> 83B	10	HDPE	No	83	1
152-2211-OEM	Carboy, HDPE, 10L, Spigot, <b>VersaCap®</b> 83B	10	HDPE	Yes	83	1
152-3111-OEM	Carboy, HDPE, 20L, <b>VersaCap®</b> 83B	20	HDPE	No	83	1
152-3211-OEM	Carboy, HDPE, 20L, Spigot, <b>VersaCap®</b> 83B	20	HDPE	Yes	83	1
152-3121-OEM	Carboy, HDPE, 20L, <b>VersaCap®</b> 120mm	20	HDPE	No	120	1
152-3221-OEM	Carboy, HDPE, 20L, Spigot <b>VersaCap®</b> 120mm	20	HDPE	Yes	120	1
152-5121-OEM	Carboy, HDPE, 40L, <b>VersaCap®</b> 120mm	40	HDPE	No	120	1
152-5221-OEM	Carboy, HDPE, 40L, Spigot <b>VersaCap®</b> 120mm	40	HDPE	Yes	120	1
152-7121-OEM	Carboy, HDPE, 75L, No Spigot <b>VersaCap®</b> 120mm	75	HDPE	No	120	1



# CliniSciences Group

## Austria

Company: CliniSciences GmbH  
Address: Sternwartestrasse 76, A-1180  
Wien - Austria  
Telephone: +43 720 115 580  
Fax: +43 720 115 577  
Email: [oesterreich@clinisciences.com](mailto:oesterreich@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Belgium

Company: CliniSciences S.R.L  
Address: Avenue Stalingrad 52, 1000  
Brussels - Belgium  
Telephone: +32 2 31 50 800  
Fax: +32 2 31 50 801  
Email: [belgium@clinisciences.com](mailto:belgium@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Denmark

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [denmark@clinisciences.com](mailto:denmark@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Finland

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [suomi@clinisciences.com](mailto:suomi@clinisciences.com)  
Web: <https://www.clinisciences.com>



## France

Company: CliniSciences S.A.S  
Address: 74 Rue des Suisses, 92000  
Nanterre- France  
Telephone: +33 9 77 40 09 09  
Fax: +33 9 77 40 10 11  
Email: [info@clinisciences.com](mailto:info@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Germany

Company: Biotrend Chemikalien GmbH  
Address: Wilhelm-Mauser-Str. 41-43,  
50827 Köln - Germany  
Telephone: +49 221 9498 320  
Fax: +49 221 9498 325  
Email: [info@biotrend.com](mailto:info@biotrend.com)  
Web: <https://www.biotrend.com>



## Iceland

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [island@clinisciences.com](mailto:island@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Ireland

Company: CliniSciences Limited  
Address: Ground Floor, 71 lower Baggot street  
Dublin D02 P593 - Ireland  
Telephone: +353 1 6971 146  
Fax: +353 1 6971 147  
Email: [ireland@clinisciences.com](mailto:ireland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Italy

Company: CliniSciences S.r.l  
Address: Via Maremmana inferiore 378  
Roma 00012 Guidonia Montecelio - Italy  
Telephone: +39 06 94 80 56 71  
Fax: +39 06 94 80 00 21  
Email: [italia@clinisciences.com](mailto:italia@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Netherlands

Company: CliniSciences B.V.  
Address: Kraijenhoffstraat 137A,  
1018RG Amsterdam, - Netherlands  
Telephone: +31 85 2082 351  
Fax: +31 85 2082 353  
Email: [nederland@clinisciences.com](mailto:nederland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Norway

Company: CliniSciences AS  
Address: c/o MerVerdi Munkerudtunet 10  
1164 Oslo - Norway  
Telephone: +47 21 988 882  
Email: [norge@clinisciences.com](mailto:norge@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Poland

Company: CliniSciences sp.Z.o.o.  
Address: ul. Rotmistrza Witolda Pileckiego 67  
lok. 200 - 02-781 Warszawa -Poland  
Telephone: +48 22 307 0535  
Fax: +48 22 307 0532  
Email: [polska@clinisciences.com](mailto:polska@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Portugal

Company: Quimigen Unipessoal LDA  
Address: Rua Almada Negreiros, Lote 5, Loja 14,  
2615-275 Alverca Do Ribatejo - Portugal  
Telephone: +351 30 8808 050  
Fax: +351 30 8808 052  
Email: [info@quimigen.com](mailto:info@quimigen.com)  
Web: <https://www.quimigen.pt>



## Spain

Company: CliniSciences Lab Solutions  
Address: C/ Hermanos del Moral 13  
(Bajo E), 28019, Madrid - Spain  
Telephone: +34 916 750 700  
Fax: +34 91 269 40 74  
Email: [espana@clinisciences.com](mailto:espana@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Sweden

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [sverige@clinisciences.com](mailto:sverige@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Switzerland

Company: CliniSciences AG  
Address: Fracht Ost Flughafen Kloten  
CH-8058 Zürich - Switzerland  
Telephone: +41 (044) 805 76 81  
Fax: +41 (044) 805 76 75  
Email: [switzerland@clinisciences.com](mailto:switzerland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## UK

Company: CliniSciences Limited  
Address: 11 Progress Business center, Whittle  
Parkway, SL1 6DQ Slough- United Kingdom  
Telephone: +44 (0)1753 866 511  
or +44 (0) 330 684 0982  
Fax: +44 (0)1753 208 899  
Email: [uk@clinisciences.com](mailto:uk@clinisciences.com)  
Web: <https://www.clinisciences.com>



## USA

Company: Biotrend Chemicals LLC  
Address: c/o Carr Riggs Ingram,  
500 Grand Boulevard, Suite 210 Miramar  
Beach, FL 32550- USA  
Telephone: +1 850 650 7790  
Fax: +1 850 650 4383  
Email: [info@biotrend-usa.com](mailto:info@biotrend-usa.com)  
Web: <https://www.biotrend-usa.com>

