A15523

Leader in Biomolecular Solutions for Life Science

VMP1 Rabbit pAb

Catalog No.: A15523



Basic Information

Observed MW 42kDa

Calculated MW 46kDa

Category Polyclonal Antibody

Applications WB,IF/ICC,ELISA

Cross-Reactivity Human,Mouse,Rat

Background

This gene encodes a transmembrane protein that plays a key regulatory role in the process of autophagy. The ectopic overexpression of the encoded protein in cultured cells triggers autophagy even under nutrient-rich conditions. This gene is overexpressed in pancreatitis affected acinar cells where the encoded protein mediates sequestration and degradation of potentially deleterious activated zymogen granules in a process termed, zymophagy.

Recommended Dilutions

Immunogen Information

WB	1:200 - 1:2000	(
IF/ICC	1:50 - 1:200	

Gene ID	Swiss Prot
81671	Q96GC9

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 300 to the C-terminus of human VMP1 (NP_112200.2).

Synonyms

EPG3; TANGO5; TMEM49; VMP1

€

Product Information


www.abclonal.com

Source Rabbit **Isotype** IgG Purification Affinity purification


Storage

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.


Validation Data



Western blot analysis of lysates from mouse kidney, using VMP1 Rabbit pAb (A15523) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 5s.



Immunofluorescence analysis of NIH-3T3 cells using VMP1 Rabbit pAb (A15523) at dilution of 1:100 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) (AS007) at 1:500 dilution. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of U-2 OS cells using VMP1 Rabbit pAb (A15523) at dilution of 1:100 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) (AS007) at 1:500 dilution. Blue: DAPI for nuclear staining.