

Rabbit IFNG (Interferon gamma) Sandwich ELISA Kit (NB-22-45289-1) NB-22-45289-1


GENERAL INFORMATION

	Kits This IFNG Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for
Description	the measurement of samples in rabbit cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Application Reactivity	Elisa Rabbit

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	9.38pg/mL
Detection	15.63~1000pg/mL
Limit	
Assay Time	3.5h
Storage	If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be
Instruction	used within 1 month, store the components separately, according to the component
	table in the manual
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.