

Product datasheet MON3242



Mouse anti-CD61, clone 2f2 (Monoclonal)

Clone no. 2f2

MONOSAN

Product name	Mouse anti-CD61, clone 2f2 (Monoclonal)
Host	Mouse
Applications	IHC-P (1:25-1:100)
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1-k
Clonality	Monoclonal
Clone number	2f2
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-CD61, clone 2f2 (Monoclonal)

Clone no. 2f2

MONOSAN

Additional info

CD61 also known as integrin beta chain beta 3 (ITGB3) is an integrin cell-surface protein associated with cellular adhesion and cell-surface mediated signaling. Immunohistochemical staining for CD61 can be useful in evaluating normal and abnormal megakaryocytes, which can aid in the identification of some hematopoietic malignancies. Anti-CD61 reactivity is also seen in platelets, osteoclasts and macrophages.

References

1. Duperray A et al. Blood. 1989 Oct; 74(5):1603-11
2. Goldman BI et al. Modern Pathology 14:589-594 (2001)
3. Thiele J et al. Virchows Archiv B Cell Pathol (1990) 58:295-302
4. Thiele J et al. Eur J Haematol 1990; 44:63-70
5. Campana D et al Leukemia. 1990 Sep; 4(9):620-4

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES