

Mouse anti-Endothelium, clone PAL-E (Monoclonal)

Clone no. PAL-E

MONOSAN

---

Product name	Mouse anti-Endothelium, clone PAL-E (Monoclonal)
Host	Mouse
Applications	IHC-Fr (1:20)
Species reactivity	human, rabbit, goat, dog, pig
Conjugate	-
Immunogen	human melanoma tissue
Isotype	IgG2a
Clonality	Monoclonal
Clone number	PAL-E
Size	1 ml
Concentration	15 ug/ml
Format	-
Storage buffer	Culture medium with 0.01M Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

## Mouse anti-Endothelium, clone PAL-E (Monoclonal)

Clone no. PAL-E

MONOSAN

**Additional info**

The antibody only stains human and various animal blood vessels with exception of arteries. It does not stain rat, mouse and chicken endothelium. The antibody is useful for immunohistochemistry on acetone fixed frozen sections. Not useful on cellsuspensions and Western blotting.

**References**

1. Jones, R.R., et al., J. Clin. Path. 39, 742-749 (1986).
2. Holden, C.A., et al., J. Immun. Meth. 91, 45-52 (1986).
3. Schlingemann, R.O., et al., Laboratory Investigation 52, 71 (1987).
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES