

# PAX6 antibody

---

## Product Description

Goat Polyclonal to PAX6.

---

## Product Information

<b>Code:</b>	NB-22-25551
<b>Clonality</b>	Polyclonal
<b>Host</b>	Goat
<b>Reactivity</b>	Human, Mouse, Rat
<b>Applications</b>	ELISA, WB
<b>Immunogen Region</b>	Internal
<b>Conjugation</b>	Unconjugated
<b>Gene ID</b>	18508 (mouse); 25509 (rat);
<b>Isotype</b>	IgG
<b>Accession Number</b>	NP_000271.1; NP_001595.2
<b>Purification</b>	Affinity purification
<b>Form</b>	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Sequence</b>	REEKLRNQRRQASN
<b>Storage</b>	Store at -20°C. Avoid repeated freeze/thaw cycles.
<b>Alternative names</b>	PAX6; paired box gene 6 (aniridia, keratitis); AN; AN2; D11S812E; MGC17209; MGDA; WAGR; paired box gene 6; paired box homeotic gene 6 (aniridia, keratitis); paired box homeotic gene-6

---