

Cell Culture Bags

Rocker Style Bioreactor

Expression Systems' new single-use cell culture bags feature a polyethylene film and are compatible with a variety of rocking platforms. Ideal for insect cell culture, these bags are leak-tested and irradiated.

- Manufactured in the USA and compatible with GE WAVE and other popular rocking platforms
- Ideal for use with the Baculovirus and Host Cell Expression System
- Standard ports include:
 - Male and female quick connects
 - Female Luer
 - Needleless sample port
- 1/8" ID x 1/4" OD and 1/4" ID x 3/8" OD tubing for welding compatibility
- Ideal for scale-up manufacturing
- Available in 10L to 25L maximum working volume sizes
- Animal component free product contact surfaces



Specifications

Film Composition	Low density polyethylene
Platform compatibility	Fits multiple rocking platforms
Shipping condition	Ambient
Storage condition	Room temperature, protect from light

For sample bags and updates— email contact info to sales@expressionsystems.com



**Expression
Systems**

Insect and Mammalian
Cell Culture Media, Systems
and Services


2537 2nd Street
Davis CA 95618

t 530 747-2035 | 877 877-7421
f 530 747-2034


Cell Culture Bags

Rocker Style Bioreactor

25L CCB
28" x 29"



10L SBS CCB
28" x 14"



P1: 1/4" ID x 3/8" OD x 39"
WITH FEMALE MPC

P2: 1/8" ID x 1/4" OD x 39"
WITH FEMALE LUER

P3: NEEDLELESS SAMPLE PORT

P4: 1/4" ID x 3/8" OD x 39"
WITH MALE MPC

P5: FILTERED AIR OUTLET

P6: FILTERED AIR INLET

Products

Expression Systems offers a variety of media formulations for insect and mammalian cell culture. Expression Systems provides tools to support the baculovirus and host cell expression system including baculovirus vectors and kits to determine baculovirus production.

	Cell Culture Bags	Dimensions		Cell Culture Media and Reagents
90-005	10L Side by Side Cell Culture Bag use 2 bags on a 25L bag platform	28" x 14"	96-001	ESF 921™ Insect Cell Culture Medium, Protein Free
90-007	10L Working Volume Cell Culture Bag	21.85" x 25.9"	99-300	ESF AF™ Insect Cell Culture Medium, Animal Free
90-006	25L Working Volume Cell Culture Bag	28" x 29"	95-006	Production Boost Additive

About Expression Systems

Headquartered in Davis, California, Expression Systems specializes in the products, services and support necessary for high quality baculovirus protein production. Our unique focus in BEVS applications helps researchers and clinicians around the world work at the highest level of consistency. At any step in the process, from gene synthesis to protein purification, Expression Systems helps achieve best-in-class results.

