

**Product Description:** CCKB antibody, human, C-term.  
CCK2  
Cholecystinin type B receptor  
Cholecystinin type 1 receptor  
G Protein-coupled receptor  
Polyclonal antibody raised against the  
C-terminus of the human CCKB receptor

**Product-No.:** S11-9032-S-R

**Immunogen:** residues 428-447 coupled to KLH (NP\_795344)  
KLH- psiaslsrslsyttistlpgp

**Host:** Rabbit

**Form:** serum with 0.02% NaN<sub>3</sub> as preservative

**Quantity:** 50 µl

**Applications:** IC, W  
W = Western blot (1:5000)  
IC = Immunocytochemistry (1:5000)

**Species cross-reactivity:** human

**Sensitivity:** 50 ng

**Storage:** 4 °C short term  
-20 °C long term

