Product datasheet MON3381



Mouse anti-Chromogranin A, clone LK2H10 (Monoclonal)

Clone no. LK2H10 MONOSAN

Product name Mouse anti-Chromogranin A, clone LK2H10 (Monoclonal)

Host Mouse

Applications IHC-P (1:100-1:500)

Species reactivity human

Conjugate -

Immunogen Unknown or proprietery to MONOSAN and/or its suppliers

lsotype lgG1-k

Clonality Monoclonal

Clone number LK2H10

Size 1 ml

Concentration n/a

Format -

Storage buffer Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON3381



Mouse anti-Chromogranin A, clone LK2H10 (Monoclonal)

Clone no. LK2H10 MONOSAN

Additional info

Immunohistochemical methods have localized chromogranin in a wide variety of endocrine tissues including the pituitary, pancreas, thyroid, and parathyroid. Neuroendocrine cells exhibit a fine granular immunoreactivity to chromogranin. It is generally accepted that the co-expression of certain keratins and chromogranin mean neuroendocrine lineage. The presence of strong chromogranin staining and absence of keratin staining should raise the possibility of paraganglioma. The co-expression of chromogranin and NSE is typical of neuroendocrine neoplasms.

References

- 1. Wilson, BS, et al., Am | Pathol; 115:458-468 (1984)
- 2 Lyda MH, Weiss LM. Hum Pathol. 31(8):980-7 (2000)
- 3. ontochristopoulous GJ et al. Dermatology.; 201(2):123-6 (2000)
- 4. Qvigstad G et al. Histochem J.; 32(9):551-6 (2000)
- 5. Bovin AM, et al. APMIS.; 110(9):658-64 (2002)

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES