

Product Description: NT1 antibody, human, C-term.
NTR1
Neurotensin 1 receptor
G Protein-coupled receptor
Polyclonal antibody raised against the
C-terminus of the human NT1 receptor

Product-No.: S12-9036-S-R

Immunogen: residues 398-418 coupled to KLH (NP_002522)
KLH- kadsvssnhtlssnatretly

Host: Rabbit

Form: serum with 0.02% NaN₃ as preservative

Quantity: 50 µl

Applications: IC, W
W = Western blot (1:5000)
IC = Immunocytochemistry (1:5000)

Species cross-reactivity: human

Sensitivity: 50 ng

Storage: 4 °C short term
-20 °C long term

