

---

## Antibody Diluent Solution

Cat # NB-23-00171-1

Size: 500ml (ready-to-use)

Cat # NB-23-00171-2

Size: 300ml as 6 x 50ml (ready-to-use)

### Intended Use

---

Antibody Diluent Solution is used for dilute primary and secondary antibodies. The solution is not suitable for those antibodies conjugated to peroxidase enzyme because  $\text{NaN}_3$  in the diluent will inhibit the activity of peroxidase.

### Components

---

Solution contains 0.05%  $\text{NaN}_3$ .

### Storage

---

Store at 2-8°C.

### Recommended Protocol

---

Refer the dilution factors on manufacture's specification sheet or insert. For example if dilution is 50x, add 1ml of the concentrated antibody into 50ml of Antibody Diluent to obtain working dilution. Neo Biotech's Antibody Diluent Solution may be used for antibodies from other source.

### Precaution

---

Wear gloves and take other necessary laboratory safety procedure.

***For Research Use Only. Not for Diagnostic or Therapeutic Use.***