Product datasheet

MONX10806



Mouse anti-Myosin Heavy Chain (developmental), clone RNMY2/9D2 (monoclonal)

Clone no. RNMy2/9		MONXtra
Product name	Mouse anti-Myosin Heavy Chain (developmental), clone RNMY2/9D (monoclonal)	2
Host	Mouse	
Applications	IHC-fr (1:10-1:40)	
Species reactivity	human, rat	
Conjugate	-	
Immunogen	Native myosin extracted from the hind limb muscle of 7 day old rat	S.
lsotype	lgG1	
Clonality	Monoclonal	
Clone number	RNMy2/9D2	
Size	1 ml	
Concentration	n/a	
Format	-	
Storage buffer	Lyophilized	
Storage until expiry date	2-8°C	

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com

Product datasheet



Mouse anti-Myosin Heavy Chain (developmental), clone RNMY2/9D2 (monoclonal)

Clone no. RNMy2/9D2

MONXtra

Additional info

Myosin is a contractile muscle specific protein composed of two heavy and four light chains. The myosin heavy chain has many isoforms which are specific for different muscles or fiber types, some of which are developmentally regulated. The range of myosin heavy chain antibodies may prove useful for investigating development of intrafusal and extrafusal muscle fibers and the course of muscle fiber regeneration. At the ultrastructural level, antibodies can reveal architectural details of the myofilament as well as the cytoplasmic and membrane sites of new myosin integration. The user is required to reconstitute the contents of the vial with the correct volume of sterile distilled water as indicated on the vial label

D C	
Reter	ences
I CICI	CIICCS

1.

2

3. 4.

5.

Davis CE et al. Neuromuscular Disorders. 1 (6): 411–421 (1991) Ecob-Prince M et al. Journal of Neurological Sciences. 91: 71–78 (1989) --

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES