

# High Speed Mini Centrifuge - N400 Series



## Features :

1. Extremely quiet operation
2. Easy-to-follow LCD display
3. Capacity: 12 x 1.5/2.0ml tubes
4. Separate short-spin key
5. Low sample heating (only 12°C after 30min, at max.speed)
6. Autoclavable rotor (121°C, 20min)
7. RMP/RCF setting as required
8. Includes standard rotor made of anodized aluminum
9. Timer can be set up to 99min
10. The small size of the unit makes it convenient for benchtop use or in a hood or in refrigerated environments

## Specifications :

Model	N400C
Cat#	NB-12-0022C
Power supply	100V or 240V 50/60Hz
Power requirement	100W
Max. rotational speed	15000 rpm/min
Relative centrifugal force (RCF)	12,100 x g
Timer	20 sec-99min
Acceleration time to max. speed	48 s
Braking time from max. speed	52 s
Dimensions in cm (WxDxH)	24x21.2x13.8



## Rotor and accessories :

1. High-strength transparent plastic cover
2. Rotor can be autoclavable (121°C, 20min)
3. Optional adapter accessory
  - 0.2ml PCR tube adapter
  - 0.5ml PCR tube adapter

Cat#	Description
NB-12-0022-01	Adapters for 0.5ml tubes
NB-12-0022-02	Adapters for 0.2ml tubes

