

Mouse anti-Liver Fatty Acid Binding Protein, clone L2B10, Biotin (Monoclonal)

Clone no. L2B10

MONOSAN

---

Product name	Mouse anti-Liver Fatty Acid Binding Protein, clone L2B10, Biotin (Monoclonal)
Host	Mouse
Applications	-
Species reactivity	-
Conjugate	Biotin
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	-
Clonality	Monoclonal
Clone number	L2B10
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	PBS with 0.1% BSA and 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

MON6023B

**MONOSAN**<sup>®</sup>

Mouse anti-Liver Fatty Acid Binding Protein, clone L2B10, Biotin (Monoclonal)

Clone no. L2B10

MONOSAN

---

Additional info

-

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES