Leader in Biomolecular Solutions for Life Science

NMUR1 Rabbit pAb

Catalog No.: A17610



Basic Information

Observed MW

53kDa

Calculated MW

47kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human, Mouse

Background

Enables neuromedin U binding activity and neuromedin U receptor activity. Involved in several processes, including activation of phospholipase C activity; chloride transport; and second-messenger-mediated signaling. Is integral component of membrane.

Recommended Dilutions

WB

1:500 - 1:2000

Immunogen Information

Gene ID 10316

Swiss Prot Q9HB89

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 360-426 of human NMUR1 (NP_006047.3).

Synonyms

FM3; FM-3; GPC-R; GPR66; NMU1R; (FM-3); NMUR1

Contact

•

www.abclonal.com

Product Information


Source Rabbit **Isotype** IgG **Purification**Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.

Validation Data



Western blot analysis of lysates from SGC-7901 cells, using NMUR1 Rabbit pAb (A17610) at 1:1000

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 90s.