

Anti-WBP2NL antibody

Cat# NB-22-47070-50 Size: 50µL

Cat# NB-22-47070-100 Size: 100µL

Cat# NB-22-47070-200 Size: 200µL

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-WBP2NL is suitable for use in Western Blot and Immunofluorescence.
Applications	WB, IF
Host/Source	Rabbit
Reactivity	Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	164684
Gene Symbol	WBP2NL
Uniprot ID	WBP2L_HUMAN
Immunogen	Recombinant protein of human WBP2NL.
Immunogen Region	
Specificity	
Immunogen sequence	



Western blot analysis of extracts of Rat spleen, using WBP2NL rabbit polyclonal antibody (NB-22-47070) at 1:1000 dilution. Secondary antibody: HRP Goat Antirabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.