

Monoclonal anti-EMI-2 (Endogenous meiotic inhibitor 2/F-box only protein 43)

NB-01-0026

Volume: 100 µl

Clones: 2EMI-2A10

Immunogen: synthetic peptide: CRRGSAKARGSKDVL corresponding to aa 613 to 627 of mouse EMI-2.

Murine Isotype: IgA, k

Specificity: mouse EMI-2.

Applications:

- ELISA
- Western Blot
- Immunocytochemistry

Recommended dilutions: 1/500-1/5000

Storage:

- Store at -20°C (up to 3 years)
- Avoid repeat freezing and thawing cycles. Does not contain any preservative