

Mouse anti-CD90, clone F15-42-1, R-PE (Monoclonal)

Clone no. F15-42-1

MONOSAN

Product name	Mouse anti-CD90, clone F15-42-1, R-PE (Monoclonal)
Host	Mouse
Applications	FC
Species reactivity	human,cynomolgus monkey
Conjugate	RPE
Immunogen	Purified human brain Thy-1.
Isotype	IgG1
Clonality	Monoclonal
Clone number	F15-42-1
Size	100 TESTS
Concentration	n/a
Format	-
Storage buffer	PBS with 0.09% sodium azide and 1% BSA, 5% sucrose
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-CD90, clone F15-42-1, R-PE (Monoclonal)

Clone no. F15-42-1

MONOSAN

Additional info

Mouse anti Human CD90 antibody, clone F15-42-1 recognizes the human CD90 cell surface antigen, a 126;25 kDa glycoprotein homologous to rat Thy1. The antigen is expressed by a subset of CD34+ve cells in the bone marrow and by prothymocytes within the thymus. CD90 is also expressed extensively within the brain. Mouse anti Human CD90 antibody, clone F15-42-1 is routinely tested in flow cytometry on the MOLT4 cell line.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES