Product datasheet MON4006



Mouse anti-Collagen IV, clone MC4-HA (Monoclonal)

Clone no. MC4-HA MONOSAN

Product name Mouse anti-Collagen IV, clone MC4-HA (Monoclonal)

Host Mouse

Applications IHC (1:20), ELISA (1:1000)

Species reactivity human

Conjugate -

Immunogen Collagen IV

Isotype IgM

Clonality Monoclonal

Clone number MC4-HA

Size 1 ml

Concentration 0.1 mg/ml

Format -

Storage buffer Purified and lyophilized

Storage until expiry date -20°C

www.monosan.com

FOR-044 22-12-2021

Revision date 09-03-2021

Product datasheet

MON4006

MONOSAN

Mouse anti-Collagen IV, clone MC4-HA (Monoclonal)

Clone no. MC4-HA MONOSAN

Additional info

Purified and lyophilized from 0.01M sodium phosphate, 0.01M sodium borate, 0.15M sodium chloride, 1% dextran, 1%mannitol, pH 8. For Reconstitution use 1 ml sterile pure water. Binds Human Collagen type IV in native and heat-denatured state. No cross-reactivity to human collagen I, II, III and V.

References

1. -

2 -

3. -

4.

5

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

 www.monosan.com
 Revision date

 FOR-044 22-12-2021
 09-03-2021