Product datasheet PS047



## Polyclonal Rabbit anti-Collagen I affinity purified

Clone no. - MONOSAN

Product name Polyclonal Rabbit anti-Collagen I affinity purified

**Host** Rabbit

Applications IHC-fr (1:20), ELISA (1:3000)

Species reactivity human

Conjugate -

Immunogen Collagen I

Isotype -

**Clonality** Polyclonal

Clone number -

Size 1 ml

**Concentration** 0.1 mg/ml

Format -

Storage buffer Purified and lyophilized

Storage until expiry date -20°C

www.monosan.com

FOR-044 22-12-2021

Revision date 09-03-2021

Product datasheet PS047

**MONOSAN** 

Polyclonal Rabbit anti-Collagen I affinity purified

Clone no. - MONOSAN

## Additional info

Purified and lyophilized from 0.01M sodium phosphate, 0.01M sodium borate, 0.15M sodium chloride, 1% dextran, 1%mannitol, pH 8. For Reconstitution use 1 ml sterile pure water. Binds human Collagen type I. No cross-reactivity to human plasma proteins and human collagen types II, III, IV and V

References

1. -

2 -

3.

4.

5

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

 www.monosan.com
 Revision date

 FOR-044 22-12-2021
 09-03-2021