

# NeoBiotech

Product Description: BRS-3 antibody, human, C-term.  
Bombesin-like receptor 3  
Bombesin receptor  
G Protein-coupled receptor  
Polyclonal antibody raised against the  
C-terminus of the human BRS-3 receptor


Product-No.: NB-49-010-50ul

Immunogen: residues 342-361 coupled to KLH (NP\_001718)  
KLH-seisvtsftgcsvkqaedrf

Host: Rabbit

Form: serum with 0.02% NaN<sub>3</sub> as preservative

Quantity: 50 µl



BRS-3 immunoreactivity in human pancreatic carcinoma.

Applications: ICH (p), IC, W  
W = Western blot (1:5000)  
IC = Immunocytochemistry (1:5000)  
ICH (p) = Immunohistochemistry of formalin-fixed, paraffin-embedded tissues (1:2500)

Species cross-reactivity: human

Sensitivity: 50 ng

Storage: 4° C short term  
-20° C long term

**Neo-Biotech**

**74 rue des Suisses – 92000 Nanterre**