# **NeoReal 1000**

## NB-12-3027-1

### NeoReal 1000 - 2 channels RealTime PCR System

- 1. The new generation of Peltier-based technology, fast ramping rate of up to 7°C/sec.
- 2. Maintenance free, long life LED illuminator and high sensitivity photoelectric detector.
- 3. Synchronous two color detection system, reducing errors caused by delay, shortening scan time.
- 4. Calibration with ROX dye is not required, guarantees extremely uniform fluorescence wells.
- 5. The software could be upgraded free for life.



### **TECHNICAL SPECIFICATIONS**

### **PERFORMANCE**

Sample Block Capacity 48 wells \*0.1ml

Reaction Volume 10-50µl (recommend 20µl)

Tubes Option Low-profile, white or clear PCR tubes or strips of tubes with optical flat cap

Healing & Cooling Technology New generation Peltier technology allow 1.000.000 cycles

Control Methods Operate via PC control or limited function on instrument

Language English

Communication Ports USE 2.0 & LAN

Display 7"Color TFT LCD and Touch Screen

Certified Non CE-IVD

#### **TEMPERATURE**

Block Temp.Range 0°C-105°C

Max.Heating Rate 7°C/sec

Max.Cooling Rate 5°C/sec

Temp.Uniformity ±0.2°C at 90°C



# **NeoReal 1000**

# NB-12-3027-1

Temp.Uniformity ±0.2°C at 90°C

Temp.Accuracy ±0.1°C (10 seconds after reach 90°C)

Display Resolution 0.1°C

Heat Lid Temp.Range 30°C-112°C

Temp.Control Mode Block & Sim-tube mode

Gradient Range 30°C-105°C

Temp.Differential Range 1°C-24°C

#### **OPTICAL MODULE**

Excitation Long life LED

Detection CCDs

Dynamic Range 1-10<sup>10</sup>

Detection Sensitivity ≥1 copy

Calibrated Dyes at Installation F1 : FAM, SYBR Green

F2: VIC, HEX, TET, JOE, CY3, NED

Fluorescence Excitation Range 300-800nm

Fluorescence Detection Range 500-800nm

Date Export Formats TXT, PDF, WORD

#### **OTHER FEATURES**

AC Power Supply 100V-240V, 50-60Hz

Consumption 400W

Net Weight 8.2kg

Dimensions (L×W×H) 320×205×380mm

## **COMPUTER** system requirements

Operating Systems Windows7, Windows10, WindowsXP

CPU Speed 2.0GHz

Display Resolution 1366×768

Memory 4.00 GB RAM

Available Hard Disk >20 GB

