



Mouse OxLDL (Oxidized Low Density Lipoprotein) Sandwich ELISA Kit

Cat # NB-22-66429 Size : 48T; 96T

GENERAL INFORMATION


Product Type	Kits
Short Description	This OxLDL Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	Polyclonal
Sensitivity	2.18pg/mL
Detection limit	6.25-400pg/mL
Assay Time	4.5hrs
Storage Instruction	Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma or other biological fluids.
Gene	
Symbol	
Uniprot ID	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.