

# <u>NeoStir OH100 -</u> <u>OverHead Stirrer</u>

NB-12-4033B


#### **NeoStir OH Elite**

Elite Series Of Overhead Stirrers

While the market is saturated with overhead stirrers, most of them fall short when it comes matching customer's expectations in terms of performance and durability. As soon as these devices encounter high-viscosity or large quantities of solutions, they falter in terms of performance which restricts the users in terms of the work that can be done with these devices. Even maintenance is also a frequent issue with such devices as they frequently break down or have to undergo frequent maintenance which puts a big question on their durability.

Neo Biotechis markedly different when it comes to our approach to designing and developing overhead stirrers. Guided by a dedicated R&D department engaged in constant research about customer needs and user-experience, we were able to design and develop our excellent range of NeoStir Elite stirrers. Featuring intelligent design and high-performance components, these devices come with unique features such as built-in overload, overcurrent and over-temperature safety, Power Intermittent mode and counter speed control to meet

even the most demanding stirring requirements.



## **Product Features**



Brushless DC motor delivers ultra quiet performance and a maintenance-free long life.



Microprocessor controlled for efficient performance even under high load.



Variable stirring speed settings upto 2000 RPM for powerful performance



Variable timer setting ranging from 1 minute to 99 hours 59 mins as well as an Infinite mode



Powerful motor delivers torque upto 100 Ncm to stir solutions with viscosity up to 70,000 mPas & 100 L volume Volume Range Upto 100L



Viscosity Upto 70,000 mpas



Smart Control Stirring feature to ensure constant stirring change even in viscosity



CW & CCW motion for efficient & thorough mixing.



## **Product Features**



Large and clear display enables easy viewing of all critical parameters



Built-in overload, overcurrent and overtemperature safety feature to ensure complete safety



In case of power failure, the Power Intermittent mode starts the device automatically as soon as power is restored and runs for the remaining left time

Safe Low Voltage Device



IP42 class compliant rating for unmatched durability



Small footprint saves valuable bench space







30,000 mpas

#### **NeoStir OH 75E**



#### NeoStir OH 100E



(100Ncm, 100L, 1300RPM) 70,000 mpas

Wide-variety of impellers available to suit all stirring needs


NeoStir OH 50 Elite		
Motor Type	Brushless DC	
Capacity (H20)	50 Liters	
Speed Range	60 - 2000 RPM	
Torque	50 Ncm	
Viscocity	30000 mPas	
Timer	1 min to 99 hr 59 mins & infinity mode	
Direction of Rotation	Reversible - CW/CCW	
Protection Class	IP 42	
Dimension (Head)	102x 249 x 242mm	
Weight	5.1 kgs	

NeoStir OH 75 Elite		
Motor Type	Brushless DC	
Capacity (H20)	75 Liters	
Speed Range	60 - 2000 RPM	
Torque	75 Ncm	
Viscocity	50000 mPas	
Timer	1 min to 99 hr 59 mins & infinity mode	
Direction of Rotation	Reversible - CW/CCW	
Protection Class	IP 42	
Dimension (Head)	102x 249 x 242mm	
Weight	5.1 kgs	

NeoStir OH 100 Elite		
Motor Type	Brushless DC	
Capacity (H20)	100 Liters	
Speed Range	60 - 1300 RPM	
Torque	100 Ncm	
Viscocity	70000 mPas	
Timer	1 min to 99 hr 59 mins & infinity mode	
Direction of Rotation	Reversible - CW/CCW	
Protection Class	IP 42	
Dimension (Head)	102x 249 x 242mm	
Weight	5.1 kgs	

## **Ordering Information**

NeoStir OH Elite		
Product Name	Cat #	
NeoStir OH 50 Elite	OHS2 - 1101-EU/UK/AU/US	
NeoStir OH 75 Elite	OHS2 - 1201-EU/UK/AU/US	
NeoStir OH 100 Elite	OHS2 - 1301-EU/UK/AU/US	



Accessories		
Cat #	Product Name	Details
OHS-C1	Key Chuck	
OHS-C2	Keyless Chuck	
OHS-C3	Collet Chuck	
OHS-CC	Chuck Cover	