

Address Avda, Universided de Coimba, sin Cancer Research Center (CLC.) Campus Miguel de Unamuno

37007 Salamanca (Spain)

Tel. / Fax: [+34) 923 294 827 E-mail: info@immunostep.com

www.immunostep.com

Anti-Human CD3/CD19

(33-2A3/A3-B1)

Fluorochrome	Reference	Size	
FITC/PE	3F119PE1-50T	50 test	

PRODUCT DESCRIPTION

Clones: 33-2A3, A3-B1

Isotype: Mouse IgG2a, Mouse IgG2a **Tested application**: flow cytometry

Immunogen: The anti-CD3 monoclonal antibody

derives from human leukocytes.

The Mouse anti-human CD19 monoclonal antibody

derives from human tonsil. **Species reactivity:** Human

Storage instruction: store in the dark at 2-8 °C

Storage buffer: aqueous buffered solution containing protein stabilizer and 0.09% sodium azide (NaN $_3$).

Recommended usage: Immunostep's CD3/CD19, is recommended for use in flow cytometry for simultaneous detection and enumeration T and B Lymphocytes. This reagent is effective for direct immunofluorescence staining of human tissue for flow cytometric analysis using I test

for 10⁶ cells.

Presentation: liquid

Source: Supernatant proceeding from an *in vitro* cell

culture of a cell hybridoma.

Purification: Affinity chromatography.

ANTIGEN DETAILS

Large description: The CD3 monoclonal antibody is directed against the CD3- antigen (T3-antigen), which is expressed on human T lymphocytes. The monoclonal antibody reacts with 80-90% human peripheral T lymphocytes and medullary thymocytes. The monoclonal antibody does not react with B-cells, monocytes, granulocytes and platelets.

The monoclonal antibody is directed against the CD19- antigen (B4-antigen), which is expressed on human B lymphocytes. The monoclonal antibody is B lineage-specific and reacts with early B-cell precursors, pre-pre-B-cells, pre-B- cells, B-cells, intermediate B-cells, mature B-cells and some plasmacytoid cells, also reacts with pre-B-cell- lines, B lymphoblastoid cell-lines and Burkitt cell- lines, and with 50% of myeloma cell-lines. Plasma cells were found to be negative. [1-4]

Please, refer to www.immunostep.com technical support for more information.

WARRANTY

Warranted only to conform to the quantity and contents stated on the label or in the product labelling at the time of delivery to the customer. Immunostep disclaims hereby other warranties. Immunostep's sole liability is limited to either the replacement of the products or refund of the purchase price.

REFERENCES

- Pezzutto A, Dörken B, Rabinovitch PS, Ledbetter JA, Moldenhauer G, Clark EA. CD19 monoclonal antibody HD37 inhibits anti-immunoglobulininduced B cell activation and proliferation. J Immunol 1987;138:2793-9.
- Sato S, Tedder TF. BC3. CD19 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 133-5.
- 3. Braylan RC, Orfao A, Borowitz MJ, Davis BH. Optimal number of reagents required to evaluate hematolymphoid neoplasias: results of an international consensus meeting. Cytometry 2001;46:23-7.

MANUFACTURED BY

Immunostep S.L



Avda. Universidad de Coimbra, s/n Cancer Research Center (CIC) Campus Miguel de Unamuno 37007 Salamanca (Spain)

Tel. (+34) 923 294 827

www.immunostep.com

Revision N° 3 Emission date: 16/02/2021 HT-H3-19-3