


## Glyco-Diosgenin

(GDN)

Cat # NB-19-0066



### Description

Glyco-diosgenin is a synthetic surfactant and detergent for extracting proteins from membranes for structure and function studies, and single-particle cryo-electron microscopy (cryoEM) studies of membrane proteins. Synthetic drop-in substitute for Digitonin.

### Product Information

<b>Codes</b>	NB-19-0066-500MG, NB-19-0066-1G, NB-19-0066-5G, NB-19-0066-25G
<b>Sizes</b>	500mg, 1g, 5g, 25g
<b>Synonyms</b>	Glyco-diosgenin, GDN
<b>CAS number</b>	1402423-29-3
<b>Formula</b>	C <sub>56</sub> H <sub>92</sub> O <sub>25</sub>
<b>Molecular Weight</b>	1165.31
<b>Purity (HPLC)</b>	>97%
<b>Note</b>	For research use only.

*For reference only  
For Research Use Only. Not for Diagnostic or Therapeutic Use.*