



Human SIRT1 (NAD-dependent protein deacetylase sirtuin-1) Sandwich ELISA Kit

Cat# NB-22-43265-1

GENERAL INFORMATION


Product Type	Kits
Short Description	This SIRT1 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	0.19ng/mL
Detection limit	0.31~20ng/mL
Assay Time	3.5h
Storage Instruction	If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual
Note	FOR RESEARCH USE ONLY (RUO)

TARGET INFORMATION

Sample Types	Serum, plasma and other biological fluids
Gene	SIRT1
Symbol	
Uniprot ID	SIR1_HUMAN



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic other medical purposes.