



Bovine CASP3 (Caspase 3) Sandwich ELISA Kit

Cat# NB-22-64210-48

Cat# NB-22-64210-96

GENERAL INFORMATION


Product Type	Kits
Short Description	This CASP3 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in bovine cell culture supernatant, serum and plasma (EDTA, citrate, heparin)
Applications	ELISA
Reactivity	Bovine

PRODUCT PROPERTIES

Detection	0.156-10ng/mL
Sensitivity	0.071ng/mL
Detection limit	0.156-10ng/mL
Assay Time	4.5 hrs
Storage Instruction	Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit
Note	FOR RESEARCH USE ONLY (RUO)

TARGET INFORMATION

Sample Types	serum, plasma, tissue homogenates or other biological fluids
Gene	
Symbol	
Uniprot ID	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic other medical purposes.

Neo Biotech

20231129

74 rue des Suisses – 92000 Nanterre