

Anti-LRP1 antibody

Cat #NB-22-3952 (100 µl)

Description

Rabbit polyclonal to LRP1.

Product Information

Model	NB-22-3952
Host	Rabbit
Reactivity	Human, Mouse
Applications	ELISA, IHC, WB
Immunogen	Synthesized peptide derived from human LRP1 around the non-phosphorylation site of S4520.
Immunogen Region	4460-4540 aa
Gene ID	4035
Gene Symbol	LRP1
Dilution range	WB 1:500-1:2000IHC 1:100-1:300ELISA 1:40000
Specificity	LRP1 Polyclonal Antibody detects endogenous levels of LRP1 protein.
Tissue Specificity	Most abundant in liver, brain and lung.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Note	For Research Use Only (RUO).
Protein Name	Prolow-density lipoprotein receptor-related protein 1 LRP-1 Alpha-2-macroglobulin receptor A2MR Apolipoprotein E receptor APOER CD antigen CD91 Low-density lipoprotein receptor-related protein 1 85-kDa subunit