


CL2A-SN-38

Cat. No.:	HY-128946
CAS No.:	1279680-68-0
Molecular Formula:	C ₇₃ H ₉₇ N ₁₁ O ₂₂
Molecular Weight:	1480.61
Target:	Drug-Linker Conjugates for ADC
Pathway:	Antibody-drug Conjugate/ADC Related
Storage:	-20°C, stored under nitrogen
	* In solvent : -80°C, 6 months; -20°C, 1 month (stored under nitrogen)



SOLVENT & SOLUBILITY

In Vitro	DMSO : 100 mg/mL (67.54 mM; Need ultrasonic)					
	Preparing Stock Solutions	Solvent Concentration	Mass	1 mg	5 mg	10 mg
		1 mM		0.6754 mL	3.3770 mL	6.7540 mL
		5 mM		0.1351 mL	0.6754 mL	1.3508 mL
		10 mM		0.0675 mL	0.3377 mL	0.6754 mL
Please refer to the solubility information to select the appropriate solvent.						
In Vivo	1. Add each solvent one by one: 10% DMSO >> 40% PEG300 >> 5% Tween-80 >> 45% saline Solubility: 2.08 mg/mL (1.40 mM); Clear solution; Need ultrasonic					
	2. Add each solvent one by one: 10% DMSO >> 90% (20% SBE-β-CD in saline) Solubility: ≥ 2.08 mg/mL (1.40 mM); Clear solution					
	3. Add each solvent one by one: 10% DMSO >> 90% corn oil Solubility: ≥ 2.08 mg/mL (1.40 mM); Clear solution					

BIOLOGICAL ACTIVITY

Description	CL2A-SN-38 is a agent-linker conjugate composed of a potent a DNA Topoisomerase I inhibitor SN-38 and a linker CL2A to make antibody agent conjugate (ADC). CL2A-SN-38 provides significant and specific antitumor effects against a range of human solid tumor types. CL2A-SN-38 uses hydrolyzable linker to deliver active agents within tumor cells and in the tumor microenvironment, resulting in bystander effects ^[1] .
IC ₅₀ & Target	Camptothecins
In Vivo	CL2A-SN-38 coupled to the anti-Trop-2 humanized antibody hRS7 (intraperitoneal injection, 0.2 or 0.4 mg/kg, twice weekly, 4 weeks) can both significantly inhibit tumor growth in mice with human colon (COLO 205) or pancreatic (Capan-1) tumour

xenografts^[2].

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

CUSTOMER VALIDATION

- Nano Today. 2024 Oct.
- J Med Chem. 2025 Jul 20.
- Pharmaceutics. 2025 Jul 25;17(8):967.
- AAPS J. 2022 May 27;24(4):70.

See more customer validations on www.MedChemExpress.com

REFERENCES

[1]. Thomas M Cardillo, et al. Humanized anti-Trop-2 IgG-SN-38 conjugate for effective treatment of diverse epithelial cancers: preclinical studies in human cancer xenograft models and monkeys. Clin Cancer Res. 2011 May 15;17(10):3157-69.

[2]. Cardillo TM, et al. Humanized anti-Trop-2 IgG-SN-38 conjugate for effective treatment of diverse epithelial cancers: preclinical studies in human cancer xenograft models and monkeys. Clin Cancer Res. 2011 May 15;17(10):3157-69.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite F, Monmouth Junction, NJ 08852, USA