



## Monoclonal anti-Cyclin L1 (CCNLI)

**NB-01-0020**

**Volume:** 100 µl

**Clones:** 3CY-2F10

**Immunogen:** synthetic peptide: (C)-SKHHGGSRSRGHGRHRR corresponding to aa 511 to 526 of human Cyclin L1.

**Murine Isotype:** IgG2a, k

**Specificity:** human, mouse and rat Cyclin L1.

### Applications:

- ELISA
- Western Blot
- Immunocytochemistry
- Immunohistochemistry

**Recommended dilutions:** 1/500-1/5000

### Storage:

- Store at -20°C (up to 3 years)
- Avoid repeat freezing and thawing cycles. Does not contain any preservative

### References:

D. Muller, R. Millon, S. The´obald, T. Hussenet, B. Wasylyk, S. du Manoir and J. Abecassis. Cyclin L1 (CCNLI) gene alterations in human head and neck squamous cell carcinoma. *British Journal of Cancer* (2006) 94, 1041 – 1044.