

Secondary Antibody Promotion 2023
sample list

Secondary Antibody Promotion 2023

Cat.no,	Description	Isotype	Clone	Size
S1010-02	Goat Anti-Mouse Ig, Human ads-FITC	IgG	Polyclonal	20µl
S1010-03	Goat Anti-Mouse Ig, Human ads-TRITC	IgG	Polyclonal	20µl
S1010-04	Goat Anti-Mouse Ig, Human ads-AP	IgG	Polyclonal	20µl
S1010-05	Goat Anti-Mouse Ig, Human ads-HRP	IgG	Polyclonal	20µl
S1010-06	Goat Anti-Mouse Ig, Human ads-BGAL	IgG	Polyclonal	20µl
S1010-07	Goat Anti-Mouse Ig, Human ads-TXRD	IgG	Polyclonal	20µl
S1010-08	Goat Anti-Mouse Ig, Human ads-BIOT	IgG	Polyclonal	20µl
S1010-09	Goat Anti-Mouse Ig, Human ads-PE	IgG	Polyclonal	20µl
S1010-30	Goat Anti-Mouse Ig, Human ads-AF488	IgG	Polyclonal	20µl
S1010-31	Goat Anti-Mouse Ig, Human ads-AF647	IgG	Polyclonal	20µl
S1010-32	Goat Anti-Mouse Ig, Human ads-AF555	IgG	Polyclonal	20µl
S1012-02	Goat F(ab) ₂ Anti-Mouse Ig, Human ads-FITC	IgG	Polyclonal	20µl
S1012-04	Goat F(ab) ₂ Anti-Mouse Ig, Human ads-AP	IgG	Polyclonal	20µl
S1012-05	Goat F(ab) ₂ Anti-Mouse Ig, Human ads-HRP	IgG	Polyclonal	20µl
S1012-08	Goat F(ab) ₂ Anti-Mouse Ig, Human ads-BIOT	IgG	Polyclonal	20µl
S1012-09	Goat F(ab) ₂ Anti-Mouse Ig, Human ads-PE	IgG	Polyclonal	20µl
S1012-11	Goat F(ab) ₂ Anti-Mouse Ig, Human ads-APC	IgG	Polyclonal	20µl
S1012-13	Goat F(ab) ₂ Anti-Mouse Ig, Human ads-SPRD	IgG	Polyclonal	20µl
S1012-15	Goat F(ab) ₂ Anti-Mouse Ig, Human ads-CY5	IgG	Polyclonal	20µl
S1012-30	Goat F(ab) ₂ Anti-Mouse Ig, Human ads-AF488	IgG	Polyclonal	20µl
S1012-31	Goat F(ab) ₂ Anti-Mouse Ig, Human ads-AF647	IgG	Polyclonal	20µl
S1012-32	Goat F(ab) ₂ Anti-Mouse Ig, Human ads-AF555	IgG	Polyclonal	20µl
S1013-02	Goat Anti-Mouse IgG Fc, Human/Bovine/Horse SP ads-FITC	IgG	Polyclonal	20µl
S1013-05	Goat Anti-Mouse IgG Fc, Human/Bovine/Horse SP ads-HRP	IgG	Polyclonal	20µl
S1015-02	Goat Anti-Mouse IgG Fab-FITC	IgG	Polyclonal	20µl
S1015-04	Goat Anti-Mouse IgG Fab-AP	IgG	Polyclonal	20µl
S1015-05	Goat Anti-Mouse IgG Fab-HRP	IgG	Polyclonal	20µl
S1015-08	Goat Anti-Mouse IgG Fab-BIOT	IgG	Polyclonal	20µl
S1015-31	Goat Anti-Mouse IgG Fab-AF647	IgG	Polyclonal	20µl
S1020-02	Goat Anti-Mouse IgM, Human ads-FITC	IgG	Polyclonal	20µl
S1020-03	Goat Anti-Mouse IgM, Human ads-TRITC	IgG	Polyclonal	20µl
S1020-04	Goat Anti-Mouse IgM, Human ads-AP	IgG	Polyclonal	20µl
S1020-05	Goat Anti-Mouse IgM, Human ads-HRP	IgG	Polyclonal	20µl
S1020-06	Goat Anti-Mouse IgM, Human ads-BGAL	IgG	Polyclonal	20µl
S1020-07	Goat Anti-Mouse IgM, Human ads-TXRD	IgG	Polyclonal	20µl
S1020-08	Goat Anti-Mouse IgM, Human ads-BIOT	IgG	Polyclonal	20µl
S1020-09	Goat Anti-Mouse IgM, Human ads-PE	IgG	Polyclonal	20µl
S1020-11	Goat Anti-Mouse IgM, Human ads-APC	IgG	Polyclonal	20µl
S1020-15	Goat Anti-Mouse IgM, Human ads-CY5	IgG	Polyclonal	20µl
S1020-17	Goat Anti-Mouse IgM, Human ads-PE/CY7	IgG	Polyclonal	20µl
S1020-30	Goat Anti-Mouse IgM, Human ads-AF488	IgG	Polyclonal	20µl
S1020-31	Goat Anti-Mouse IgM, Human ads-AF647	IgG	Polyclonal	20µl
S1020-32	Goat Anti-Mouse IgM, Human ads-AF555	IgG	Polyclonal	20µl
S1021-02	Goat Anti-Mouse IgM-FITC	IgG	Polyclonal	20µl
S1021-03	Goat Anti-Mouse IgM-TRITC	IgG	Polyclonal	20µl
S1021-04	Goat Anti-Mouse IgM-AP	IgG	Polyclonal	20µl
S1021-05	Goat Anti-Mouse IgM-HRP	IgG	Polyclonal	20µl
S1021-07	Goat Anti-Mouse IgM-TXRD	IgG	Polyclonal	20µl
S1021-08	Goat Anti-Mouse IgM-BIOT	IgG	Polyclonal	20µl
S1021-09	Goat Anti-Mouse IgM-PE	IgG	Polyclonal	20µl
S1021-30	Goat Anti-Mouse IgM-AF488	IgG	Polyclonal	20µl
S1021-31	Goat Anti-Mouse IgM-AF647	IgG	Polyclonal	20µl
S1021-32	Goat Anti-Mouse IgM-AF555	IgG	Polyclonal	20µl
S1022-02	Goat F(ab) ₂ Anti-Mouse IgM, Human ads-FITC	IgG	Polyclonal	20µl
S1022-05	Goat F(ab) ₂ Anti-Mouse IgM, Human ads-HRP	IgG	Polyclonal	20µl
S1022-08	Goat F(ab) ₂ Anti-Mouse IgM, Human ads-BIOT	IgG	Polyclonal	20µl
S1022-09	Goat F(ab) ₂ Anti-Mouse IgM, Human ads-PE	IgG	Polyclonal	20µl
S1022-30	Goat F(ab) ₂ Anti-Mouse IgM, Human ads-AF488	IgG	Polyclonal	20µl
S1022-31	Goat F(ab) ₂ Anti-Mouse IgM, Human ads-AF647	IgG	Polyclonal	20µl
S1022-32	Goat F(ab) ₂ Anti-Mouse IgM, Human ads-AF555	IgG	Polyclonal	20µl
S1030-02	Goat Anti-Mouse IgG, Human ads-FITC	IgG	Polyclonal	20µl
S1030-03	Goat Anti-Mouse IgG, Human ads-TRITC	IgG	Polyclonal	20µl
S1030-04	Goat Anti-Mouse IgG, Human ads-AP	IgG	Polyclonal	20µl
S1030-05	Goat Anti-Mouse IgG, Human ads-HRP	IgG	Polyclonal	20µl
S1030-06	Goat Anti-Mouse IgG, Human ads-BGAL	IgG	Polyclonal	20µl
S1030-07	Goat Anti-Mouse IgG, Human ads-TXRD	IgG	Polyclonal	20µl
S1030-08	Goat Anti-Mouse IgG, Human ads-BIOT	IgG	Polyclonal	20µl

Secondary Antibody Promotion 2023
sample list

Cat.no.	Description	Isotype	Clone	Size
S1030-09	Goat Anti-Mouse IgG, Human ads-PE	IgG	Polyclonal	20µl
S1030-15	Goat Anti-Mouse IgG, Human ads-CY5	IgG	Polyclonal	20µl
S1030-30	Goat Anti-Mouse IgG, Human ads-AF488	IgG	Polyclonal	20µl
S1030-31	Goat Anti-Mouse IgG, Human ads-AF647	IgG	Polyclonal	20µl
S1030-32	Goat Anti-Mouse IgG, Human ads-AF555	IgG	Polyclonal	20µl
S1031-02	Goat Anti-Mouse IgG(H+L), Human ads-FITC	IgG	Polyclonal	20µl
S1031-03	Goat Anti-Mouse IgG(H+L), Human ads-TRITC	IgG	Polyclonal	20µl
S1031-04	Goat Anti-Mouse IgG(H+L), Human ads-AP	IgG	Polyclonal	20µl
S1031-05	Goat Anti-Mouse IgG(H+L), Human ads-HRP	IgG	Polyclonal	20µl
S1031-07	Goat Anti-Mouse IgG(H+L), Human ads-TXRD	IgG	Polyclonal	20µl
S1031-08	Goat Anti-Mouse IgG(H+L), Human ads-BIOT	IgG	Polyclonal	20µl
S1031-09	Goat Anti-Mouse IgG(H+L), Human ads-PE	IgG	Polyclonal	20µl
S1031-11	Goat Anti-Mouse IgG(H+L), Human ads-APC	IgG	Polyclonal	20µl
S1031-15	Goat Anti-Mouse IgG(H+L), Human ads-CY5	IgG	Polyclonal	20µl
S1031-30	Goat Anti-Mouse IgG(H+L), Human ads-AF488	IgG	Polyclonal	20µl
S1031-31	Goat Anti-Mouse IgG(H+L), Human ads-AF647	IgG	Polyclonal	20µl
S1031-32	Goat Anti-Mouse IgG(H+L), Human ads-AF555	IgG	Polyclonal	20µl
S1032-02	Goat F(ab) ₂ Anti-Mouse IgG(H+L), Human ads-FITC	IgG	Polyclonal	20µl
S1032-04	Goat F(ab) ₂ Anti-Mouse IgG(H+L), Human ads-AP	IgG	Polyclonal	20µl
S1032-05	Goat F(ab) ₂ Anti-Mouse IgG(H+L), Human ads-HRP	IgG	Polyclonal	20µl
S1032-08	Goat F(ab) ₂ Anti-Mouse IgG(H+L), Human ads-BIOT	IgG	Polyclonal	20µl
S1032-09	Goat F(ab) ₂ Anti-Mouse IgG(H+L), Human ads-PE	IgG	Polyclonal	20µl
S1032-11	Goat F(ab) ₂ Anti-Mouse IgG(H+L), Human ads-APC	IgG	Polyclonal	20µl
S1032-13	Goat F(ab) ₂ Anti-Mouse IgG(H+L), Human ads-SPRD	IgG	Polyclonal	20µl
S1032-15	Goat F(ab) ₂ Anti-Mouse IgG(H+L), Human ads-CY5	IgG	Polyclonal	20µl
S1032-30	Goat F(ab) ₂ Anti-Mouse IgG(H+L), Human ads-AF488	IgG	Polyclonal	20µl
S1032-32	Goat F(ab) ₂ Anti-Mouse IgG(H+L), Human ads-AF555	IgG	Polyclonal	20µl
S1033-02	Goat Anti-Mouse IgG Fc-FITC	IgG	Polyclonal	20µl
S1033-04	Goat Anti-Mouse IgG Fc-AP	IgG	Polyclonal	20µl
S1033-05	Goat Anti-Mouse IgG Fc-HRP	IgG	Polyclonal	20µl
S1033-08	Goat Anti-Mouse IgG Fc-BIOT	IgG	Polyclonal	20µl
S1033-30	Goat Anti-Mouse IgG Fc-AF488	IgG	Polyclonal	20µl
S1033-31	Goat Anti-Mouse IgG Fc-AF647	IgG	Polyclonal	20µl
S1034-02	Goat Anti-Mouse IgG(H+L), Rat ads-FITC	IgG	Polyclonal	20µl
S1034-03	Goat Anti-Mouse IgG(H+L), Rat ads-TRITC	IgG	Polyclonal	20µl
S1034-04	Goat Anti-Mouse IgG(H+L), Rat ads-AP	IgG	Polyclonal	20µl
S1034-05	Goat Anti-Mouse IgG(H+L), Rat ads-HRP	IgG	Polyclonal	20µl
S1034-07	Goat Anti-Mouse IgG(H+L), Rat ads-TXRD	IgG	Polyclonal	20µl
S1034-08	Goat Anti-Mouse IgG(H+L), Rat ads-BIOT	IgG	Polyclonal	20µl
S1034-09	Goat Anti-Mouse IgG(H+L), Rat ads-PE	IgG	Polyclonal	20µl
S1034-32	Goat Anti-Mouse IgG(H+L), Rat ads-AF555	IgG	Polyclonal	20µl
S1036-02	Goat Anti-Mouse IgG(H+L)-FITC	IgG	Polyclonal	20µl
S1036-04	Goat Anti-Mouse IgG(H+L)-AP	IgG	Polyclonal	20µl
S1036-05	Goat Anti-Mouse IgG(H+L)-HRP	IgG	Polyclonal	20µl
S1036-08	Goat Anti-Mouse IgG(H+L)-BIOT	IgG	Polyclonal	20µl
S1036-09	Goat Anti-Mouse IgG(H+L)-PE	IgG	Polyclonal	20µl
S1036-30	Goat Anti-Mouse IgG(H+L)-AF488	IgG	Polyclonal	20µl
S1036-31	Goat Anti-Mouse IgG(H+L)-AF647	IgG	Polyclonal	20µl
S1036-32	Goat Anti-Mouse IgG(H+L)-AF555	IgG	Polyclonal	20µl
S1037-02	Goat Anti-Mouse IgG(H+L), Human/Bovine/Horse SP ads-FITC	IgG	Polyclonal	20µl
S1037-04	Goat Anti-Mouse IgG(H+L), Human/Bovine/Horse SP ads-AP	IgG	Polyclonal	20µl
S1037-05	Goat Anti-Mouse IgG(H+L), Human/Bovine/Horse SP ads-HRP	IgG	Polyclonal	20µl
S1037-08	Goat Anti-Mouse IgG(H+L), Human/Bovine/Horse SP ads-BIOT	IgG	Polyclonal	20µl
S1037-32	Goat Anti-Mouse IgG(H+L), Human/Bovine/Horse SP ads-AF555	IgG	Polyclonal	20µl
S1038-02	Goat Anti-Mouse IgG(H+L), Multi-Species SP ads-FITC	IgG	Polyclonal	20µl
S1038-04	Goat Anti-Mouse IgG(H+L), Multi-Species SP ads-AP	IgG	Polyclonal	20µl
S1038-05	Goat Anti-Mouse IgG(H+L), Multi-Species SP ads-HRP	IgG	Polyclonal	20µl
S1038-08	Goat Anti-Mouse IgG(H+L), Multi-Species SP ads-BIOT	IgG	Polyclonal	20µl
S1038-09	Goat Anti-Mouse IgG(H+L), Multi-Species SP ads-PE	IgG	Polyclonal	20µl
S1038-30	Goat Anti-Mouse IgG(H+L), Multi-Species SP ads-AF488	IgG	Polyclonal	20µl
S1038-31	Goat Anti-Mouse IgG(H+L), Multi-Species SP ads-AF647	IgG	Polyclonal	20µl
S1038-32	Goat Anti-Mouse IgG(H+L), Multi-Species SP ads-AF555	IgG	Polyclonal	20µl
S1050-02	Goat Anti-Mouse Kappa-FITC	IgG	Polyclonal	20µl
S1050-03	Goat Anti-Mouse Kappa-TRITC	IgG	Polyclonal	20µl
S1050-04	Goat Anti-Mouse Kappa-AP	IgG	Polyclonal	20µl
S1050-05	Goat Anti-Mouse Kappa-HRP	IgG	Polyclonal	20µl
S1050-07	Goat Anti-Mouse Kappa-TXRD	IgG	Polyclonal	20µl
S1050-08	Goat Anti-Mouse Kappa-BIOT	IgG	Polyclonal	20µl
S1050-09	Goat Anti-Mouse Kappa-PE	IgG	Polyclonal	20µl
S1050-31	Goat Anti-Mouse Kappa-AF647	IgG	Polyclonal	20µl
S1050-32	Goat Anti-Mouse Kappa-AF555	IgG	Polyclonal	20µl
S1052-05	Goat F(ab) ₂ Anti-Mouse Kappa-HRP	IgG	Polyclonal	20µl
S1052-08	Goat F(ab) ₂ Anti-Mouse Kappa-BIOT	IgG	Polyclonal	20µl
S1060-02	Goat Anti-Mouse Lambda-FITC	IgG	Polyclonal	20µl
S1060-03	Goat Anti-Mouse Lambda-TRITC	IgG	Polyclonal	20µl

Secondary Antibody Promotion 2023
sample list

Cat.no.	Description	Isotype	Clone	Size
S1060-04	Goat Anti-Mouse Lambda-AP	IgG	Polyclonal	20µl
S1060-05	Goat Anti-Mouse Lambda-HRP	IgG	Polyclonal	20µl
S1060-07	Goat Anti-Mouse Lambda-TXRD	IgG	Polyclonal	20µl
S1060-08	Goat Anti-Mouse Lambda-BIOT	IgG	Polyclonal	20µl
S1060-09	Goat Anti-Mouse Lambda-PE	IgG	Polyclonal	20µl
S1060-32	Goat Anti-Mouse Lambda-AF555	IgG	Polyclonal	20µl
S1070-02	Goat Anti-Mouse IgG ₁ , Human ads-FITC	IgG	Polyclonal	20µl
S1070-03	Goat Anti-Mouse IgG ₁ , Human ads-TRITC	IgG	Polyclonal	20µl
S1070-04	Goat Anti-Mouse IgG ₁ , Human ads-AP	IgG	Polyclonal	20µl
S1070-05	Goat Anti-Mouse IgG ₁ , Human ads-HRP	IgG	Polyclonal	20µl
S1070-07	Goat Anti-Mouse IgG ₁ , Human ads-TXRD	IgG	Polyclonal	20µl
S1070-08	Goat Anti-Mouse IgG ₁ , Human ads-BIOT	IgG	Polyclonal	20µl
S1070-09	Goat Anti-Mouse IgG ₁ , Human ads-PE	IgG	Polyclonal	20µl
S1070-11	Goat Anti-Mouse IgG ₁ , Human ads-APC	IgG	Polyclonal	20µl
S1070-15	Goat Anti-Mouse IgG ₁ , Human ads-CY5	IgG	Polyclonal	20µl
S1070-19	Goat Anti-Mouse IgG ₁ , Human ads-APC/CY7	IgG	Polyclonal	20µl
S1070-30	Goat Anti-Mouse IgG ₁ , Human ads-AF488	IgG	Polyclonal	20µl
S1070-31	Goat Anti-Mouse IgG ₁ , Human ads-AF647	IgG	Polyclonal	20µl
S1070-32	Goat Anti-Mouse IgG ₁ , Human ads-AF555	IgG	Polyclonal	20µl
S1071-02	Goat Anti-Mouse IgG ₁ -FITC	IgG	Polyclonal	20µl
S1071-04	Goat Anti-Mouse IgG ₁ -AP	IgG	Polyclonal	20µl
S1071-05	Goat Anti-Mouse IgG ₁ -HRP	IgG	Polyclonal	20µl
S1071-08	Goat Anti-Mouse IgG ₁ -BIOT	IgG	Polyclonal	20µl
S1071-30	Goat Anti-Mouse IgG ₁ -AF488	IgG	Polyclonal	20µl
S1071-31	Goat Anti-Mouse IgG ₁ -AF647	IgG	Polyclonal	20µl
S1071-32	Goat Anti-Mouse IgG ₁ -AF555	IgG	Polyclonal	20µl
S1072-02	Goat F(ab) ₂ Anti-Mouse IgG ₁ , Human ads-FITC	IgG	Polyclonal	20µl
S1072-08	Goat F(ab) ₂ Anti-Mouse IgG ₁ , Human ads-BIOT	IgG	Polyclonal	20µl
S1072-09	Goat F(ab) ₂ Anti-Mouse IgG ₁ , Human ads-PE	IgG	Polyclonal	20µl
S1072-30	Goat F(ab) ₂ Anti-Mouse IgG ₁ , Human ads-AF488	IgG	Polyclonal	20µl
S1072-31	Goat F(ab) ₂ Anti-Mouse IgG ₁ , Human ads-AF647	IgG	Polyclonal	20µl
S1072-32	Goat F(ab) ₂ Anti-Mouse IgG ₁ , Human ads-AF555	IgG	Polyclonal	20µl
S1073-02	Goat Anti-Mouse IgG ₁ , Human/Bovine/Horse SP ads-FITC	IgG	Polyclonal	20µl
S1073-04	Goat Anti-Mouse IgG ₁ , Human/Bovine/Horse SP ads-AP	IgG	Polyclonal	20µl
S1073-05	Goat Anti-Mouse IgG ₁ , Human/Bovine/Horse SP ads-HRP	IgG	Polyclonal	20µl
S1073-08	Goat Anti-Mouse IgG ₁ , Human/Bovine/Horse SP ads-BIOT	IgG	Polyclonal	20µl
S1073-30	Goat Anti-Mouse IgG ₁ , Human/Bovine/Horse SP ads-AF488	IgG	Polyclonal	20µl
S1073-31	Goat Anti-Mouse IgG ₁ , Human/Bovine/Horse SP ads-AF647	IgG	Polyclonal	20µl
S1073-32	Goat Anti-Mouse IgG ₁ , Human/Bovine/Horse SP ads-AF555	IgG	Polyclonal	20µl
S1077-02	Goat Anti-Mouse IgG _{2c} , Human/Bovine/Horse SP ads-FITC	IgG	Polyclonal	20µl
S1077-04	Goat Anti-Mouse IgG _{2c} , Human/Bovine/Horse SP ads-AP	IgG	Polyclonal	20µl
S1077-05	Goat Anti-Mouse IgG _{2c} , Human/Bovine/Horse SP ads-HRP	IgG	Polyclonal	20µl
S1077-08	Goat Anti-Mouse IgG _{2c} , Human/Bovine/Horse SP ads-BIOT	IgG	Polyclonal	20µl
S1077-30	Goat Anti-Mouse IgG _{2c} , Human/Bovine/Horse SP ads-AF488	IgG	Polyclonal	20µl
S1077-31	Goat Anti-Mouse IgG _{2c} , Human/Bovine/Horse SP ads-AF647	IgG	Polyclonal	20µl
S1077-32	Goat Anti-Mouse IgG _{2c} , Human/Bovine/Horse SP ads-AF555	IgG	Polyclonal	20µl
S1078-02	Goat Anti-Mouse IgG _{2c} -FITC	IgG	Polyclonal	20µl
S1078-04	Goat Anti-Mouse IgG _{2c} -AP	IgG	Polyclonal	20µl
S1078-05	Goat Anti-Mouse IgG _{2c} -HRP	IgG	Polyclonal	20µl
S1078-08	Goat Anti-Mouse IgG _{2c} -BIOT	IgG	Polyclonal	20µl
S1079-02	Goat Anti-Mouse IgG _{2c} , Human ads-FITC	IgG	Polyclonal	20µl
S1079-04	Goat Anti-Mouse IgG _{2c} , Human ads-AP	IgG	Polyclonal	20µl
S1079-05	Goat Anti-Mouse IgG _{2c} , Human ads-HRP	IgG	Polyclonal	20µl
S1079-08	Goat Anti-Mouse IgG _{2c} , Human ads-BIOT	IgG	Polyclonal	20µl
S1079-09	Goat Anti-Mouse IgG _{2c} , Human ads-PE	IgG	Polyclonal	20µl
S1079-11	Goat Anti-Mouse IgG _{2c} , Human ads-APC	IgG	Polyclonal	20µl
S1079-30	Goat Anti-Mouse IgG _{2c} , Human ads-AF488	IgG	Polyclonal	20µl
S1079-31	Goat Anti-Mouse IgG _{2c} , Human ads-AF647	IgG	Polyclonal	20µl
S1079-32	Goat Anti-Mouse IgG _{2c} , Human ads-AF555	IgG	Polyclonal	20µl
S1080-02	Goat Anti-Mouse IgG _{2a} , Human ads-FITC	IgG	Polyclonal	20µl
S1080-03	Goat Anti-Mouse IgG _{2a} , Human ads-TRITC	IgG	Polyclonal	20µl
S1080-04	Goat Anti-Mouse IgG _{2a} , Human ads-AP	IgG	Polyclonal	20µl
S1080-05	Goat Anti-Mouse IgG _{2a} , Human ads-HRP	IgG	Polyclonal	20µl
S1080-07	Goat Anti-Mouse IgG _{2a} , Human ads-TXRD	IgG	Polyclonal	20µl
S1080-08	Goat Anti-Mouse IgG _{2a} , Human ads-BIOT	IgG	Polyclonal	20µl
S1080-09	Goat Anti-Mouse IgG _{2a} , Human ads-PE	IgG	Polyclonal	20µl
S1080-10	Goat Anti-Mouse IgG _{2a} , Human ads-PE/TXRD	IgG	Polyclonal	20µl
S1080-11	Goat Anti-Mouse IgG _{2a} , Human ads-APC	IgG	Polyclonal	20µl
S1080-13	Goat Anti-Mouse IgG _{2a} , Human ads-SPRD	IgG	Polyclonal	20µl
S1080-15	Goat Anti-Mouse IgG _{2a} , Human ads-CY5	IgG	Polyclonal	20µl
S1080-17	Goat Anti-Mouse IgG _{2a} , Human ads-PE/CY7	IgG	Polyclonal	20µl
S1080-19	Goat Anti-Mouse IgG _{2a} , Human ads-APC/CY7	IgG	Polyclonal	20µl
S1080-30	Goat Anti-Mouse IgG _{2a} , Human ads-AF488	IgG	Polyclonal	20µl
S1080-31	Goat Anti-Mouse IgG _{2a} , Human ads-AF647	IgG	Polyclonal	20µl
S1080-32	Goat Anti-Mouse IgG _{2a} , Human ads-AF555	IgG	Polyclonal	20µl

Secondary Antibody Promotion 2023
sample list

Cat.no.	Description	Isotype	Clone	Size
S1081-02	Goat Anti-Mouse IgG _{2a} -FITC	IgG	Polyclonal	20µl
S1081-04	Goat Anti-Mouse IgG _{2a} -AP	IgG	Polyclonal	20µl
S1081-05	Goat Anti-Mouse IgG _{2a} -HRP	IgG	Polyclonal	20µl
S1081-08	Goat Anti-Mouse IgG _{2a} -BIOT	IgG	Polyclonal	20µl
S1081-30	Goat Anti-Mouse IgG _{2a} -AF488	IgG	Polyclonal	20µl
S1081-31	Goat Anti-Mouse IgG _{2a} -AF647	IgG	Polyclonal	20µl
S1081-32	Goat Anti-Mouse IgG _{2a} -AF555	IgG	Polyclonal	20µl
S1082-02	Goat F(ab) ₂ Anti-Mouse IgG _{2a} , Human ads-FITC	IgG	Polyclonal	20µl
S1082-08	Goat F(ab) ₂ Anti-Mouse IgG _{2a} , Human ads-BIOT	IgG	Polyclonal	20µl
S1082-09	Goat F(ab) ₂ Anti-Mouse IgG _{2a} , Human ads-PE	IgG	Polyclonal	20µl
S1082-30	Goat F(ab) ₂ Anti-Mouse IgG _{2a} , Human ads-AF488	IgG	Polyclonal	20µl
S1082-31	Goat F(ab) ₂ Anti-Mouse IgG _{2a} , Human ads-AF647	IgG	Polyclonal	20µl
S1082-32	Goat F(ab) ₂ Anti-Mouse IgG _{2a} , Human ads-AF555	IgG	Polyclonal	20µl
S1083-02	Goat Anti-Mouse IgG _{2a} , Human/Bovine/Horse SP ads-FITC	IgG	Polyclonal	20µl
S1083-04	Goat Anti-Mouse IgG _{2a} , Human/Bovine/Horse SP ads-AP	IgG	Polyclonal	20µl
S1083-05	Goat Anti-Mouse IgG _{2a} , Human/Bovine/Horse SP ads-HRP	IgG	Polyclonal	20µl
S1083-08	Goat Anti-Mouse IgG _{2a} , Human/Bovine/Horse SP ads-BIOT	IgG	Polyclonal	20µl
S1090-02	Goat Anti-Mouse IgG _{2b} , Human ads-FITC	IgG	Polyclonal	20µl
S1090-03	Goat Anti-Mouse IgG _{2b} , Human ads-TRITC	IgG	Polyclonal	20µl
S1090-04	Goat Anti-Mouse IgG _{2b} , Human ads-AP	IgG	Polyclonal	20µl
S1090-05	Goat Anti-Mouse IgG _{2b} , Human ads-HRP	IgG	Polyclonal	20µl
S1090-07	Goat Anti-Mouse IgG _{2b} , Human ads-TXRD	IgG	Polyclonal	20µl
S1090-08	Goat Anti-Mouse IgG _{2b} , Human ads-BIOT	IgG	Polyclonal	20µl
S1090-09	Goat Anti-Mouse IgG _{2b} , Human ads-PE	IgG	Polyclonal	20µl
S1090-11	Goat Anti-Mouse IgG _{2b} , Human ads-APC	IgG	Polyclonal	20µl
S1090-13	Goat Anti-Mouse IgG _{2b} , Human ads-SPRD	IgG	Polyclonal	20µl
S1090-15	Goat Anti-Mouse IgG _{2b} , Human ads-CY5	IgG	Polyclonal	20µl
S1090-17	Goat Anti-Mouse IgG _{2b} , Human ads-PE/CY7	IgG	Polyclonal	20µl
S1090-19	Goat Anti-Mouse IgG _{2b} , Human ads-APC/CY7	IgG	Polyclonal	20µl
S1090-30	Goat Anti-Mouse IgG _{2b} , Human ads-AF488	IgG	Polyclonal	20µl
S1090-31	Goat Anti-Mouse IgG _{2b} , Human ads-AF647	IgG	Polyclonal	20µl
S1090-32	Goat Anti-Mouse IgG _{2b} , Human ads-AF555	IgG	Polyclonal	20µl
S1091-02	Goat Anti-Mouse IgG _{2b} -FITC	IgG	Polyclonal	20µl
S1091-04	Goat Anti-Mouse IgG _{2b} -AP	IgG	Polyclonal	20µl
S1091-05	Goat Anti-Mouse IgG _{2b} -HRP	IgG	Polyclonal	20µl
S1091-08	Goat Anti-Mouse IgG _{2b} -BIOT	IgG	Polyclonal	20µl
S1091-30	Goat Anti-Mouse IgG _{2b} -AF488	IgG	Polyclonal	20µl
S1091-31	Goat Anti-Mouse IgG _{2b} -AF647	IgG	Polyclonal	20µl
S1091-32	Goat Anti-Mouse IgG _{2b} -AF555	IgG	Polyclonal	20µl
S1092-02	Goat F(ab) ₂ Anti-Mouse IgG _{2b} , Human ads-FITC	IgG	Polyclonal	20µl
S1092-08	Goat F(ab) ₂ Anti-Mouse IgG _{2b} , Human ads-BIOT	IgG	Polyclonal	20µl
S1092-09	Goat F(ab) ₂ Anti-Mouse IgG _{2b} , Human ads-PE	IgG	Polyclonal	20µl
S1092-30	Goat F(ab) ₂ Anti-Mouse IgG _{2b} , Human ads-AF488	IgG	Polyclonal	20µl
S1092-31	Goat F(ab) ₂ Anti-Mouse IgG _{2b} , Human ads-AF647	IgG	Polyclonal	20µl
S1093-02	Goat Anti-Mouse IgG _{2b} , Human/Bovine/Horse SP ads-FITC	IgG	Polyclonal	20µl
S1093-04	Goat Anti-Mouse IgG _{2b} , Human/Bovine/Horse SP ads-AP	IgG	Polyclonal	20µl
S1093-05	Goat Anti-Mouse IgG _{2b} , Human/Bovine/Horse SP ads-HRP	IgG	Polyclonal	20µl
S1093-08	Goat Anti-Mouse IgG _{2b} , Human/Bovine/Horse SP ads-BIOT	IgG	Polyclonal	20µl
S1100-02	Goat Anti-Mouse IgG ₃ , Human ads-FITC	IgG	Polyclonal	20µl
S1100-03	Goat Anti-Mouse IgG ₃ , Human ads-TRITC	IgG	Polyclonal	20µl
S1100-04	Goat Anti-Mouse IgG ₃ , Human ads-AP	IgG	Polyclonal	20µl
S1100-05	Goat Anti-Mouse IgG ₃ , Human ads-HRP	IgG	Polyclonal	20µl
S1100-07	Goat Anti-Mouse IgG ₃ , Human ads-TXRD	IgG	Polyclonal	20µl
S1100-08	Goat Anti-Mouse IgG ₃ , Human ads-BIOT	IgG	Polyclonal	20µl
S1100-09	Goat Anti-Mouse IgG ₃ , Human ads-PE	IgG	Polyclonal	20µl
S1100-13	Goat Anti-Mouse IgG ₃ , Human ads-SPRD	IgG	Polyclonal	20µl
S1100-15	Goat Anti-Mouse IgG ₃ , Human ads-CY5	IgG	Polyclonal	20µl
S1100-17	Goat Anti-Mouse IgG ₃ , Human ads-PE/CY7	IgG	Polyclonal	20µl
S1100-19	Goat Anti-Mouse IgG ₃ , Human ads-APC/CY7	IgG	Polyclonal	20µl
S1100-30	Goat Anti-Mouse IgG ₃ , Human ads-AF488	IgG	Polyclonal	20µl
S1100-31	Goat Anti-Mouse IgG ₃ , Human ads-AF647	IgG	Polyclonal	20µl
S1100-32	Goat Anti-Mouse IgG ₃ , Human ads-AF555	IgG	Polyclonal	20µl
S1101-02	Goat Anti-Mouse IgG ₃ -FITC	IgG	Polyclonal	20µl
S1101-05	Goat Anti-Mouse IgG ₃ -HRP	IgG	Polyclonal	20µl
S1101-30	Goat Anti-Mouse IgG ₃ -AF488	IgG	Polyclonal	20µl
S1101-31	Goat Anti-Mouse IgG ₃ -AF647	IgG	Polyclonal	20µl
S1101-32	Goat Anti-Mouse IgG ₃ -AF555	IgG	Polyclonal	20µl
S1102-02	Goat F(ab) ₂ Anti-Mouse IgG ₃ , Human ads-FITC	IgG	Polyclonal	20µl
S1102-08	Goat F(ab) ₂ Anti-Mouse IgG ₃ , Human ads-BIOT	IgG	Polyclonal	20µl
S1102-09	Goat F(ab) ₂ Anti-Mouse IgG ₃ , Human ads-PE	IgG	Polyclonal	20µl
S1103-02	Goat Anti-Mouse IgG ₃ , Human/Bovine/Horse SP ads-FITC	IgG	Polyclonal	20µl
S1103-04	Goat Anti-Mouse IgG ₃ , Human/Bovine/Horse SP ads-AP	IgG	Polyclonal	20µl
S1103-05	Goat Anti-Mouse IgG ₃ , Human/Bovine/Horse SP ads-HRP	IgG	Polyclonal	20µl
S1103-08	Goat Anti-Mouse IgG ₃ , Human/Bovine/Horse SP ads-BIOT	IgG	Polyclonal	20µl

Secondary Antibody Promotion 2023
sample list

Cat.no.	Description	Isotype	Clone	Size
S1139-04	Rat Anti-Mouse IgM-AP	IgG ₁	SB73a	20µl
S1139-05	Rat Anti-Mouse IgM-HRP	IgG ₁	SB73a	20µl
S1140-02	Rat Anti-Mouse IgM-FITC	IgG _{2b}	1B4B1	20µl
S1140-04	Rat Anti-Mouse IgM-AP	IgG _{2b}	1B4B1	20µl
S1140-05	Rat Anti-Mouse IgM-HRP	IgG _{2b}	1B4B1	20µl
S1140-08	Rat Anti-Mouse IgM-BIOT	IgG _{2b}	1B4B1	20µl
S1140-09	Rat Anti-Mouse IgM-PE	IgG _{2b}	1B4B1	20µl
S1140-11	Rat Anti-Mouse IgM-APC	IgG _{2b}	1B4B1	20µl
S1140-17	Rat Anti-Mouse IgM-PE/CY7	IgG _{2b}	1B4B1	20µl
S1140-30	Rat Anti-Mouse IgM-AF488	IgG _{2b}	1B4B1	20µl
S1140-31	Rat Anti-Mouse IgM-AF647	IgG _{2b}	1B4B1	20µl
S1140-32	Rat Anti-Mouse IgM-AF555	IgG _{2b}	1B4B1	20µl
S1144-02	Rat Anti-Mouse IgG ₁ -FITC	IgG _{2b}	SB77e	20µl
S1144-04	Rat Anti-Mouse IgG ₁ -AP	IgG _{2b}	SB77e	20µl
S1144-05	Rat Anti-Mouse IgG ₁ -HRP	IgG _{2b}	SB77e	20µl
S1144-08	Rat Anti-Mouse IgG ₁ -BIOT	IgG _{2b}	SB77e	20µl
S1144-09	Rat Anti-Mouse IgG ₁ -PE	IgG _{2b}	SB77e	20µl
S1144-30	Rat Anti-Mouse IgG ₁ -AF488	IgG _{2b}	SB77e	20µl
S1144-31	Rat Anti-Mouse IgG ₁ -AF647	IgG _{2b}	SB77e	20µl
S1144-32	Rat Anti-Mouse IgG ₁ -AF555	IgG _{2b}	SB77e	20µl
S1155-02	Rat Anti-Mouse IgG _{2a} -FITC	IgG ₁	SB84a	20µl
S1155-04	Rat Anti-Mouse IgG _{2a} -AP	IgG ₁	SB84a	20µl
S1155-05	Rat Anti-Mouse IgG _{2a} -HRP	IgG ₁	SB84a	20µl
S1155-08	Rat Anti-Mouse IgG _{2a} -BIOT	IgG ₁	SB84a	20µl
S1155-09	Rat Anti-Mouse IgG _{2a} -PE	IgG ₁	SB84a	20µl
S1155-30	Rat Anti-Mouse IgG _{2a} -AF488	IgG ₁	SB84a	20µl
S1155-31	Rat Anti-Mouse IgG _{2a} -AF647	IgG ₁	SB84a	20µl
S1155-32	Rat Anti-Mouse IgG _{2a} -AF555	IgG ₁	SB84a	20µl
S1170-02	Rat Anti-Mouse Kappa-FITC	IgG ₁	187.1	20µl
S1170-04	Rat Anti-Mouse Kappa-AP	IgG ₁	187.1	20µl
S1170-05	Rat Anti-Mouse Kappa-HRP	IgG ₁	187.1	20µl
S1170-08	Rat Anti-Mouse Kappa-BIOT	IgG ₁	187.1	20µl
S1170-09	Rat Anti-Mouse Kappa-PE	IgG ₁	187.1	20µl
S1175-02	Rat Anti-Mouse Lambda-FITC	IgG _{2b}	JC5-1	20µl
S1175-04	Rat Anti-Mouse Lambda-AP	IgG _{2b}	JC5-1	20µl
S1175-05	Rat Anti-Mouse Lambda-HRP	IgG _{2b}	JC5-1	20µl
S1175-08	Rat Anti-Mouse Lambda-BIOT	IgG _{2b}	JC5-1	20µl
S1175-09	Rat Anti-Mouse Lambda-PE	IgG _{2b}	JC5-1	20µl
S1180-02	Rat Anti-Mouse Kappa-FITC	IgG ₁	H139-52.1	20µl
S1180-04	Rat Anti-Mouse Kappa-AP	IgG ₁	H139-52.1	20µl
S1180-05	Rat Anti-Mouse Kappa-HRP	IgG ₁	H139-52.1	20µl
S1180-08	Rat Anti-Mouse Kappa-BIOT	IgG ₁	H139-52.1	20µl
S1180-09	Rat Anti-Mouse Kappa-PE	IgG ₁	H139-52.1	20µl
S1186-02	Rat Anti-Mouse IgG _{2b} -FITC	IgG _{2b}	SB74g	20µl
S1186-04	Rat Anti-Mouse IgG _{2b} -AP	IgG _{2b}	SB74g	20µl
S1186-05	Rat Anti-Mouse IgG _{2b} -HRP	IgG _{2b}	SB74g	20µl
S1186-08	Rat Anti-Mouse IgG _{2b} -BIOT	IgG _{2b}	SB74g	20µl
S1186-09	Rat Anti-Mouse IgG _{2b} -PE	IgG _{2b}	SB74g	20µl
S1191-02	Rat Anti-Mouse IgG ₃ -FITC	IgG _{2a}	SB76b	20µl
S1191-04	Rat Anti-Mouse IgG ₃ -AP	IgG _{2a}	SB76b	20µl
S1191-05	Rat Anti-Mouse IgG ₃ -HRP	IgG _{2a}	SB76b	20µl
S1191-08	Rat Anti-Mouse IgG ₃ -BIOT	IgG _{2a}	SB76b	20µl
S1191-09	Rat Anti-Mouse IgG ₃ -PE	IgG _{2a}	SB76b	20µl
S1930-05	Mouse Anti-Hamster IgG ₃ -HRP	IgG ₁	SB139b	20µl
S1930-08	Mouse Anti-Hamster IgG ₃ -BIOT	IgG ₁	SB139b	20µl
S1935-02	Mouse Anti-Hamster IgG ₂ /IgG ₃ -FITC	IgG ₁	SB139e	20µl
S1935-05	Mouse Anti-Hamster IgG ₂ /IgG ₃ -HRP	IgG ₁	SB139e	20µl
S1935-08	Mouse Anti-Hamster IgG ₂ /IgG ₃ -BIOT	IgG ₁	SB139e	20µl
S1935-09	Mouse Anti-Hamster IgG ₂ /IgG ₃ -PE	IgG ₁	SB139e	20µl
S1940-02	Mouse Anti-Hamster IgG ₁ -FITC	IgG _{2b}	SB139a	20µl
S1940-05	Mouse Anti-Hamster IgG ₁ -HRP	IgG _{2b}	SB139a	20µl
S1940-08	Mouse Anti-Hamster IgG ₁ -BIOT	IgG _{2b}	SB139a	20µl
S1940-09	Mouse Anti-Hamster IgG ₁ -PE	IgG _{2b}	SB139a	20µl
S2010-02	Goat Anti-Human Ig-FITC	IgG	Polyclonal	20µl
S2010-03	Goat Anti-Human Ig-TRITC	IgG	Polyclonal	20µl
S2010-04	Goat Anti-Human Ig-AP	IgG	Polyclonal	20µl
S2010-05	Goat Anti-Human Ig-HRP	IgG	Polyclonal	20µl
S2010-06	Goat Anti-Human Ig-BGAL	IgG	Polyclonal	20µl
S2010-07	Goat Anti-Human Ig-TXRD	IgG	Polyclonal	20µl
S2010-08	Goat Anti-Human Ig-BIOT	IgG	Polyclonal	20µl
S2010-09	Goat Anti-Human Ig-PE	IgG	Polyclonal	20µl
S2010-30	Goat Anti-Human Ig-AF488	IgG	Polyclonal	20µl
S2010-31	Goat Anti-Human Ig-AF647	IgG	Polyclonal	20µl
S2010-32	Goat Anti-Human Ig-AF555	IgG	Polyclonal	20µl
S2012-02	Goat F(ab) ₂ Anti-Human Ig-FITC	IgG	Polyclonal	20µl

Secondary Antibody Promotion 2023
sample list

Cat.no.	Description	Isotype	Clone	Size
S2012-04	Goat F(ab) ₂ Anti-Human Ig-AP	IgG	Polyclonal	20µl
S2012-05	Goat F(ab) ₂ Anti-Human Ig-HRP	IgG	Polyclonal	20µl
S2012-08	Goat F(ab) ₂ Anti-Human Ig-BIOT	IgG	Polyclonal	20µl
S2012-09	Goat F(ab) ₂ Anti-Human Ig-PE	IgG	Polyclonal	20µl
S2012-30	Goat F(ab) ₂ Anti-Human Ig-AF488	IgG	Polyclonal	20µl
S2012-32	Goat F(ab) ₂ Anti-Human Ig-AF555	IgG	Polyclonal	20µl
S2014-02	Goat Anti-Human IgG Fc, Multi-Species SP ads-FITC	IgG	Polyclonal	20µl
S2014-04	Goat Anti-Human IgG Fc, Multi-Species SP ads-AP	IgG	Polyclonal	20µl
S2014-05	Goat Anti-Human IgG Fc, Multi-Species SP ads-HRP	IgG	Polyclonal	20µl
S2014-08	Goat Anti-Human IgG Fc, Multi-Species SP ads-BIOT	IgG	Polyclonal	20µl
S2014-09	Goat Anti-Human IgG Fc, Multi-Species SP ads-PE	IgG	Polyclonal	20µl
S2014-30	Goat Anti-Human IgG Fc, Multi-Species SP ads-AF488	IgG	Polyclonal	20µl
S2014-31	Goat Anti-Human IgG Fc, Multi-Species SP ads-AF647	IgG	Polyclonal	20µl
S2014-32	Goat Anti-Human IgG Fc, Multi-Species SP ads-AF555	IgG	Polyclonal	20µl
S2015-02	Goat Anti-Human IgG(H+L)-FITC	IgG	Polyclonal	20µl
S2015-04	Goat Anti-Human IgG(H+L)-AP	IgG	Polyclonal	20µl
S2015-05	Goat Anti-Human IgG(H+L)-HRP	IgG	Polyclonal	20µl
S2015-08	Goat Anti-Human IgG(H+L)-BIOT	IgG	Polyclonal	20µl
S2015-30	Goat Anti-Human IgG(H+L)-AF488	IgG	Polyclonal	20µl
S2015-31	Goat Anti-Human IgG(H+L)-AF647	IgG	Polyclonal	20µl
S2015-32	Goat Anti-Human IgG(H+L)-AF555	IgG	Polyclonal	20µl
S2016-02	Goat Anti-Human IgG(H+L), Mouse/Bovine/Horse SP ads-FITC	IgG	Polyclonal	20µl
S2016-04	Goat Anti-Human IgG(H+L), Mouse/Bovine/Horse SP ads-AP	IgG	Polyclonal	20µl
S2016-05	Goat Anti-Human IgG(H+L), Mouse/Bovine/Horse SP ads-HRP	IgG	Polyclonal	20µl
S2016-08	Goat Anti-Human IgG(H+L), Mouse/Bovine/Horse SP ads-BIOT	IgG	Polyclonal	20µl
S2016-31	Goat Anti-Human IgG(H+L), Mouse/Bovine/Horse SP ads-AF647	IgG	Polyclonal	20µl
S2016-32	Goat Anti-Human IgG(H+L), Mouse/Bovine/Horse SP ads-AF555	IgG	Polyclonal	20µl
S2020-02	Goat Anti-Human IgM-FITC	IgG	Polyclonal	20µl
S2020-03	Goat Anti-Human IgM-TRITC	IgG	Polyclonal	20µl
S2020-04	Goat Anti-Human IgM-AP	IgG	Polyclonal	20µl
S2020-05	Goat Anti-Human IgM-HRP	IgG	Polyclonal	20µl
S2020-07	Goat Anti-Human IgM-TXRD	IgG	Polyclonal	20µl
S2020-08	Goat Anti-Human IgM-BIOT	IgG	Polyclonal	20µl
S2020-09	Goat Anti-Human IgM-PE	IgG	Polyclonal	20µl
S2020-30	Goat Anti-Human IgM-AF488	IgG	Polyclonal	20µl
S2020-31	Goat Anti-Human IgM-AF647	IgG	Polyclonal	20µl
S2020-32	Goat Anti-Human IgM-AF555	IgG	Polyclonal	20µl
S2022-02	Goat F(ab) ₂ Anti-Human IgM-FITC	IgG	Polyclonal	20µl
S2022-04	Goat F(ab) ₂ Anti-Human IgM-AP	IgG	Polyclonal	20µl
S2022-05	Goat F(ab) ₂ Anti-Human IgM-HRP	IgG	Polyclonal	20µl
S2022-08	Goat F(ab) ₂ Anti-Human IgM-BIOT	IgG	Polyclonal	20µl
S2022-09	Goat F(ab) ₂ Anti-Human IgM-PE	IgG	Polyclonal	20µl
S2022-30	Goat F(ab) ₂ Anti-Human IgM-AF488	IgG	Polyclonal	20µl
S2022-31	Goat F(ab) ₂ Anti-Human IgM-AF647	IgG	Polyclonal	20µl
S2022-32	Goat F(ab) ₂ Anti-Human IgM-AF555	IgG	Polyclonal	20µl
S2023-02	Goat Anti-Human IgM, Mouse/Bovine/Horse SP ads-FITC	IgG	Polyclonal	20µl
S2023-04	Goat Anti-Human IgM, Mouse/Bovine/Horse SP ads-AP	IgG	Polyclonal	20µl
S2023-05	Goat Anti-Human IgM, Mouse/Bovine/Horse SP ads-HRP	IgG	Polyclonal	20µl
S2023-08	Goat Anti-Human IgM, Mouse/Bovine/Horse SP ads-BIOT	IgG	Polyclonal	20µl
S2023-31	Goat Anti-Human IgM, Mouse/Bovine/Horse SP ads-AF647	IgG	Polyclonal	20µl
S2040-02	Goat Anti-Human IgG-FITC	IgG	Polyclonal	20µl
S2040-03	Goat Anti-Human IgG-TRITC	IgG	Polyclonal	20µl
S2040-04	Goat Anti-Human IgG-AP	IgG	Polyclonal	20µl
S2040-05	Goat Anti-Human IgG-HRP	IgG	Polyclonal	20µl
S2040-06	Goat Anti-Human IgG-BGAL	IgG	Polyclonal	20µl
S2040-07	Goat Anti-Human IgG-TXRD	IgG	Polyclonal	20µl
S2040-08	Goat Anti-Human IgG-BIOT	IgG	Polyclonal	20µl
S2040-09	Goat Anti-Human IgG-PE	IgG	Polyclonal	20µl
S2040-30	Goat Anti-Human IgG-AF488	IgG	Polyclonal	20µl
S2040-31	Goat Anti-Human IgG-AF647	IgG	Polyclonal	20µl
S2040-32	Goat Anti-Human IgG-AF555	IgG	Polyclonal	20µl
S2041-02	Goat Fab Anti-Human IgG-FITC	IgG	Polyclonal	20µl
S2041-04	Goat Fab Anti-Human IgG-AP	IgG	Polyclonal	20µl
S2041-05	Goat Fab Anti-Human IgG-HRP	IgG	Polyclonal	20µl
S2041-08	Goat Fab Anti-Human IgG-BIOT	IgG	Polyclonal	20µl
S2042-02	Goat F(ab) ₂ Anti-Human IgG-FITC	IgG	Polyclonal	20µl
S2042-04	Goat F(ab) ₂ Anti-Human IgG-AP	IgG	Polyclonal	20µl
S2042-05	Goat F(ab) ₂ Anti-Human IgG-HRP	IgG	Polyclonal	20µl
S2042-08	Goat F(ab) ₂ Anti-Human IgG-BIOT	IgG	Polyclonal	20µl
S2042-09	Goat F(ab) ₂ Anti-Human IgG-PE	IgG	Polyclonal	20µl
S2042-30	Goat F(ab) ₂ Anti-Human IgG-AF488	IgG	Polyclonal	20µl
S2042-31	Goat F(ab) ₂ Anti-Human IgG-AF647	IgG	Polyclonal	20µl
S2042-32	Goat F(ab) ₂ Anti-Human IgG-AF555	IgG	Polyclonal	20µl
S2043-02	Goat F(ab) ₂ Anti-Human IgG, Mouse ads-FITC	IgG	Polyclonal	20µl
S2043-05	Goat F(ab) ₂ Anti-Human IgG, Mouse ads-HRP	IgG	Polyclonal	20µl

Secondary Antibody Promotion 2023
sample list

Cat.no.	Description	Isotype	Clone	Size
S2043-08	Goat F(ab) ₂ Anti-Human IgG, Mouse ads-BIOT	IgG	Polyclonal	20µl
S2043-09	Goat F(ab) ₂ Anti-Human IgG, Mouse ads-PE	IgG	Polyclonal	20µl
S2043-30	Goat F(ab) ₂ Anti-Human IgG, Mouse ads-AF488	IgG	Polyclonal	20µl
S2043-31	Goat F(ab) ₂ Anti-Human IgG, Mouse ads-AF647	IgG	Polyclonal	20µl
S2043-32	Goat F(ab) ₂ Anti-Human IgG, Mouse ads-AF555	IgG	Polyclonal	20µl
S2044-02	Goat Anti-Human IgG, Mouse ads-FITC	IgG	Polyclonal	20µl
S2044-04	Goat Anti-Human IgG, Mouse ads-AP	IgG	Polyclonal	20µl
S2044-05	Goat Anti-Human IgG, Mouse ads-HRP	IgG	Polyclonal	20µl
S2044-08	Goat Anti-Human IgG, Mouse ads-BIOT	IgG	Polyclonal	20µl
S2044-09	Goat Anti-Human IgG, Mouse ads-PE	IgG	Polyclonal	20µl
S2044-30	Goat Anti-Human IgG, Mouse ads-AF488	IgG	Polyclonal	20µl
S2044-31	Goat Anti-Human IgG, Mouse ads-AF647	IgG	Polyclonal	20µl
S2044-32	Goat Anti-Human IgG, Mouse ads-AF555	IgG	Polyclonal	20µl
S2045-02	Goat Anti-Human IgG, Mouse/Bovine/Horse SP ads-FITC	IgG	Polyclonal	20µl
S2045-04	Goat Anti-Human IgG, Mouse/Bovine/Horse SP ads-AP	IgG	Polyclonal	20µl
S2045-05	Goat Anti-Human IgG, Mouse/Bovine/Horse SP ads-HRP	IgG	Polyclonal	20µl
S2045-08	Goat Anti-Human IgG, Mouse/Bovine/Horse SP ads-BIOT	IgG	Polyclonal	20µl
S2045-30	Goat Anti-Human IgG, Mouse/Bovine/Horse SP ads-AF488	IgG	Polyclonal	20µl
S2045-31	Goat Anti-Human IgG, Mouse/Bovine/Horse SP ads-AF647	IgG	Polyclonal	20µl
S2045-32	Goat Anti-Human IgG, Mouse/Bovine/Horse SP ads-AF555	IgG	Polyclonal	20µl
S2046-02	Goat Anti-Human IgG Fd-FITC	IgG	Polyclonal	20µl
S2046-04	Goat Anti-Human IgG Fd-AP	IgG	Polyclonal	20µl
S2046-05	Goat Anti-Human IgG Fd-HRP	IgG	Polyclonal	20µl
S2046-08	Goat Anti-Human IgG Fd-BIOT	IgG	Polyclonal	20µl
S2046-31	Goat Anti-Human IgG Fd-AF647	IgG	Polyclonal	20µl
S2047-04	Goat Anti-Human Ig Fc-AP	IgG	Polyclonal	20µl
S2047-05	Goat Anti-Human Ig Fc-HRP	IgG	Polyclonal	20µl
S2047-08	Goat Anti-Human Ig Fc-BIOT	IgG	Polyclonal	20µl
S2047-30	Goat Anti-Human Ig Fc-AF488	IgG	Polyclonal	20µl
S2048-02	Goat Anti-Human IgG Fc-FITC	IgG	Polyclonal	20µl
S2048-04	Goat Anti-Human IgG Fc-AP	IgG	Polyclonal	20µl
S2048-05	Goat Anti-Human IgG Fc-HRP	IgG	Polyclonal	20µl
S2048-08	Goat Anti-Human IgG Fc-BIOT	IgG	Polyclonal	20µl
S2048-09	Goat Anti-Human IgG Fc-PE	IgG	Polyclonal	20µl
S2048-30	Goat Anti-Human IgG Fc-AF488	IgG	Polyclonal	20µl
S2048-31	Goat Anti-Human IgG Fc-AF647	IgG	Polyclonal	20µl
S2048-32	Goat Anti-Human IgG Fc-AF555	IgG	Polyclonal	20µl
S2049-02	Goat Anti-Human IgG, Monkey ads-FITC	IgG	Polyclonal	20µl
S2049-04	Goat Anti-Human IgG, Monkey ads-AP	IgG	Polyclonal	20µl
S2049-05	Goat Anti-Human IgG, Monkey ads-HRP	IgG	Polyclonal	20µl
S2049-08	Goat Anti-Human IgG, Monkey ads-BIOT	IgG	Polyclonal	20µl
S2049-30	Goat Anti-Human IgG, Monkey ads-AF488	IgG	Polyclonal	20µl
S2049-31	Goat Anti-Human IgG, Monkey ads-AF647	IgG	Polyclonal	20µl
S2049-32	Goat Anti-Human IgG, Monkey ads-AF555	IgG	Polyclonal	20µl
S2060-02	Goat Anti-Human Kappa-FITC	IgG	Polyclonal	20µl
S2060-03	Goat Anti-Human Kappa-TRITC	IgG	Polyclonal	20µl
S2060-04	Goat Anti-Human Kappa-AP	IgG	Polyclonal	20µl
S2060-05	Goat Anti-Human Kappa-HRP	IgG	Polyclonal	20µl
S2060-06	Goat Anti-Human Kappa-BGAL	IgG	Polyclonal	20µl
S2060-07	Goat Anti-Human Kappa-TXRD	IgG	Polyclonal	20µl
S2060-08	Goat Anti-Human Kappa-BIOT	IgG	Polyclonal	20µl
S2060-09	Goat Anti-Human Kappa-PE	IgG	Polyclonal	20µl
S2060-30	Goat Anti-Human Kappa-AF488	IgG	Polyclonal	20µl
S2060-31	Goat Anti-Human Kappa-AF647	IgG	Polyclonal	20µl
S2060-32	Goat Anti-Human Kappa-AF555	IgG	Polyclonal	20µl
S2061-02	Goat Anti-Human Kappa, Mouse ads-FITC	IgG	Polyclonal	20µl
S2061-03	Goat Anti-Human Kappa, Mouse ads-TRITC	IgG	Polyclonal	20µl
S2061-04	Goat Anti-Human Kappa, Mouse ads-AP	IgG	Polyclonal	20µl
S2061-05	Goat Anti-Human Kappa, Mouse ads-HRP	IgG	Polyclonal	20µl
S2061-07	Goat Anti-Human Kappa, Mouse ads-TXRD	IgG	Polyclonal	20µl
S2061-08	Goat Anti-Human Kappa, Mouse ads-BIOT	IgG	Polyclonal	20µl
S2061-09	Goat Anti-Human Kappa, Mouse ads-PE	IgG	Polyclonal	20µl
S2061-31	Goat Anti-Human Kappa, Mouse ads-AF647	IgG	Polyclonal	20µl
S2062-02	Goat F(ab) ₂ Anti-Human Kappa-FITC	IgG	Polyclonal	20µl
S2062-04	Goat F(ab) ₂ Anti-Human Kappa-AP	IgG	Polyclonal	20µl
S2062-05	Goat F(ab) ₂ Anti-Human Kappa-HRP	IgG	Polyclonal	20µl
S2062-08	Goat F(ab) ₂ Anti-Human Kappa-BIOT	IgG	Polyclonal	20µl
S2062-09	Goat F(ab) ₂ Anti-Human Kappa-PE	IgG	Polyclonal	20µl
S2062-30	Goat F(ab) ₂ Anti-Human Kappa-AF488	IgG	Polyclonal	20µl
S2062-31	Goat F(ab) ₂ Anti-Human Kappa-AF647	IgG	Polyclonal	20µl
S2062-32	Goat F(ab) ₂ Anti-Human Kappa-AF555	IgG	Polyclonal	20µl
S2063-02	Goat F(ab) ₂ Anti-Human Kappa, Mouse ads-FITC	IgG	Polyclonal	20µl
S2063-04	Goat F(ab) ₂ Anti-Human Kappa, Mouse ads-AP	IgG	Polyclonal	20µl
S2063-08	Goat F(ab) ₂ Anti-Human Kappa, Mouse ads-BIOT	IgG	Polyclonal	20µl
S2063-09	Goat F(ab) ₂ Anti-Human Kappa, Mouse ads-PE	IgG	Polyclonal	20µl

Secondary Antibody Promotion 2023
sample list

Cat.no.	Description	Isotype	Clone	Size
S2063-31	Goat F(ab) ₂ Anti-Human Kappa, Mouse ads-AF647	IgG	Polyclonal	20µl
S2064-02	Goat Anti-Human Kappa, Monkey ads-FITC	IgG	Polyclonal	20µl
S2064-04	Goat Anti-Human Kappa, Monkey ads-AP	IgG	Polyclonal	20µl
S2064-05	Goat Anti-Human Kappa, Monkey ads-HRP	IgG	Polyclonal	20µl
S2064-08	Goat Anti-Human Kappa, Monkey ads-BIOT	IgG	Polyclonal	20µl
S2064-30	Goat Anti-Human Kappa, Monkey ads-AF488	IgG	Polyclonal	20µl
S2064-31	Goat Anti-Human Kappa, Monkey ads-AF647	IgG	Polyclonal	20µl
S2070-02	Goat Anti-Human Lambda-FITC	IgG	Polyclonal	20µl
S2070-03	Goat Anti-Human Lambda-TRITC	IgG	Polyclonal	20µl
S2070-04	Goat Anti-Human Lambda-AP	IgG	Polyclonal	20µl
S2070-05	Goat Anti-Human Lambda-HRP	IgG	Polyclonal	20µl
S2070-07	Goat Anti-Human Lambda-TXRD	IgG	Polyclonal	20µl
S2070-08	Goat Anti-Human Lambda-BIOT	IgG	Polyclonal	20µl
S2070-09	Goat Anti-Human Lambda-PE	IgG	Polyclonal	20µl
S2070-30	Goat Anti-Human Lambda-AF488	IgG	Polyclonal	20µl
S2070-31	Goat Anti-Human Lambda-AF647	IgG	Polyclonal	20µl
S2070-32	Goat Anti-Human Lambda-AF555	IgG	Polyclonal	20µl
S2071-02	Goat Anti-Human Lambda, Mouse ads-FITC	IgG	Polyclonal	20µl
S2071-04	Goat Anti-Human Lambda, Mouse ads-AP	IgG	Polyclonal	20µl
S2071-05	Goat Anti-Human Lambda, Mouse ads-HRP	IgG	Polyclonal	20µl
S2071-08	Goat Anti-Human Lambda, Mouse ads-BIOT	IgG	Polyclonal	20µl
S2071-09	Goat Anti-Human Lambda, Mouse ads-PE	IgG	Polyclonal	20µl
S2071-31	Goat Anti-Human Lambda, Mouse ads-AF647	IgG	Polyclonal	20µl
S2072-02	Goat F(ab) ₂ Anti-Human Lambda-FITC	IgG	Polyclonal	20µl
S2072-04	Goat F(ab) ₂ Anti-Human Lambda-AP	IgG	Polyclonal	20µl
S2072-05	Goat F(ab) ₂ Anti-Human Lambda-HRP	IgG	Polyclonal	20µl
S2072-08	Goat F(ab) ₂ Anti-Human Lambda-BIOT	IgG	Polyclonal	20µl
S2072-09	Goat F(ab) ₂ Anti-Human Lambda-PE	IgG	Polyclonal	20µl
S2072-31	Goat F(ab) ₂ Anti-Human Lambda-AF647	IgG	Polyclonal	20µl
S2073-02	Goat F(ab) ₂ Anti-Human Lambda, Mouse ads-FITC	IgG	Polyclonal	20µl
S2073-04	Goat F(ab) ₂ Anti-Human Lambda, Mouse ads-AP	IgG	Polyclonal	20µl
S2073-05	Goat F(ab) ₂ Anti-Human Lambda, Mouse ads-HRP	IgG	Polyclonal	20µl
S2073-08	Goat F(ab) ₂ Anti-Human Lambda, Mouse ads-BIOT	IgG	Polyclonal	20µl
S2073-09	Goat F(ab) ₂ Anti-Human Lambda, Mouse ads-PE	IgG	Polyclonal	20µl
S2073-30	Goat F(ab) ₂ Anti-Human Lambda, Mouse ads-AF488	IgG	Polyclonal	20µl
S2081-04	Goat Anti-Human IgG Fc, Mouse/Bovine/Horse SP ads-AP	IgG	Polyclonal	20µl
S2081-05	Goat Anti-Human IgG Fc, Mouse/Bovine/Horse SP ads-HRP	IgG	Polyclonal	20µl
S2081-08	Goat Anti-Human IgG Fc, Mouse/Bovine/Horse SP ads-BIOT	IgG	Polyclonal	20µl
S2081-30	Goat Anti-Human IgG Fc, Mouse/Bovine/Horse SP ads-AF488	IgG	Polyclonal	20µl
S2081-31	Goat Anti-Human IgG Fc, Mouse/Bovine/Horse SP ads-AF647	IgG	Polyclonal	20µl
S2085-02	Goat Anti-Human Ig Fab-FITC	IgG	Polyclonal	20µl
S2085-04	Goat Anti-Human Ig Fab-AP	IgG	Polyclonal	20µl
S2085-05	Goat Anti-Human Ig Fab-HRP	IgG	Polyclonal	20µl
S2085-08	Goat Anti-Human Ig Fab-BIOT	IgG	Polyclonal	20µl
S2087-02	Goat Anti-Human IgG(H+L), Multi-Species SP ads-FITC	IgG	Polyclonal	20µl
S2087-04	Goat Anti-Human IgG(H+L), Multi-Species SP ads-AP	IgG	Polyclonal	20µl
S2087-05	Goat Anti-Human IgG(H+L), Multi-Species SP ads-HRP	IgG	Polyclonal	20µl
S2087-08	Goat Anti-Human IgG(H+L), Multi-Species SP ads-BIOT	IgG	Polyclonal	20µl
S3000-02	Goat Fab Anti-Rat Ig, Mouse ads-FITC	IgG	Polyclonal	20µl
S3000-05	Goat Fab Anti-Rat Ig, Mouse ads-HRP	IgG	Polyclonal	20µl
S3000-08	Goat Fab Anti-Rat Ig, Mouse ads-BIOT	IgG	Polyclonal	20µl
S3010-02	Goat Anti-Rat Ig, Mouse ads-FITC	IgG	Polyclonal	20µl
S3010-03	Goat Anti-Rat Ig, Mouse ads-TRITC	IgG	Polyclonal	20µl
S3010-04	Goat Anti-Rat Ig, Mouse ads-AP	IgG	Polyclonal	20µl
S3010-05	Goat Anti-Rat Ig, Mouse ads-HRP	IgG	Polyclonal	20µl
S3010-07	Goat Anti-Rat Ig, Mouse ads-TXRD	IgG	Polyclonal	20µl
S3010-08	Goat Anti-Rat Ig, Mouse ads-BIOT	IgG	Polyclonal	20µl
S3010-09	Goat Anti-Rat Ig, Mouse ads-PE	IgG	Polyclonal	20µl
S3010-11	Goat Anti-Rat Ig, Mouse ads-APC	IgG	Polyclonal	20µl
S3010-31	Goat Anti-Rat Ig, Mouse ads-AF647	IgG	Polyclonal	20µl
S3010-32	Goat Anti-Rat Ig, Mouse ads-AF555	IgG	Polyclonal	20µl
S3020-02	Goat Anti-Rat IgM-FITC	IgG	Polyclonal	20µl
S3020-03	Goat Anti-Rat IgM-TRITC	IgG	Polyclonal	20µl
S3020-04	Goat Anti-Rat IgM-AP	IgG	Polyclonal	20µl
S3020-05	Goat Anti-Rat IgM-HRP	IgG	Polyclonal	20µl
S3020-07	Goat Anti-Rat IgM-TXRD	IgG	Polyclonal	20µl
S3020-08	Goat Anti-Rat IgM-BIOT	IgG	Polyclonal	20µl
S3020-09	Goat Anti-Rat IgM-PE	IgG	Polyclonal	20µl
S3020-30	Goat Anti-Rat IgM-AF488	IgG	Polyclonal	20µl
S3020-31	Goat Anti-Rat IgM-AF647	IgG	Polyclonal	20µl
S3020-32	Goat Anti-Rat IgM-AF555	IgG	Polyclonal	20µl
S3030-02	Goat Anti-Rat IgG-FITC	IgG	Polyclonal	20µl
S3030-03	Goat Anti-Rat IgG-TRITC	IgG	Polyclonal	20µl
S3030-04	Goat Anti-Rat IgG-AP	IgG	Polyclonal	20µl

Secondary Antibody Promotion 2023
sample list

Cat.no.	Description	Isotype	Clone	Size
S3030-05	Goat Anti-Rat IgG-HRP	IgG	Polyclonal	20µl
S3030-06	Goat Anti-Rat IgG-BGAL	IgG	Polyclonal	20µl
S3030-08	Goat Anti-Rat IgG-BIOT	IgG	Polyclonal	20µl
S3030-09	Goat Anti-Rat IgG-PE	IgG	Polyclonal	20µl
S3030-30	Goat Anti-Rat IgG-AF488	IgG	Polyclonal	20µl
S3030-31	Goat Anti-Rat IgG-AF647	IgG	Polyclonal	20µl
S3030-32	Goat Anti-Rat IgG-AF555	IgG	Polyclonal	20µl
S3032-05	Goat F(ab) ₂ Anti-Rat IgG-HRP	IgG	Polyclonal	20µl
S3050-02	Goat Anti-Rat IgG(H+L), Mouse ads-FITC	IgG	Polyclonal	20µl
S3050-03	Goat Anti-Rat IgG(H+L), Mouse ads-TRITC	IgG	Polyclonal	20µl
S3050-04	Goat Anti-Rat IgG(H+L), Mouse ads-AP	IgG	Polyclonal	20µl
S3050-05	Goat Anti-Rat IgG(H+L), Mouse ads-HRP	IgG	Polyclonal	20µl
S3050-06	Goat Anti-Rat IgG(H+L), Mouse ads-BGAL	IgG	Polyclonal	20µl
S3050-07	Goat Anti-Rat IgG(H+L), Mouse ads-TXRD	IgG	Polyclonal	20µl
S3050-08	Goat Anti-Rat IgG(H+L), Mouse ads-BIOT	IgG	Polyclonal	20µl
S3050-09	Goat Anti-Rat IgG(H+L), Mouse ads-PE	IgG	Polyclonal	20µl
S3050-09	Goat Anti-Rat IgG(H+L), Mouse ads-PE	IgG	Polyclonal	20µl
S3050-30	Goat Anti-Rat IgG(H+L), Mouse ads-AF488	IgG	Polyclonal	20µl
S3050-31	Goat Anti-Rat IgG(H+L), Mouse ads-AF647	IgG	Polyclonal	20µl
S3050-32	Goat Anti-Rat IgG(H+L), Mouse ads-AF555	IgG	Polyclonal	20µl
S3051-02	Goat Anti-Rat IgG(H+L), Mouse/Human ads-FITC	IgG	Polyclonal	20µl
S3051-05	Goat Anti-Rat IgG(H+L), Mouse/Human ads-HRP	IgG	Polyclonal	20µl
S3051-08	Goat Anti-Rat IgG(H+L), Mouse/Human ads-BIOT	IgG	Polyclonal	20µl
S3052-02	Goat F(ab) ₂ Anti-Rat IgG(H+L), Mouse ads-FITC	IgG	Polyclonal	20µl
S3052-07	Goat F(ab) ₂ Anti-Rat IgG(H+L), Mouse ads-TXRD	IgG	Polyclonal	20µl
S3052-08	Goat F(ab) ₂ Anti-Rat IgG(H+L), Mouse ads-BIOT	IgG	Polyclonal	20µl
S3052-09	Goat F(ab) ₂ Anti-Rat IgG(H+L), Mouse ads-PE	IgG	Polyclonal	20µl
S3052-10	Goat F(ab) ₂ Anti-Rat IgG(H+L), Mouse ads-PE/TXRD	IgG	Polyclonal	20µl
S3052-13	Goat F(ab) ₂ Anti-Rat IgG(H+L), Mouse ads-SPRD	IgG	Polyclonal	20µl
S3052-30	Goat F(ab) ₂ Anti-Rat IgG(H+L), Mouse ads-AF488	IgG	Polyclonal	20µl
S3052-31	Goat F(ab) ₂ Anti-Rat IgG(H+L), Mouse ads-AF647	IgG	Polyclonal	20µl
S3053-02	Mouse Anti-Rat IgG(H+L), Mouse SP ads-FITC	IgG	Polyclonal	20µl
S3053-08	Mouse Anti-Rat IgG(H+L), Mouse SP ads-BIOT	IgG	Polyclonal	20µl
S3055-31	Goat Anti-Rat IgG(H+L), Multi-Species SP ads-AF647	IgG	Polyclonal	20µl
S3060-02	Mouse Anti-Rat IgG ₁ -FITC	IgG ₁	G17E7	20µl
S3060-05	Mouse Anti-Rat IgG ₁ -HRP	IgG ₁	G17E7	20µl
S3060-08	Mouse Anti-Rat IgG ₁ -BIOT	IgG ₁	G17E7	20µl
S3060-09	Mouse Anti-Rat IgG ₁ -PE	IgG ₁	G17E7	20µl
S3060-15	Mouse Anti-Rat IgG ₁ -CY5	IgG ₁	G17E7	20µl
S3061-02	Mouse Anti-Rat IgG ₁ -FITC	IgG _{2b}	G11C5	20µl
S3061-04	Mouse Anti-Rat IgG ₁ -AP	IgG _{2b}	G11C5	20µl
S3061-05	Mouse Anti-Rat IgG ₁ -HRP	IgG _{2b}	G11C5	20µl
S3061-08	Mouse Anti-Rat IgG ₁ -BIOT	IgG _{2b}	G11C5	20µl
S3061-09	Mouse Anti-Rat IgG ₁ -PE	IgG _{2b}	G11C5	20µl
S3061-15	Mouse Anti-Rat IgG ₁ -CY5	IgG _{2b}	G11C5	20µl
S3061-30	Mouse Anti-Rat IgG ₁ -AF488	IgG _{2b}	G11C5	20µl
S3061-31	Mouse Anti-Rat IgG ₁ -AF647	IgG _{2b}	G11C5	20µl
S3061-32	Mouse Anti-Rat IgG ₁ -AF555	IgG _{2b}	G11C5	20µl
S3065-02	Mouse Anti-Rat IgG _{2a} -FITC	IgG ₁	2A8F4	20µl
S3065-04	Mouse Anti-Rat IgG _{2a} -AP	IgG ₁	2A8F4	20µl
S3065-05	Mouse Anti-Rat IgG _{2a} -HRP	IgG ₁	2A8F4	20µl
S3065-08	Mouse Anti-Rat IgG _{2a} -BIOT	IgG ₁	2A8F4	20µl
S3065-09	Mouse Anti-Rat IgG _{2a} -PE	IgG ₁	2A8F4	20µl
S3065-15	Mouse Anti-Rat IgG _{2a} -CY5	IgG ₁	2A8F4	20µl
S3065-30	Mouse Anti-Rat IgG _{2a} -AF488	IgG ₁	2A8F4	20µl
S3065-31	Mouse Anti-Rat IgG _{2a} -AF647	IgG ₁	2A8F4	20µl
S3065-32	Mouse Anti-Rat IgG _{2a} -AF555	IgG ₁	2A8F4	20µl
S3070-02	Mouse Anti-Rat IgG _{2b} -FITC	IgG _{2b}	2B10A8	20µl
S3070-04	Mouse Anti-Rat IgG _{2b} -AP	IgG _{2b}	2B10A8	20µl
S3070-05	Mouse Anti-Rat IgG _{2b} -HRP	IgG _{2b}	2B10A8	20µl
S3070-09	Mouse Anti-Rat IgG _{2b} -PE	IgG _{2b}	2B10A8	20µl
S3070-15	Mouse Anti-Rat IgG _{2b} -CY5	IgG _{2b}	2B10A8	20µl
S3070-28	Mouse Anti-Rat IgG _{2b} -BIMA	IgG _{2b}	2B10A8	20µl
S3070-30	Mouse Anti-Rat IgG _{2b} -AF488	IgG _{2b}	2B10A8	20µl
S3070-31	Mouse Anti-Rat IgG _{2b} -AF647	IgG _{2b}	2B10A8	20µl
S3070-32	Mouse Anti-Rat IgG _{2b} -AF555	IgG _{2b}	2B10A8	20µl
S3075-02	Mouse Anti-Rat IgG _{2c} -FITC	IgG _{2a}	2C8F1	20µl
S3075-04	Mouse Anti-Rat IgG _{2c} -AP	IgG _{2a}	2C8F1	20µl
S3075-05	Mouse Anti-Rat IgG _{2c} -HRP	IgG _{2a}	2C8F1	20µl
S3075-08	Mouse Anti-Rat IgG _{2c} -BIOT	IgG _{2a}	2C8F1	20µl
S3075-09	Mouse Anti-Rat IgG _{2c} -PE	IgG _{2a}	2C8F1	20µl
S3075-15	Mouse Anti-Rat IgG _{2c} -CY5	IgG _{2a}	2C8F1	20µl
S3080-02	Mouse Anti-Rat IgM-FITC	IgG ₁	M2A1	20µl
S3080-05	Mouse Anti-Rat IgM-HRP	IgG ₁	M2A1	20µl
S3080-09	Mouse Anti-Rat IgM-PE	IgG ₁	M2A1	20µl

Secondary Antibody Promotion 2023
sample list

Cat.no.	Description	Isotype	Clone	Size
S3080-15	Mouse Anti-Rat IgM-CY5	IgG ₁	M2A1	20µl
S3080-30	Mouse Anti-Rat IgM-AF488	IgG ₁	M2A1	20µl
S3080-31	Mouse Anti-Rat IgM-AF647	IgG ₁	M2A1	20µl
S3080-32	Mouse Anti-Rat IgM-AF555	IgG ₁	M2A1	20µl
S3090-02	Mouse Anti-Rat Kappa-FITC	IgG _{2a}	K4F5	20µl
S3090-04	Mouse Anti-Rat Kappa-AP	IgG _{2a}	K4F5	20µl
S3090-05	Mouse Anti-Rat Kappa-HRP	IgG _{2a}	K4F5	20µl
S3090-08	Mouse Anti-Rat Kappa-BIOT	IgG _{2a}	K4F5	20µl
S3090-09	Mouse Anti-Rat Kappa-PE	IgG _{2a}	K4F5	20µl
S4010-02	Goat Anti-Rabbit Ig, Human ads-FITC	IgG	Polyclonal	20µl
S4010-03	Goat Anti-Rabbit Ig, Human ads-TRITC	IgG	Polyclonal	20µl
S4010-04	Goat Anti-Rabbit Ig, Human ads-AP	IgG	Polyclonal	20µl
S4010-05	Goat Anti-Rabbit Ig, Human ads-HRP	IgG	Polyclonal	20µl
S4010-06	Goat Anti-Rabbit Ig, Human ads-BGAL	IgG	Polyclonal	20µl
S4010-07	Goat Anti-Rabbit Ig, Human ads-TXRD	IgG	Polyclonal	20µl
S4010-08	Goat Anti-Rabbit Ig, Human ads-BIOT	IgG	Polyclonal	20µl
S4010-09	Goat Anti-Rabbit Ig, Human ads-PE	IgG	Polyclonal	20µl
S4010-09	Goat Anti-Rabbit Ig, Human ads-PE	IgG	Polyclonal	20µl
S4010-30	Goat Anti-Rabbit Ig, Human ads-AF488	IgG	Polyclonal	20µl
S4010-31	Goat Anti-Rabbit Ig, Human ads-AF647	IgG	Polyclonal	20µl
S4010-32	Goat Anti-Rabbit Ig, Human ads-AF555	IgG	Polyclonal	20µl
S4020-02	Goat Anti-Rabbit IgM-FITC	IgG	Polyclonal	20µl
S4020-03	Goat Anti-Rabbit IgM-TRITC	IgG	Polyclonal	20µl
S4020-04	Goat Anti-Rabbit IgM-AP	IgG	Polyclonal	20µl
S4020-05	Goat Anti-Rabbit IgM-HRP	IgG	Polyclonal	20µl
S4020-08	Goat Anti-Rabbit IgM-BIOT	IgG	Polyclonal	20µl
S4020-09	Goat Anti-Rabbit IgM-PE	IgG	Polyclonal	20µl
S4020-31	Goat Anti-Rabbit IgM-AF647	IgG	Polyclonal	20µl
S4020-32	Goat Anti-Rabbit IgM-AF555	IgG	Polyclonal	20µl
S4021-05	Goat Anti-Rabbit IgM, Multi-Species SP ads-HRP	IgG	Polyclonal	20µl
S4030-02	Goat Anti-Rabbit IgG-FITC	IgG	Polyclonal	20µl
S4030-03	Goat Anti-Rabbit IgG-TRITC	IgG	Polyclonal	20µl
S4030-04	Goat Anti-Rabbit IgG-AP	IgG	Polyclonal	20µl
S4030-05	Goat Anti-Rabbit IgG-HRP	IgG	Polyclonal	20µl
S4030-08	Goat Anti-Rabbit IgG-BIOT	IgG	Polyclonal	20µl
S4030-09	Goat Anti-Rabbit IgG-PE	IgG	Polyclonal	20µl
S4030-30	Goat Anti-Rabbit IgG-AF488	IgG	Polyclonal	20µl
S4030-31	Goat Anti-Rabbit IgG-AF647	IgG	Polyclonal	20µl
S4030-32	Goat Anti-Rabbit IgG-AF555	IgG	Polyclonal	20µl
S4040-02	Goat Anti-Rabbit IgG Fab-FITC	IgG	Polyclonal	20µl
S4040-05	Goat Anti-Rabbit IgG Fab-HRP	IgG	Polyclonal	20µl
S4040-08	Goat Anti-Rabbit IgG Fab-BIOT	IgG	Polyclonal	20µl
S4041-02	Goat Anti-Rabbit IgG Fc-FITC	IgG	Polyclonal	20µl
S4041-05	Goat Anti-Rabbit IgG Fc-HRP	IgG	Polyclonal	20µl
S4041-08	Goat Anti-Rabbit IgG Fc-BIOT	IgG	Polyclonal	20µl
S4049-04	Goat Anti-Rabbit IgG(H+L), Mouse/Rat/Human ads-AP	IgG	Polyclonal	20µl
S4049-05	Goat Anti-Rabbit IgG(H+L), Mouse/Rat/Human ads-HRP	IgG	Polyclonal	20µl
S4049-32	Goat Anti-Rabbit IgG(H+L), Mouse/Rat/Human ads-AF555	IgG	Polyclonal	20µl
S4050-02	Goat Anti-Rabbit IgG(H+L), Mouse/Human ads-FITC	IgG	Polyclonal	20µl
S4050-03	Goat Anti-Rabbit IgG(H+L), Mouse/Human ads-TRITC	IgG	Polyclonal	20µl
S4050-04	Goat Anti-Rabbit IgG(H+L), Mouse/Human ads-AP	IgG	Polyclonal	20µl
S4050-05	Goat Anti-Rabbit IgG(H+L), Mouse/Human ads-HRP	IgG	Polyclonal	20µl
S4050-07	Goat Anti-Rabbit IgG(H+L), Mouse/Human ads-TXRD	IgG	Polyclonal	20µl
S4050-08	Goat Anti-Rabbit IgG(H+L), Mouse/Human ads-BIOT	IgG	Polyclonal	20µl
S4050-09	Goat Anti-Rabbit IgG(H+L), Mouse/Human ads-PE	IgG	Polyclonal	20µl
S4050-11	Goat Anti-Rabbit IgG(H+L), Mouse/Human ads-APC	IgG	Polyclonal	20µl
S4050-13	Goat Anti-Rabbit IgG(H+L), Mouse/Human ads-SPRD	IgG	Polyclonal	20µl
S4050-30	Goat Anti-Rabbit IgG(H+L), Mouse/Human ads-AF488	IgG	Polyclonal	20µl
S4050-31	Goat Anti-Rabbit IgG(H+L), Mouse/Human ads-AF647	IgG	Polyclonal	20µl
S4050-32	Goat Anti-Rabbit IgG(H+L), Mouse/Human ads-AF555	IgG	Polyclonal	20µl
S4051-02	Goat Anti-Rabbit IgG(H+L), Human ads-FITC	IgG	Polyclonal	20µl
S4051-04	Goat Anti-Rabbit IgG(H+L), Human ads-AP	IgG	Polyclonal	20µl
S4051-05	Goat Anti-Rabbit IgG(H+L), Human ads-HRP	IgG	Polyclonal	20µl
S4051-08	Goat Anti-Rabbit IgG(H+L), Human ads-BIOT	IgG	Polyclonal	20µl
S4051-30	Goat Anti-Rabbit IgG(H+L), Human ads-AF488	IgG	Polyclonal	20µl
S4051-32	Goat Anti-Rabbit IgG(H+L), Human ads-AF555	IgG	Polyclonal	20µl
S4052-02	Goat F(ab) ₂ Anti-Rabbit IgG(H+L), Mouse/Human ads-FITC	IgG	Polyclonal	20µl
S4052-04	Goat F(ab) ₂ Anti-Rabbit IgG(H+L), Mouse/Human ads-AP	IgG	Polyclonal	20µl
S4052-05	Goat F(ab) ₂ Anti-Rabbit IgG(H+L), Mouse/Human ads-HRP	IgG	Polyclonal	20µl
S4052-08	Goat F(ab) ₂ Anti-Rabbit IgG(H+L), Mouse/Human ads-BIOT	IgG	Polyclonal	20µl
S4052-09	Goat F(ab) ₂ Anti-Rabbit IgG(H+L), Mouse/Human ads-PE	IgG	Polyclonal	20µl
S4052-31	Goat F(ab) ₂ Anti-Rabbit IgG(H+L), Mouse/Human ads-AF647	IgG	Polyclonal	20µl
S4052-32	Goat F(ab) ₂ Anti-Rabbit IgG(H+L), Mouse/Human ads-AF555	IgG	Polyclonal	20µl
S4054-02	Goat F(ab) ₂ Anti-Rabbit IgG(H+L), Mouse/Rat/Human ads-FITC	IgG	Polyclonal	20µl

Secondary Antibody Promotion 2023
sample list

Cat.no.	Description	Isotype	Clone	Size
S4054-05	Goat F(ab) ₂ Anti-Rabbit IgG(H+L), Mouse/Rat/Human ads-HRP	IgG	Polyclonal	20µl
S4054-08	Goat F(ab) ₂ Anti-Rabbit IgG(H+L), Mouse/Rat/Human ads-BIOT	IgG	Polyclonal	20µl
S4054-30	Goat F(ab) ₂ Anti-Rabbit IgG(H+L), Mouse/Rat/Human ads-AF488	IgG	Polyclonal	20µl
S4054-32	Goat F(ab) ₂ Anti-Rabbit IgG(H+L), Mouse/Rat/Human ads-AF555	IgG	Polyclonal	20µl
S4055-02	Goat Anti-Rabbit IgG(H+L)-FITC	IgG	Polyclonal	20µl
S4055-04	Goat Anti-Rabbit IgG(H+L)-AP	IgG	Polyclonal	20µl
S4055-05	Goat Anti-Rabbit IgG(H+L)-HRP	IgG	Polyclonal	20µl
S4055-08	Goat Anti-Rabbit IgG(H+L)-BIOT	IgG	Polyclonal	20µl
S4055-30	Goat Anti-Rabbit IgG(H+L)-AF488	IgG	Polyclonal	20µl
S4055-31	Goat Anti-Rabbit IgG(H+L)-AF647	IgG	Polyclonal	20µl
S4055-32	Goat Anti-Rabbit IgG(H+L)-AF555	IgG	Polyclonal	20µl
S4058-02	Goat Anti-Rabbit IgG Fab, Mouse/Human SP ads-FITC	IgG	Polyclonal	20µl
S4058-04	Goat Anti-Rabbit IgG Fab, Mouse/Human SP ads-AP	IgG	Polyclonal	20µl
S4058-05	Goat Anti-Rabbit IgG Fab, Mouse/Human SP ads-HRP	IgG	Polyclonal	20µl
S4058-08	Goat Anti-Rabbit IgG Fab, Mouse/Human SP ads-BIOT	IgG	Polyclonal	20µl
S4060-04	Mouse Anti-Rabbit Light Chain-AP	IgG _{2b}	SB62a	20µl
S4060-05	Mouse Anti-Rabbit Light Chain-HRP	IgG _{2b}	SB62a	20µl
S4065-05	Rat Anti-Rabbit Ig-HRP	IgG _{2a}	SB87a	20µl
S4065-08	Rat Anti-Rabbit Ig-BIOT	IgG _{2a}	SB87a	20µl
S4065-09	Rat Anti-Rabbit Ig-PE	IgG _{2a}	SB87a	20µl
S4065-13	Rat Anti-Rabbit Ig-SPRD	IgG _{2a}	SB87a	20µl
S4090-02	Mouse Anti-Rabbit IgG-FITC	IgG ₁	2A9	20µl
S4090-04	Mouse Anti-Rabbit IgG-AP	IgG ₁	2A9	20µl
S4090-05	Mouse Anti-Rabbit IgG-HRP	IgG ₁	2A9	20µl
S4090-08	Mouse Anti-Rabbit IgG-BIOT	IgG ₁	2A9	20µl
S4090-09	Mouse Anti-Rabbit IgG-PE	IgG ₁	2A9	20µl
S4090-13	Mouse Anti-Rabbit IgG-SPRD	IgG ₁	2A9	20µl
S4090-32	Mouse Anti-Rabbit IgG-AF555	IgG ₁	2A9	20µl
S4091-02	Mouse Anti-Rabbit IgG(H+L), Mouse/Human/Goat/Sheep SP ads-FITC	IgG	Polyclonal	20µl
S4091-08	Mouse Anti-Rabbit IgG(H+L), Mouse/Human/Goat/Sheep SP ads-BIOT	IgG	Polyclonal	20µl
S4700-02	Mouse Anti-Monkey IgG-FITC	IgG ₁	SB108a	20µl
S4700-04	Mouse Anti-Monkey IgG-AP	IgG ₁	SB108a	20µl
S4700-05	Mouse Anti-Monkey IgG-HRP	IgG ₁	SB108a	20µl
S4700-08	Mouse Anti-Monkey IgG-BIOT	IgG ₁	SB108a	20µl
S4700-09	Mouse Anti-Monkey IgG-PE	IgG ₁	SB108a	20µl
S4700-30	Mouse Anti-Monkey IgG-AF488	IgG ₁	SB108a	20µl
S4700-31	Mouse Anti-Monkey IgG-AF647	IgG ₁	SB108a	20µl
S4700-32	Mouse Anti-Monkey IgG-AF555	IgG ₁	SB108a	20µl
S5820-04	Mouse Anti-Llama IgM-AP	IgG _{2b}	SB189b	20µl
S5820-05	Mouse Anti-Llama IgM-HRP	IgG _{2b}	SB189b	20µl
S5820-08	Mouse Anti-Llama IgM-BIOT	IgG _{2b}	SB189b	20µl
S5870-04	Mouse Anti-Llama IgG ₁ -AP	IgG ₁	SB182L	20µl
S5870-05	Mouse Anti-Llama IgG ₁ -HRP	IgG ₁	SB182L	20µl
S5870-08	Mouse Anti-Llama IgG ₁ -BIOT	IgG ₁	SB182L	20µl
S5880-04	Mouse Anti-Llama IgG ₂ /IgG ₃ -AP	IgG _{2b}	SB183c	20µl
S5880-05	Mouse Anti-Llama IgG ₂ /IgG ₃ -HRP	IgG _{2b}	SB183c	20µl
S5880-08	Mouse Anti-Llama IgG ₂ /IgG ₃ -BIOT	IgG _{2b}	SB183c	20µl
S5880-09	Mouse Anti-Llama IgG ₂ /IgG ₃ -PE	IgG _{2b}	SB183c	20µl
S5880-11	Mouse Anti-Llama IgG ₂ /IgG ₃ -APC	IgG _{2b}	SB183c	20µl
S5880-31	Mouse Anti-Llama IgG ₂ /IgG ₃ -AF647	IgG _{2b}	SB183c	20µl
S6000-02	Rabbit F(ab) ₂ Anti-Human IgG(H+L)-FITC	IgG	Polyclonal	20µl
S6000-04	Rabbit F(ab) ₂ Anti-Human IgG(H+L)-AP	IgG	Polyclonal	20µl
S6000-05	Rabbit F(ab) ₂ Anti-Human IgG(H+L)-HRP	IgG	Polyclonal	20µl
S6000-30	Rabbit F(ab) ₂ Anti-Human IgG(H+L)-AF488	IgG	Polyclonal	20µl
S6005-02	Rabbit F(ab) ₂ Anti-Human IgG(H+L), Mouse ads-FITC	IgG	Polyclonal	20µl
S6005-04	Rabbit F(ab) ₂ Anti-Human IgG(H+L), Mouse ads-AP	IgG	Polyclonal	20µl
S6005-05	Rabbit F(ab) ₂ Anti-Human IgG(H+L), Mouse ads-HRP	IgG	Polyclonal	20µl
S6010-02	Rabbit F(ab) ₂ Anti-Sheep IgG(H+L)-FITC	IgG	Polyclonal	20µl
S6010-04	Rabbit F(ab) ₂ Anti-Sheep IgG(H+L)-AP	IgG	Polyclonal	20µl
S6010-05	Rabbit F(ab) ₂ Anti-Sheep IgG(H+L)-HRP	IgG	Polyclonal	20µl
S6016-02	Rabbit F(ab) ₂ Anti-Sheep IgG(H+L), Human SP ads-FITC	IgG	Polyclonal	20µl
S6016-04	Rabbit F(ab) ₂ Anti-Sheep IgG(H+L), Human SP ads-AP	IgG	Polyclonal	20µl
S6016-05	Rabbit F(ab) ₂ Anti-Sheep IgG(H+L), Human SP ads-HRP	IgG	Polyclonal	20µl
S6020-02	Rabbit F(ab) ₂ Anti-Goat IgG(H+L)-FITC	IgG	Polyclonal	20µl
S6020-04	Rabbit F(ab) ₂ Anti-Goat IgG(H+L)-AP	IgG	Polyclonal	20µl
S6020-05	Rabbit F(ab) ₂ Anti-Goat IgG(H+L)-HRP	IgG	Polyclonal	20µl
S6022-05	Rabbit Anti-Goat IgG Fab-HRP	IgG	Polyclonal	20µl
S6022-08	Rabbit Anti-Goat IgG Fab-BIOT	IgG	Polyclonal	20µl
S6026-02	Rabbit F(ab) ₂ Anti-Goat IgG(H+L), Human SP ads-FITC	IgG	Polyclonal	20µl
S6026-04	Rabbit F(ab) ₂ Anti-Goat IgG(H+L), Human SP ads-AP	IgG	Polyclonal	20µl
S6026-05	Rabbit F(ab) ₂ Anti-Goat IgG(H+L), Human SP ads-HRP	IgG	Polyclonal	20µl
S6030-02	Goat Anti-Bovine IgG(H+L)-FITC	IgG	Polyclonal	20µl
S6030-04	Goat Anti-Bovine IgG(H+L)-AP	IgG	Polyclonal	20µl
S6030-05	Goat Anti-Bovine IgG(H+L)-HRP	IgG	Polyclonal	20µl
S6030-08	Goat Anti-Bovine IgG(H+L)-BIOT	IgG	Polyclonal	20µl

Secondary Antibody Promotion 2023
sample list

Cat.no.	Description	Isotype	Clone	Size
S6040-02	Goat Anti-Equine IgG(H+L)-FITC	IgG	Polyclonal	20µl
S6040-04	Goat Anti-Equine IgG(H+L)-AP	IgG	Polyclonal	20µl
S6040-05	Goat Anti-Equine IgG(H+L)-HRP	IgG	Polyclonal	20µl
S6040-08	Goat Anti-Equine IgG(H+L)-BIOT	IgG	Polyclonal	20µl
S6040-09	Goat Anti-Equine IgG(H+L)-PE	IgG	Polyclonal	20µl
S6045-02	Goat Anti-Llama IgG(H+L)-FITC	IgG	Polyclonal	20µl
S6045-04	Goat Anti-Llama IgG(H+L)-AP	IgG	Polyclonal	20µl
S6045-05	Goat Anti-Llama IgG(H+L)-HRP	IgG	Polyclonal	20µl
S6045-08	Goat Anti-Llama IgG(H+L)-BIOT	IgG	Polyclonal	20µl
S6045-30	Goat Anti-Llama IgG(H+L)-AF488	IgG	Polyclonal	20µl
S6045-31	Goat Anti-Llama IgG(H+L)-AF647	IgG	Polyclonal	20µl
S6050-02	Goat Anti-Porcine IgG(H+L)-FITC	IgG	Polyclonal	20µl
S6050-04	Goat Anti-Porcine IgG(H+L)-AP	IgG	Polyclonal	20µl
S6050-05	Goat Anti-Porcine IgG(H+L)-HRP	IgG	Polyclonal	20µl
S6050-08	Goat Anti-Porcine IgG(H+L)-BIOT	IgG	Polyclonal	20µl
S6050-09	Goat Anti-Porcine IgG(H+L)-PE	IgG	Polyclonal	20µl
S6050-30	Goat Anti-Porcine IgG(H+L)-AF488	IgG	Polyclonal	20µl
S6060-02	Goat Anti-Hamster IgG(H+L)-FITC	IgG	Polyclonal	20µl
S6060-04	Goat Anti-Hamster IgG(H+L)-AP	IgG	Polyclonal	20µl
S6060-05	Goat Anti-Hamster IgG(H+L)-HRP	IgG	Polyclonal	20µl
S6060-08	Goat Anti-Hamster IgG(H+L)-BIOT	IgG	Polyclonal	20µl
S6061-02	Goat Anti-Hamster IgG(H+L), Mouse/Rat ads-FITC	IgG	Polyclonal	20µl
S6061-04	Goat Anti-Hamster IgG(H+L), Mouse/Rat ads-AP	IgG	Polyclonal	20µl
S6061-05	Goat Anti-Hamster IgG(H+L), Mouse/Rat ads-HRP	IgG	Polyclonal	20µl
S6061-08	Goat Anti-Hamster IgG(H+L), Mouse/Rat ads-BIOT	IgG	Polyclonal	20µl
S6061-09	Goat Anti-Hamster IgG(H+L), Mouse/Rat ads-PE	IgG	Polyclonal	20µl
S6061-11	Goat Anti-Hamster IgG(H+L), Mouse/Rat ads-APC	IgG	Polyclonal	20µl
S6061-17	Goat Anti-Hamster IgG(H+L), Mouse/Rat ads-PE/CY7	IgG	Polyclonal	20µl
S6062-02	Goat F(ab) ₂ Anti-Hamster IgG(H+L), Mouse/Rat ads-FITC	IgG	Polyclonal	20µl
S6062-08	Goat F(ab) ₂ Anti-Hamster IgG(H+L), Mouse/Rat ads-BIOT	IgG	Polyclonal	20µl
S6062-09	Goat F(ab) ₂ Anti-Hamster IgG(H+L), Mouse/Rat ads-PE	IgG	Polyclonal	20µl
S6070-02	Goat Anti-Canine IgG(H+L)-FITC	IgG	Polyclonal	20µl
S6070-04	Goat Anti-Canine IgG(H+L)-AP	IgG	Polyclonal	20µl
S6070-05	Goat Anti-Canine IgG(H+L)-HRP	IgG	Polyclonal	20µl
S6070-08	Goat Anti-Canine IgG(H+L)-BIOT	IgG	Polyclonal	20µl
S6070-30	Goat Anti-Canine IgG(H+L)-AF488	IgG	Polyclonal	20µl
S6080-02	Goat Anti-Feline IgG(H+L)-FITC	IgG	Polyclonal	20µl
S6080-04	Goat Anti-Feline IgG(H+L)-AP	IgG	Polyclonal	20µl
S6080-05	Goat Anti-Feline IgG(H+L)-HRP	IgG	Polyclonal	20µl
S6090-02	Goat Anti-Guinea Pig IgG(H+L)-FITC	IgG	Polyclonal	20µl
S6090-04	Goat Anti-Guinea Pig IgG(H+L)-AP	IgG	Polyclonal	20µl
S6090-05	Goat Anti-Guinea Pig IgG(H+L)-HRP	IgG	Polyclonal	20µl
S6090-08	Goat Anti-Guinea Pig IgG(H+L)-BIOT	IgG	Polyclonal	20µl
S6100-02	Goat Anti-Chicken IgY(H+L)-FITC	IgG	Polyclonal	20µl
S6100-04	Goat Anti-Chicken IgY(H+L)-AP	IgG	Polyclonal	20µl
S6100-05	Goat Anti-Chicken IgY(H+L)-HRP	IgG	Polyclonal	20µl
S6100-08	Goat Anti-Chicken IgY(H+L)-BIOT	IgG	Polyclonal	20µl
S6100-09	Goat Anti-Chicken IgY(H+L)-PE	IgG	Polyclonal	20µl
S6100-30	Goat Anti-Chicken IgY(H+L)-AF488	IgG	Polyclonal	20µl
S6105-02	Goat Anti-Chicken IgY(H+L), Multi-Species SP ads-FITC	IgG	Polyclonal	20µl
S6105-05	Goat Anti-Chicken IgY(H+L), Multi-Species SP ads-HRP	IgG	Polyclonal	20µl
S6105-08	Goat Anti-Chicken IgY(H+L), Multi-Species SP ads-BIOT	IgG	Polyclonal	20µl
S6110-02	Goat Anti-Turkey IgY(H+L)-FITC	IgG	Polyclonal	20µl
S6110-04	Goat Anti-Turkey IgY(H+L)-AP	IgG	Polyclonal	20µl
S6110-05	Goat Anti-Turkey IgY(H+L)-HRP	IgG	Polyclonal	20µl
S6110-08	Goat Anti-Turkey IgY(H+L)-BIOT	IgG	Polyclonal	20µl
S6110-09	Goat Anti-Turkey IgY(H+L)-PE	IgG	Polyclonal	20µl
S6120-02	Rabbit F(ab) ₂ Anti-Mouse IgG(H+L)-FITC	IgG	Polyclonal	20µl
S6120-04	Rabbit F(ab) ₂ Anti-Mouse IgG(H+L)-AP	IgG	Polyclonal	20µl
S6120-05	Rabbit F(ab) ₂ Anti-Mouse IgG(H+L)-HRP	IgG	Polyclonal	20µl
S6120-08	Rabbit F(ab) ₂ Anti-Mouse IgG(H+L)-BIOT	IgG	Polyclonal	20µl
S6125-02	Rabbit F(ab) ₂ Anti-Mouse IgG(H+L), Human ads-FITC	IgG	Polyclonal	20µl
S6125-04	Rabbit F(ab) ₂ Anti-Mouse IgG(H+L), Human ads-AP	IgG	Polyclonal	20µl
S6125-05	Rabbit F(ab) ₂ Anti-Mouse IgG(H+L), Human ads-HRP	IgG	Polyclonal	20µl
S6125-08	Rabbit F(ab) ₂ Anti-Mouse IgG(H+L), Human ads-BIOT	IgG	Polyclonal	20µl
S6130-02	Rabbit F(ab) ₂ Anti-Rat IgG(H+L)-FITC	IgG	Polyclonal	20µl
S6130-03	Rabbit F(ab) ₂ Anti-Rat IgG(H+L)-TRITC	IgG	Polyclonal	20µl
S6130-04	Rabbit F(ab) ₂ Anti-Rat IgG(H+L)-AP	IgG	Polyclonal	20µl
S6130-05	Rabbit F(ab) ₂ Anti-Rat IgG(H+L)-HRP	IgG	Polyclonal	20µl
S6130-08	Rabbit F(ab) ₂ Anti-Rat IgG(H+L)-BIOT	IgG	Polyclonal	20µl
S6135-02	Rabbit F(ab) ₂ Anti-Rat IgG(H+L), Human ads-FITC	IgG	Polyclonal	20µl
S6135-04	Rabbit F(ab) ₂ Anti-Rat IgG(H+L), Human ads-AP	IgG	Polyclonal	20µl
S6135-05	Rabbit F(ab) ₂ Anti-Rat IgG(H+L), Human ads-HRP	IgG	Polyclonal	20µl
S6140-02	Rabbit Anti-Human IgG(H+L)-FITC	IgG	Polyclonal	20µl
S6140-04	Rabbit Anti-Human IgG(H+L)-AP	IgG	Polyclonal	20µl

Secondary Antibody Promotion 2023
sample list

Cat.no.	Description	Isotype	Clone	Size
S6140-05	Rabbit Anti-Human IgG(H+L)-HRP	IgG	Polyclonal	20µl
S6140-08	Rabbit Anti-Human IgG(H+L)-BIOT	IgG	Polyclonal	20µl
S6140-09	Rabbit Anti-Human IgG(H+L)-PE	IgG	Polyclonal	20µl
S6145-02	Rabbit Anti-Human IgG(H+L), Mouse ads-FITC	IgG	Polyclonal	20µl
S6145-04	Rabbit Anti-Human IgG(H+L), Mouse ads-AP	IgG	Polyclonal	20µl
S6145-05	Rabbit Anti-Human IgG(H+L), Mouse ads-HRP	IgG	Polyclonal	20µl
S6145-08	Rabbit Anti-Human IgG(H+L), Mouse ads-BIOT	IgG	Polyclonal	20µl
S6150-02	Rabbit Anti-Sheep IgG(H+L)-FITC	IgG	Polyclonal	20µl
S6150-03	Rabbit Anti-Sheep IgG(H+L)-TRITC	IgG	Polyclonal	20µl
S6150-04	Rabbit Anti-Sheep IgG(H+L)-AP	IgG	Polyclonal	20µl
S6150-05	Rabbit Anti-Sheep IgG(H+L)-HRP	IgG	Polyclonal	20µl
S6150-08	Rabbit Anti-Sheep IgG(H+L)-BIOT	IgG	Polyclonal	20µl
S6156-02	Rabbit Anti-Sheep IgG(H+L), Human SP ads-FITC	IgG	Polyclonal	20µl
S6156-04	Rabbit Anti-Sheep IgG(H+L), Human SP ads-AP	IgG	Polyclonal	20µl
S6156-05	Rabbit Anti-Sheep IgG(H+L), Human SP ads-HRP	IgG	Polyclonal	20µl
S6156-08	Rabbit Anti-Sheep IgG(H+L), Human SP ads-BIOT	IgG	Polyclonal	20µl
S6157-02	Mouse Anti-Goat IgG Fc-FITC	IgG _{2b}	SB115h	20µl
S6157-05	Mouse Anti-Goat IgG Fc-HRP	IgG _{2b}	SB115h	20µl
S6157-08	Mouse Anti-Goat IgG Fc-BIOT	IgG _{2b}	SB115h	20µl
S6157-09	Mouse Anti-Goat IgG Fc-PE	IgG _{2b}	SB115h	20µl
S6158-02	Mouse Anti-Goat IgG Fc-FITC	IgG _{2a}	SB115d	20µl
S6158-05	Mouse Anti-Goat IgG Fc-HRP	IgG _{2a}	SB115d	20µl
S6158-08	Mouse Anti-Goat IgG Fc-BIOT	IgG _{2a}	SB115d	20µl
S6158-30	Mouse Anti-Goat IgG Fc-AF488	IgG _{2a}	SB115d	20µl
S6158-31	Mouse Anti-Goat IgG Fc-AF647	IgG _{2a}	SB115d	20µl
S6158-32	Mouse Anti-Goat IgG Fc-AF555	IgG _{2a}	SB115d	20µl
S6160-02	Rabbit Anti-Goat IgG(H+L)-FITC	IgG	Polyclonal	20µl
S6160-03	Rabbit Anti-Goat IgG(H+L)-TRITC	IgG	Polyclonal	20µl
S6160-04	Rabbit Anti-Goat IgG(H+L)-AP	IgG	Polyclonal	20µl
S6160-05	Rabbit Anti-Goat IgG(H+L)-HRP	IgG	Polyclonal	20µl
S6160-08	Rabbit Anti-Goat IgG(H+L)-BIOT	IgG	Polyclonal	20µl
S6160-30	Rabbit Anti-Goat IgG(H+L)-AF488	IgG	Polyclonal	20µl
S6160-31	Rabbit Anti-Goat IgG(H+L)-AF647	IgG	Polyclonal	20µl
S6160-32	Rabbit Anti-Goat IgG(H+L)-AF555	IgG	Polyclonal	20µl
S6163-02	Rabbit Anti-Goat IgG Fc-FITC	IgG	Polyclonal	20µl
S6163-05	Rabbit Anti-Goat IgG Fc-HRP	IgG	Polyclonal	20µl
S6164-02	Rabbit Anti-Goat IgG(H+L), Human SP ads-FITC	IgG	Polyclonal	20µl
S6164-04	Rabbit Anti-Goat IgG(H+L), Human SP ads-AP	IgG	Polyclonal	20µl
S6164-05	Rabbit Anti-Goat IgG(H+L), Human SP ads-HRP	IgG	Polyclonal	20µl
S6164-08	Rabbit Anti-Goat IgG(H+L), Human SP ads-BIOT	IgG	Polyclonal	20µl
S6170-02	Rabbit Anti-Mouse IgG(H+L)-FITC	IgG	Polyclonal	20µl
S6170-04	Rabbit Anti-Mouse IgG(H+L)-AP	IgG	Polyclonal	20µl
S6170-05	Rabbit Anti-Mouse IgG(H+L)-HRP	IgG	Polyclonal	20µl
S6170-08	Rabbit Anti-Mouse IgG(H+L)-BIOT	IgG	Polyclonal	20µl
S6170-30	Rabbit Anti-Mouse IgG(H+L)-AF488	IgG	Polyclonal	20µl
S6170-31	Rabbit Anti-Mouse IgG(H+L)-AF647	IgG	Polyclonal	20µl
S6175-02	Rabbit Anti-Mouse IgG(H+L), Human ads-FITC	IgG	Polyclonal	20µl
S6175-04	Rabbit Anti-Mouse IgG(H+L), Human ads-AP	IgG	Polyclonal	20µl
S6175-05	Rabbit Anti-Mouse IgG(H+L), Human ads-HRP	IgG	Polyclonal	20µl
S6175-08	Rabbit Anti-Mouse IgG(H+L), Human ads-BIOT	IgG	Polyclonal	20µl
S6175-09	Rabbit Anti-Mouse IgG(H+L), Human ads-PE	IgG	Polyclonal	20µl
S6175-31	Rabbit Anti-Mouse IgG(H+L), Human ads-AF647	IgG	Polyclonal	20µl
S6180-02	Rabbit Anti-Rat IgG(H+L)-FITC	IgG	Polyclonal	20µl
S6180-04	Rabbit Anti-Rat IgG(H+L)-AP	IgG	Polyclonal	20µl
S6180-05	Rabbit Anti-Rat IgG(H+L)-HRP	IgG	Polyclonal	20µl
S6180-08	Rabbit Anti-Rat IgG(H+L)-BIOT	IgG	Polyclonal	20µl
S6185-02	Rabbit Anti-Rat IgG(H+L), Human ads-FITC	IgG	Polyclonal	20µl
S6185-04	Rabbit Anti-Rat IgG(H+L), Human ads-AP	IgG	Polyclonal	20µl
S6185-05	Rabbit Anti-Rat IgG(H+L), Human ads-HRP	IgG	Polyclonal	20µl
S6185-08	Rabbit Anti-Rat IgG(H+L), Human ads-BIOT	IgG	Polyclonal	20µl
S6200-02	Goat Anti-Rhesus IgG(H+L)-FITC	IgG	Polyclonal	20µl
S6200-03	Goat Anti-Rhesus IgG(H+L)-TRITC	IgG	Polyclonal	20µl
S6200-04	Goat Anti-Rhesus IgG(H+L)-AP	IgG	Polyclonal	20µl
S6200-05	Goat Anti-Rhesus IgG(H+L)-HRP	IgG	Polyclonal	20µl
S6200-08	Goat Anti-Rhesus IgG(H+L)-BIOT	IgG	Polyclonal	20µl
S6200-09	Goat Anti-Rhesus IgG(H+L)-PE	IgG	Polyclonal	20µl
S6200-31	Goat Anti-Rhesus IgG(H+L)-AF647	IgG	Polyclonal	20µl
S6211-05	Rabbit Anti-Hamster IgG(H+L), Mouse ads-HRP	IgG	Polyclonal	20µl
S6211-09	Rabbit Anti-Hamster IgG(H+L), Mouse ads-PE	IgG	Polyclonal	20µl
S6215-02	Rabbit Anti-Hamster IgG(H+L), Mouse/Rat SP ads-FITC	IgG	Polyclonal	20µl
S6215-05	Rabbit Anti-Hamster IgG(H+L), Mouse/Rat SP ads-HRP	IgG	Polyclonal	20µl
S6215-08	Rabbit Anti-Hamster IgG(H+L), Mouse/Rat SP ads-BIOT	IgG	Polyclonal	20µl
S6215-09	Rabbit Anti-Hamster IgG(H+L), Mouse/Rat SP ads-PE	IgG	Polyclonal	20µl
S6300-02	Swine Anti-Goat IgG(H+L), Human/Rat/Mouse SP ads-FITC	IgG	Polyclonal	20µl
S6300-04	Swine Anti-Goat IgG(H+L), Human/Rat/Mouse SP ads-AP	IgG	Polyclonal	20µl

Secondary Antibody Promotion 2023
sample list

Cat.no.	Description	Isotype	Clone	Size
S6300-05	Swine Anti-Goat IgG(H+L), Human/Rat/Mouse SP ads-HRP	IgG	Polyclonal	20µl
S6300-08	Swine Anti-Goat IgG(H+L), Human/Rat/Mouse SP ads-BIOT	IgG	Polyclonal	20µl
S6300-09	Swine Anti-Goat IgG(H+L), Human/Rat/Mouse SP ads-PE	IgG	Polyclonal	20µl
S6300-31	Swine Anti-Goat IgG(H+L), Human/Rat/Mouse SP ads-AF647	IgG	Polyclonal	20µl
S6300-32	Swine Anti-Goat IgG(H+L), Human/Rat/Mouse SP ads-AF555	IgG	Polyclonal	20µl
S6310-05	Swine Anti-Rabbit IgG(H+L), Mouse/Rat/Human/Goat SP ads-HRP	IgG	Polyclonal	20µl
S6311-02	Swine Anti-Rabbit IgG(H+L)-FITC	IgG	Polyclonal	20µl
S6311-04	Swine Anti-Rabbit IgG(H+L)-AP	IgG	Polyclonal	20µl
S6311-05	Swine Anti-Rabbit IgG(H+L)-HRP	IgG	Polyclonal	20µl
S6312-04	Swine Anti-Rabbit IgG(H+L), Multi-Species SP ads-AP	IgG	Polyclonal	20µl
S6312-05	Swine Anti-Rabbit IgG(H+L), Multi-Species SP ads-HRP	IgG	Polyclonal	20µl
S6312-08	Swine Anti-Rabbit IgG(H+L), Multi-Species SP ads-BIOT	IgG	Polyclonal	20µl
S6400-04	Sheep Anti-FITC-AP	IgG	Polyclonal	20µl
S6400-05	Sheep Anti-FITC-HRP	IgG	Polyclonal	20µl
S6400-08	Sheep Anti-FITC-BIOT	IgG	Polyclonal	20µl
S6400-09	Sheep Anti-FITC-PE	IgG	Polyclonal	20µl
S6402-04	Sheep F(ab) ₂ Anti-FITC-AP	IgG	Polyclonal	20µl
S6402-05	Sheep F(ab) ₂ Anti-FITC-HRP	IgG	Polyclonal	20µl
S6404-02	Mouse Anti-Biotin-FITC	IgG ₁	SB58c	20µl
S6404-04	Mouse Anti-Biotin-AP	IgG ₁	SB58c	20µl
S6404-05	Mouse Anti-Biotin-HRP	IgG ₁	SB58c	20µl
S6404-09	Mouse Anti-Biotin-PE	IgG ₁	SB58c	20µl
S6410-02	Donkey Anti-Mouse IgG(H+L), Human SP ads-FITC	IgG	Polyclonal	20µl
S6410-04	Donkey Anti-Mouse IgG(H+L), Human SP ads-AP	IgG	Polyclonal	20µl
S6410-05	Donkey Anti-Mouse IgG(H+L), Human SP ads-HRP	IgG	Polyclonal	20µl
S6410-08	Donkey Anti-Mouse IgG(H+L), Human SP ads-BIOT	IgG	Polyclonal	20µl
S6410-09	Donkey Anti-Mouse IgG(H+L), Human SP ads-PE	IgG	Polyclonal	20µl
S6410-15	Donkey Anti-Mouse IgG(H+L), Human SP ads-CY5	IgG	Polyclonal	20µl
S6410-30	Donkey Anti-Mouse IgG(H+L), Human SP ads-AF488	IgG	Polyclonal	20µl
S6410-31	Donkey Anti-Mouse IgG(H+L), Human SP ads-AF647	IgG	Polyclonal	20µl
S6410-32	Donkey Anti-Mouse IgG(H+L), Human SP ads-AF555	IgG	Polyclonal	20µl
S6411-02	Donkey Anti-Mouse IgG(H+L)-FITC	IgG	Polyclonal	20µl
S6411-04	Donkey Anti-Mouse IgG(H+L)-AP	IgG	Polyclonal	20µl
S6411-05	Donkey Anti-Mouse IgG(H+L)-HRP	IgG	Polyclonal	20µl
S6411-08	Donkey Anti-Mouse IgG(H+L)-BIOT	IgG	Polyclonal	20µl
S6411-30	Donkey Anti-Mouse IgG(H+L)-AF488	IgG	Polyclonal	20µl
S6411-31	Donkey Anti-Mouse IgG(H+L)-AF647	IgG	Polyclonal	20µl
S6411-32	Donkey Anti-Mouse IgG(H+L)-AF555	IgG	Polyclonal	20µl
S6412-02	Donkey F(ab) ₂ Anti-Mouse IgG(H+L)-FITC	IgG	Polyclonal	20µl
S6412-05	Donkey F(ab) ₂ Anti-Mouse IgG(H+L)-HRP	IgG	Polyclonal	20µl
S6412-08	Donkey F(ab) ₂ Anti-Mouse IgG(H+L)-BIOT	IgG	Polyclonal	20µl
S6412-30	Donkey F(ab) ₂ Anti-Mouse IgG(H+L)-AF488	IgG	Polyclonal	20µl
S6412-31	Donkey F(ab) ₂ Anti-Mouse IgG(H+L)-AF647	IgG	Polyclonal	20µl
S6412-32	Donkey F(ab) ₂ Anti-Mouse IgG(H+L)-AF555	IgG	Polyclonal	20µl
S6415-02	Donkey Anti-Mouse IgG(H+L), Multi-Species SP ads-FITC	IgG	Polyclonal	20µl
S6415-04	Donkey Anti-Mouse IgG(H+L), Multi-Species SP ads-AP	IgG	Polyclonal	20µl
S6415-05	Donkey Anti-Mouse IgG(H+L), Multi-Species SP ads-HRP	IgG	Polyclonal	20µl
S6415-08	Donkey Anti-Mouse IgG(H+L), Multi-Species SP ads-BIOT	IgG	Polyclonal	20µl
S6415-15	Donkey Anti-Mouse IgG(H+L), Multi-Species SP ads-CY5	IgG	Polyclonal	20µl
S6415-30	Donkey Anti-Mouse IgG(H+L), Multi-Species SP ads-AF488	IgG	Polyclonal	20µl
S6415-31	Donkey Anti-Mouse IgG(H+L), Multi-Species SP ads-AF647	IgG	Polyclonal	20µl
S6415-32	Donkey Anti-Mouse IgG(H+L), Multi-Species SP ads-AF555	IgG	Polyclonal	20µl
S6420-02	Donkey Anti-Goat IgG(H+L), Mouse/Rat SP ads-FITC	IgG	Polyclonal	20µl
S6420-04	Donkey Anti-Goat IgG(H+L), Mouse/Rat SP ads-AP	IgG	Polyclonal	20µl
S6420-05	Donkey Anti-Goat IgG(H+L), Mouse/Rat SP ads-HRP	IgG	Polyclonal	20µl
S6420-08	Donkey Anti-Goat IgG(H+L), Mouse/Rat SP ads-BIOT	IgG	Polyclonal	20µl
S6420-09	Donkey Anti-Goat IgG(H+L), Mouse/Rat SP ads-PE	IgG	Polyclonal	20µl
S6420-15	Donkey Anti-Goat IgG(H+L), Mouse/Rat SP ads-CY5	IgG	Polyclonal	20µl
S6420-30	Donkey Anti-Goat IgG(H+L), Mouse/Rat SP ads-AF488	IgG	Polyclonal	20µl
S6420-31	Donkey Anti-Goat IgG(H+L), Mouse/Rat SP ads-AF647	IgG	Polyclonal	20µl
S6420-32	Donkey Anti-Goat IgG(H+L), Mouse/Rat SP ads-AF555	IgG	Polyclonal	20µl
S6421-02	Donkey F(ab) ₂ Anti-Goat IgG(H+L), Mouse/Rat SP ads-FITC	IgG	Polyclonal	20µl
S6421-04	Donkey F(ab) ₂ Anti-Goat IgG(H+L), Mouse/Rat SP ads-AP	IgG	Polyclonal	20µl
S6421-05	Donkey F(ab) ₂ Anti-Goat IgG(H+L), Mouse/Rat SP ads-HRP	IgG	Polyclonal	20µl
S6421-08	Donkey F(ab) ₂ Anti-Goat IgG(H+L), Mouse/Rat SP ads-BIOT	IgG	Polyclonal	20µl
S6421-30	Donkey F(ab) ₂ Anti-Goat IgG(H+L), Mouse/Rat SP ads-AF488	IgG	Polyclonal	20µl
S6421-31	Donkey F(ab) ₂ Anti-Goat IgG(H+L), Mouse/Rat SP ads-AF647	IgG	Polyclonal	20µl
S6421-32	Donkey F(ab) ₂ Anti-Goat IgG(H+L), Mouse/Rat SP ads-AF555	IgG	Polyclonal	20µl
S6425-04	Donkey Anti-Goat IgG(H+L), Multi-Species SP ads-AP	IgG	Polyclonal	20µl
S6425-05	Donkey Anti-Goat IgG(H+L), Multi-Species SP ads-HRP	IgG	Polyclonal	20µl
S6425-08	Donkey Anti-Goat IgG(H+L), Multi-Species SP ads-BIOT	IgG	Polyclonal	20µl
S6425-15	Donkey Anti-Goat IgG(H+L), Multi-Species SP ads-CY5	IgG	Polyclonal	20µl
S6425-30	Donkey Anti-Goat IgG(H+L), Multi-Species SP ads-AF488	IgG	Polyclonal	20µl
S6425-31	Donkey Anti-Goat IgG(H+L), Multi-Species SP ads-AF647	IgG	Polyclonal	20µl
S6425-32	Donkey Anti-Goat IgG(H+L), Multi-Species SP ads-AF555	IgG	Polyclonal	20µl

Secondary Antibody Promotion 2023
sample list

Cat.no.	Description	Isotype	Clone	Size
S6430-02	Donkey Anti-Rat IgG(H+L), Mouse SP ads-FITC	IgG	Polyclonal	20µl
S6430-04	Donkey Anti-Rat IgG(H+L), Mouse SP ads-AP	IgG	Polyclonal	20µl
S6430-05	Donkey Anti-Rat IgG(H+L), Mouse SP ads-HRP	IgG	Polyclonal	20µl
S6430-08	Donkey Anti-Rat IgG(H+L), Mouse SP ads-BIOT	IgG	Polyclonal	20µl
S6430-09	Donkey Anti-Rat IgG(H+L), Mouse SP ads-PE	IgG	Polyclonal	20µl
S6430-30	Donkey Anti-Rat IgG(H+L), Mouse SP ads-AF488	IgG	Polyclonal	20µl
S6430-31	Donkey Anti-Rat IgG(H+L), Mouse SP ads-AF647	IgG	Polyclonal	20µl
S6430-32	Donkey Anti-Rat IgG(H+L), Mouse SP ads-AF555	IgG	Polyclonal	20µl
S6431-04	Donkey F(ab) ₂ Anti-Rat IgG(H+L), Mouse SP ads-AP	IgG	Polyclonal	20µl
S6431-05	Donkey F(ab) ₂ Anti-Rat IgG(H+L), Mouse SP ads-HRP	IgG	Polyclonal	20µl
S6431-08	Donkey F(ab) ₂ Anti-Rat IgG(H+L), Mouse SP ads-BIOT	IgG	Polyclonal	20µl
S6431-30	Donkey F(ab) ₂ Anti-Rat IgG(H+L), Mouse SP ads-AF488	IgG	Polyclonal	20µl
S6431-31	Donkey F(ab) ₂ Anti-Rat IgG(H+L), Mouse SP ads-AF647	IgG	Polyclonal	20µl
S6431-32	Donkey F(ab) ₂ Anti-Rat IgG(H+L), Mouse SP ads-AF555	IgG	Polyclonal	20µl
S6440-02	Donkey Anti-Rabbit IgG(H+L), Mouse/Rat/Human SP ads-FITC	IgG	Polyclonal	20µl
S6440-04	Donkey Anti-Rabbit IgG(H+L), Mouse/Rat/Human SP ads-AP	IgG	Polyclonal	20µl
S6440-05	Donkey Anti-Rabbit IgG(H+L), Mouse/Rat/Human SP ads-HRP	IgG	Polyclonal	20µl
S6440-08	Donkey Anti-Rabbit IgG(H+L), Mouse/Rat/Human SP ads-BIOT	IgG	Polyclonal	20µl
S6440-09	Donkey Anti-Rabbit IgG(H+L), Mouse/Rat/Human SP ads-PE	IgG	Polyclonal	20µl
S6440-15	Donkey Anti-Rabbit IgG(H+L), Mouse/Rat/Human SP ads-CY5	IgG	Polyclonal	20µl
S6440-30	Donkey Anti-Rabbit IgG(H+L), Mouse/Rat/Human SP ads-AF488	IgG	Polyclonal	20µl
S6440-31	Donkey Anti-Rabbit IgG(H+L), Mouse/Rat/Human SP ads-AF647	IgG	Polyclonal	20µl
S6440-32	Donkey Anti-Rabbit IgG(H+L), Mouse/Rat/Human SP ads-AF555	IgG	Polyclonal	20µl
S6441-02	Donkey Anti-Rabbit IgG(H+L)-FITC	IgG	Polyclonal	20µl
S6441-04	Donkey Anti-Rabbit IgG(H+L)-AP	IgG	Polyclonal	20µl
S6441-05	Donkey Anti-Rabbit IgG(H+L)-HRP	IgG	Polyclonal	20µl
S6441-08	Donkey Anti-Rabbit IgG(H+L)-BIOT	IgG	Polyclonal	20µl
S6441-30	Donkey Anti-Rabbit IgG(H+L)-AF488	IgG	Polyclonal	20µl
S6441-31	Donkey Anti-Rabbit IgG(H+L)-AF647	IgG	Polyclonal	20µl
S6441-32	Donkey Anti-Rabbit IgG(H+L)-AF555	IgG	Polyclonal	20µl
S6442-04	Donkey Anti-Rabbit IgG(H+L), Multi-Species SP ads-AP	IgG	Polyclonal	20µl
S6442-05	Donkey Anti-Rabbit IgG(H+L), Multi-Species SP ads-HRP	IgG	Polyclonal	20µl
S6442-08	Donkey Anti-Rabbit IgG(H+L), Multi-Species SP ads-BIOT	IgG	Polyclonal	20µl
S6442-15	Donkey Anti-Rabbit IgG(H+L), Multi-Species SP ads-CY5	IgG	Polyclonal	20µl
S6442-30	Donkey Anti-Rabbit IgG(H+L), Multi-Species SP ads-AF488	IgG	Polyclonal	20µl
S6442-31	Donkey Anti-Rabbit IgG(H+L), Multi-Species SP ads-AF647	IgG	Polyclonal	20µl
S6442-32	Donkey Anti-Rabbit IgG(H+L), Multi-Species SP ads-AF555	IgG	Polyclonal	20µl
S6444-04	Donkey F(ab) ₂ Anti-Rabbit IgG(H+L), Multi-Species SP ads-AP	IgG	Polyclonal	20µl
S6444-05	Donkey F(ab) ₂ Anti-Rabbit IgG(H+L), Multi-Species SP ads-HRP	IgG	Polyclonal	20µl
S6444-08	Donkey F(ab) ₂ Anti-Rabbit IgG(H+L), Multi-Species SP ads-BIOT	IgG	Polyclonal	20µl
S6444-30	Donkey F(ab) ₂ Anti-Rabbit IgG(H+L), Multi-Species SP ads-AF488	IgG	Polyclonal	20µl
S6444-31	Donkey F(ab) ₂ Anti-Rabbit IgG(H+L), Multi-Species SP ads-AF647	IgG	Polyclonal	20µl
S6444-32	Donkey F(ab) ₂ Anti-Rabbit IgG(H+L), Multi-Species SP ads-AF555	IgG	Polyclonal	20µl
S6446-04	Donkey F(ab) ₂ Anti-Rabbit IgG(H+L), Mouse/Rat/Human SP ads-AP	IgG	Polyclonal	20µl
S6446-05	Donkey F(ab) ₂ Anti-Rabbit IgG(H+L), Mouse/Rat/Human SP ads-HRP	IgG	Polyclonal	20µl
S6446-08	Donkey F(ab) ₂ Anti-Rabbit IgG(H+L), Mouse/Rat/Human SP ads-BIOT	IgG	Polyclonal	20µl
S6446-30	Donkey F(ab) ₂ Anti-Rabbit IgG(H+L), Mouse/Rat/Human SP ads-AF488	IgG	Polyclonal	20µl
S6446-31	Donkey F(ab) ₂ Anti-Rabbit IgG(H+L), Mouse/Rat/Human SP ads-AF647	IgG	Polyclonal	20µl
S6446-32	Donkey F(ab) ₂ Anti-Rabbit IgG(H+L), Mouse/Rat/Human SP ads-AF555	IgG	Polyclonal	20µl
S6460-04	Donkey Anti-Goat IgG Fc-AP	IgG	Polyclonal	20µl
S6460-05	Donkey Anti-Goat IgG Fc-HRP	IgG	Polyclonal	20µl
S6460-08	Donkey Anti-Goat IgG Fc-BIOT	IgG	Polyclonal	20µl
S6460-30	Donkey Anti-Goat IgG Fc-AF488	IgG	Polyclonal	20µl
S6460-31	Donkey Anti-Goat IgG Fc-AF647	IgG	Polyclonal	20µl
S6460-32	Donkey Anti-Goat IgG Fc-AF555	IgG	Polyclonal	20µl
S6800-02	Goat Anti-Donkey IgG(H+L)-FITC	IgG	Polyclonal	20µl
S6800-05	Goat Anti-Donkey IgG(H+L)-HRP	IgG	Polyclonal	20µl
S6800-08	Goat Anti-Donkey IgG(H+L)-BIOT	IgG	Polyclonal	20µl
S6900-04	Sheep Anti-Rat IgG(H+L)-AP	IgG	Polyclonal	20µl
S6900-05	Sheep Anti-Rat IgG(H+L)-HRP	IgG	Polyclonal	20µl
S6900-08	Sheep Anti-Rat IgG(H+L)-BIOT	IgG	Polyclonal	20µl
S7105-02	Streptavidin-FITC		Polyclonal	20µl
S7105-03	Streptavidin-TRITC		Polyclonal	20µl
S7105-04	Streptavidin-AP		Polyclonal	20µl
S7105-05	Streptavidin-HRP		Polyclonal	20µl
S7105-06	Streptavidin-BGAL		Polyclonal	20µl
S7105-07	Streptavidin-TXRD		Polyclonal	20µl
S7105-09	Streptavidin-PE		Polyclonal	20µl
S7105-10	Streptavidin-PE/TXRD		Polyclonal	20µl
S7105-11	Streptavidin-APC		Polyclonal	20µl
S7105-12	Streptavidin-CY3		Polyclonal	20µl
S7105-13	Streptavidin-SPRD		Polyclonal	20µl
S7105-15	Streptavidin-CY5		Polyclonal	20µl
S7105-17	Streptavidin-PE/CY7		Polyclonal	20µl
S7105-19	Streptavidin-APC/CY7		Polyclonal	20µl

Secondary Antibody Promotion 2023
sample list

Cat.no.	Description	Isotype	Clone	Size
S7200-02	Neutralite Avidin-FITC		Polyclonal	20µl
S7200-03	Neutralite Avidin-TRITC		Polyclonal	20µl
S7200-04	Neutralite Avidin-AP		Polyclonal	20µl
S7200-05	Neutralite Avidin-HRP		Polyclonal	20µl
S7200-06	Neutralite Avidin-BGAL		Polyclonal	20µl
S7200-07	Neutralite Avidin-TXRD		Polyclonal	20µl
S7200-09	Neutralite Avidin-PE		Polyclonal	20µl
S7200-11	Neutralite Avidin-APC		Polyclonal	20µl
S7200-12	Neutralite Avidin-CY3		Polyclonal	20µl
S7200-13	Neutralite Avidin-SPRD		Polyclonal	20µl
S7200-15	Neutralite Avidin-CY5		Polyclonal	20µl
S9020-02	Mouse Anti-Human IgM-FITC	IgG ₁	SA-DA4	20µl
S9020-04	Mouse Anti-Human IgM-AP	IgG ₁	SA-DA4	20µl
S9020-05	Mouse Anti-Human IgM-HRP	IgG ₁	SA-DA4	20µl
S9020-08	Mouse Anti-Human IgM-BIOT	IgG ₁	SA-DA4	20µl
S9020-09	Mouse Anti-Human IgM-PE	IgG ₁	SA-DA4	20µl
S9020-11	Mouse Anti-Human IgM-APC	IgG ₁	SA-DA4	20µl
S9020-30	Mouse Anti-Human IgM-AF488	IgG ₁	SA-DA4	20µl
S9020-31	Mouse Anti-Human IgM-AF647	IgG ₁	SA-DA4	20µl
S9020-32	Mouse Anti-Human IgM-AF555	IgG ₁	SA-DA4	20µl
S9022-02	Mouse Anti-Human IgM-FITC	IgG ₃	UHB	20µl
S9022-04	Mouse Anti-Human IgM-AP	IgG ₃	UHB	20µl
S9022-05	Mouse Anti-Human IgM-HRP	IgG ₃	UHB	20µl
S9022-08	Mouse Anti-Human IgM-BIOT	IgG ₃	UHB	20µl
S9022-09	Mouse Anti-Human IgM-PE	IgG ₃	UHB	20µl
S9022-13	Mouse Anti-Human IgM-SPRD	IgG ₃	UHB	20µl
S9040-02	Mouse Anti-Human IgG Fc-FITC	IgG ₁	JDC-10	20µl
S9040-04	Mouse Anti-Human IgG Fc-AP	IgG ₁	JDC-10	20µl
S9040-05	Mouse Anti-Human IgG Fc-HRP	IgG ₁	JDC-10	20µl
S9040-08	Mouse Anti-Human IgG Fc-BIOT	IgG ₁	JDC-10	20µl
S9040-09	Mouse Anti-Human IgG Fc-PE	IgG ₁	JDC-10	20µl
S9040-30	Mouse Anti-Human IgG Fc-AF488	IgG ₁	JDC-10	20µl
S9040-31	Mouse Anti-Human IgG Fc-AF647	IgG ₁	JDC-10	20µl
S9040-32	Mouse Anti-Human IgG Fc-AF555	IgG ₁	JDC-10	20µl
S9042-02	Mouse Anti-Human IgG Fc-FITC	IgG _{2b}	H2	20µl
S9042-04	Mouse Anti-Human IgG Fc-AP	IgG _{2b}	H2	20µl
S9042-05	Mouse Anti-Human IgG Fc-HRP	IgG _{2b}	H2	20µl
S9042-08	Mouse Anti-Human IgG Fc-BIOT	IgG _{2b}	H2	20µl
S9042-09	Mouse Anti-Human IgG Fc-PE	IgG _{2b}	H2	20µl
S9042-11	Mouse Anti-Human IgG Fc-APC	IgG _{2b}	H2	20µl
S9042-13	Mouse Anti-Human IgG Fc-SPRD	IgG _{2b}	H2	20µl
S9042-15	Mouse Anti-Human IgG Fc-CY5	IgG _{2b}	H2	20µl
S9042-30	Mouse Anti-Human IgG Fc-AF488	IgG _{2b}	H2	20µl
S9042-31	Mouse Anti-Human IgG Fc-AF647	IgG _{2b}	H2	20µl
S9042-32	Mouse Anti-Human IgG Fc-AF555	IgG _{2b}	H2	20µl
S9180-02	Mouse Anti-Human Lambda-FITC	IgG ₁	JDC-12	20µl
S9180-04	Mouse Anti-Human Lambda-AP	IgG ₁	JDC-12	20µl
S9180-05	Mouse Anti-Human Lambda-HRP	IgG ₁	JDC-12	20µl
S9180-08	Mouse Anti-Human Lambda-BIOT	IgG ₁	JDC-12	20µl
S9180-09	Mouse Anti-Human Lambda-PE	IgG ₁	JDC-12	20µl
S9180-30	Mouse Anti-Human Lambda-AF488	IgG ₁	JDC-12	20µl
S9230-02	Mouse Anti-Human Kappa-FITC	IgG ₁	SB81a	20µl
S9230-04	Mouse Anti-Human Kappa-AP	IgG ₁	SB81a	20µl
S9230-05	Mouse Anti-Human Kappa-HRP	IgG ₁	SB81a	20µl
S9230-08	Mouse Anti-Human Kappa-BIOT	IgG ₁	SB81a	20µl
S9230-09	Mouse Anti-Human Kappa-PE	IgG ₁	SB81a	20µl
S9230-11	Mouse Anti-Human Kappa-APC	IgG ₁	SB81a	20µl
S9230-30	Mouse Anti-Human Kappa-AF488	IgG ₁	SB81a	20µl
S9230-31	Mouse Anti-Human Kappa-AF647	IgG ₁	SB81a	20µl
S9230-32	Mouse Anti-Human Kappa-AF555	IgG ₁	SB81a	20µl